

# Remote Controlled R/F System FLEXAVISION SF package Cassette (CR) SFD system



Remote Controlled R/F System  $FLEXAVISION^{^{\rm TM}}$  SF package Cassette (CR) SFD system

# **The Multi-Purpose System**

As examinations become increasingly diversified, R/F systems used typically for gastrointestinal examinations must now provide expanded functionality and improved performance to accommodate a wider range of examinations. The "FLEX" in "FLEXAVISION" stands for "flexibility". True to its name, FLEXAVISION is an R/F system that meets the needs of the modern age by responding to the need for multi-purpose examinations.

.

4

7

8

8

9


10

Entire System Designed to Ensure High Image Quality Digital Network Compatibility Dose Management A Variety of Options

Easy-to-Operate, High-Throughput System

Built-in Flexibility Expands the Scope of Examinations

2 FLEXAVISION SF Package



# Easy to Operate, High-Throughput System

FLEXAVISION's wide range of functions supports chest imaging plus many other types of examinations, ensuring easy performance of cassette radiography.

System operation is designed to be as simple as possible and allows easy selection of the optimum conditions for each patient's body type. The functionality of this system consistently produces images of the high quality, regardless of the operator or patient.





#### Simple Setting of R/F Table

#### Imaging Chain Extension Function

Extending the imaging chain 150 cm for thoracic radiography expands the system's imaging capability to include low magnification chest radiography. The radiography conditions automatically change to compensasate for the change in source to object distance.

A selection of cassette sizes, from  $18 \times 24$  cm to  $35 \times 43$  cm, can be used to perform a wide range of examinations, from thoracic to orthopedic. Also, a xenon-filled phototimer (option) with up to 4 spot fields and minimal radiation scattering ensure that X-ray efficiency is high and radiography is stable.




3

#### Positioning and Setting of Radiography Conditions

#### "Memory Shot" Function Ensures Optimal Radiography Conditions

For spot film radiography performed during a fluoroscopy examination, the optimum radiography conditions for radiography of the region of interest, is instantly and automatically set from the required fluoroscopic exposure condition. Additionally, the systems APR function allows the presetting of up to 12 exposure conditions. This allows for quick and optimal imaging of common procedures.



## Simple Cassette Mounting

#### Wide Selection of Cassette Sizes

#### Front-Loading Design

A simple cassette placement mechanism results in a quick and reliable cassette insertion and removal process at the SFD.

The cassette automatically retracts out of the FOV allowing fluoroscopy to be performed.






lsers can create rograms for any rocedure, based on 2 types of radiography echnique. Radiography under Optimum Conditions and Timing

#### **Cassette Spot Filming Function**

Advanced technology in the Spot Film Device gives precise radiography timing. Cassette movement from the "park" position to the "exposure" position is almost instantaneous, allowing capture of the region of interest.



# **Built-in Flexibility Expands Scope of Examinations**

The key concept of FLEXAVISION is flexibility. FLEXAVISION can accommodate a wide range of radiology applications. These applications include chest imaging, general radiography, orthopaedic imaging, urology procedures and gastrointestinal examinations.

Advanced functions incorporated into the system design give the operator a system that is quick and easy to use for general radiography through to the more specialized procedures.

#### **Spot Filming Function**

#### Multiple Cassette Sizes Facility

During orthopaedic examinations, radiography can be performed immediately when fluoroscopy identifies the desired projection, region or anatomy. We recommend that this function, which illustrates the merits of an R/F system, is used for orthopedic examinations.

A selection of cassette sizes, from 18 X 24 cm to 35 X 43 cm, can be used to handle a broad range of examinations, from thoracic to orthopedic. This system's phototimer(option) ensures consistent radiography and use of the 35 X 43 cm cassette allows capture of a large field in a single image.



Sub-divisional Spot Filming

90°

regions.

Examination range: 133 cm (when using  $35 \times 43$  cm cassette)



Wide coverage ensures a large examination area





**Bedside Switches** 


Tabletop and imaging chain operation switches are provided on the tabletop to allow examinations and positioning while caring for the patient.

High-speed Sub-divisional spot filming function makes it easy to perform accurately timed left/right 2-frame radiography of the esophagus and other

(Top/bottom 2-frame imaging and 4-frame imaging are available as options.)








Allows control of contrast media flow and performance of orthopaedic examinations

Lateral table movement:



Allows easy imaging of peripheral anatomy without the need to shift the patient on the tabletop

## **Entire System Designed to Ensure High Image Quality**

Since developing the world's first remote-control R/F system, Shimadzu has, with continued to provide cutting-edge technology FLEXAVISION has been packed with the functions and performance required by clinical facilities. Its components, including the R/F table, X-ray generator, X-ray tube, and image intensifier, were all developed in-house to create a unified system that produces high-quality images.

#### Easy-to-Use Operating Console (Intercom is an option.)





Designed considering the space restrictions of an operating room, a tabletop-type console that fits on a standard work desk is also available.




#### "Memory Shot" Function

Automatically and instantly set the optimum radiography conditions for the region of interest and operating procedure from the fluoroscopy conditions.

# **Digital Network Compatibility**

Improve examination efficiency by storing digital images acquired via the CR system using an electronic patient record system equipped with an image filing function and use them for observation and creating reports.





## **Dose Management**

FLEXAVISION not only provides high-definition images with the optimal image quality for each examination. It also effectively reduces the total exposure dose in pediatric, gynecological and other examinations where low dose exposures are required.

#### Lower Dose with a Removable Grid

The FLEXAVISION grid can be inserted or removed to suit the radiography application. The grid can easily be removed for pediatric, obstetrics and gynecological examinations when the radiation dose to the patient must be kept to a minimum.


#### Area Dosimeter

addon console.

An integrated area dosimeter can be fitted within collimator. The exposurer dose reading are displayed on the

#### Automatic BH Filters Switch to Suit the Examination

Three beam hardening (BH) filters are provided as standard to efficiently remove unnecessary soft X-rays that do not contribute to image quality. The optimal BH filter is automatically selected to suit the examination, so image quality is increased while exposure dose to the patient is reduced.






\* Physical DAP meter is also available as an option.

#### **Polygonal Aperture Reduces** Excessive Exposure

FLEXAVISION has polygonal aperture feature. It reduces exposure levels by blocking the unnecessary X-rays at the four corners of the field of view.



Blocks excessive exposure Polygonal aperture Conventional rectangular aperture





### ■ For General Radiographic / Gastro-intestinal studies

#### X-ray tube 180° swing unit

180° rotation of the X-ray tube/collimator easily and effectively accommodates chest examinations using a bucky stand.



#### X-ray tube swing-out unit

G.I. examination kit \*

Compression cone & Barium cup holder Allows rotating X-ray tube 37°/ 90° to perform chest radiography with a leader stand



#### Rolling step

Rolling step has the left-and-right rotational function It helps an elderly person to

rotate left and right smoothly, and helps a workflow at mass examination.



#### Compression band

Mattress

Oblique projection unit \*

Table elevation function \*, \*\* Lateral cassette holder



Compression band

## For Urogenital / Endoscopic studies

#### Drain bag

It is for receiving excrement during the examination of the urethra.



Leg supports One pair is mounted on the tabletop to support legs.



Auxiliary tabletop

The tabletop extension kit to mount at the head side .

It is useful for flexible positioning in ERCP or lower extremities' studies etc.

Dimension : 65cm×50cm(W×L)



Endoscope support Elbow Rest for Operator in urologic procedures



## Other Options

#### 2nd -tube option with CH-200M

Combination with CH-200M allows general radiographic techniques such as lateral and decubitus radiography to be performed on the R/F table.



#### Foot switch

In examination room, it is possible to control radiography and fluoroscopy by foot.





#### 2nd -tube option with FH-series

Combination with FH-series and a radiography stand enables the simple and safe performance of a range of procedures from chest and abdomen radiography to general radiography.



- Local console
- Monitor cart
- Maximum allowable load UP kit
- Direct phototimer
- Intercom for Operation Console
- DAP meter



Local Cons

\*Indicates options installed at the factory. \*\* Indicates an option by upon inquiry

#### Note 1)

Some options cannot be compatible with some system configurations. Please inquire for the details to our sales representative

Label Description : REMOTE-CONTROLLED R/F SYSTEM FLEXAVISION

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. We maintain a global network of sales, service, technical support and applications centers on six continents, and have established long-term relationships with a host of highly trained distributors located in over 100 countries. For information about Shimadzu, and to contact your local office, please visit our website at www.shimadzu.com



#### Shimadzu Corporation

#### Headquarters

1, Nishinokyo-Kuwabara-cho, Nakagyo-ku, Kyoto 604-8511, Japan https://www.shimadzu.com/med/



TÜVRheinland CERTIFIED

Shimadzu Corporation Medical Systems Division has been certified by TÜV Rheinland as a manufacturer of medical systems in compliance with ISO9001:2015 Quality Management Systems and ISO13485:2016 Medical Devices Quality Management Systems.

Remarks:

- Every value in this catalogue is a standard value, and it may vary a little from the actual at each site.
- The appearances and specifications are subject to change for reasons of improvement without notice.
- Items and components in the photos may include optional items. Please confirm with your sales representative for details.
- Certain configurations may not be available pending regulatory clearance. Contact your Shimadzu representative for information on specific configurations.
- Before operating this system, you should first thoroughly review the Instruction
- Manual.