

Quantitative Analysis of Flavonols in Tea Leaves

In collaboration with the National Agriculture and Food Research Organization,



Shimadzu Corporation has been developing a simple, quick and accurate method of analyzing functional components in agricultural and food products. This report introduces a quantitative method for flavonols analysis in tea leaves and presents the results obtained in two kinds of them.

Flavonols, a kind of polyphenols, are classified into flavonoids. Generally, to determine the content of flavonols, the glycosides are hydrolyzed to provide only the aglycone form rate. In this report, an analysis method of flavonols shown in Table 1 was developed to determine the contents of glycosides and aglycone without hydrolysis.

M. Kawashima

Table 1 Target Compounds

Compound
Kaempferol
Quercetin
Isoquercitrin (quercetin 3-glucoside)
Hyperoside (quercetin 3-galactoside)
Rutin (quercetin 3-rutinoside)
Myricetin
Myricetin 3-glucoside
Myricetin 3-galactoside

Sample Pretreatment

The extraction was performed following the conditions determined in the reference of pre-existing method¹⁾. The workflow is shown in Fig.1. The pretreatment was conducted by extracting crushed tea leaves with 80% MeOH solution and a dilution by 10 times in water.




Fig. 1 Pretreatment Workflow

Analytical Conditions

The analytical conditions were determined in the reference of the pre-existing method^{1), 2)}. The analytical conditions are shown in Table 2.

Table 2 Analytical Conditions

System	: Nexera™ X3
Column	: Shim-pack™ GIST C18 (150 mm × 4.6 mm I.D., 3 μm P/N : 227-30011-07)
Mobile phases	: A) 0.1% Formic acid in H ₂ O B) Acetonitrile
Gradient Program	: B conc. 15% (0-14.00 min) - 95% (22.01-24.00 min) -15% (24.01-30.00 min)
Flow rate	: 1.0 mL/min
Column Temp.	: 40 °C
Injection volume	: 10 μL
Detection	: UV 370 nm

Analysis Results of Standards

The linearities were determined by the standards analysis. Fig. 2 shows a representative chromatogram and Table 3 shows the dynamic range and the coefficients of determination. Good linearities were obtained with a coefficient of determination (R²) ≥ 0.997 for all compounds.




Fig. 2 Chromatogram of Standard Sample

Table 3 Linear range and Coefficient of determination (R²)

Compound	Linear range (μg/mL)	Coefficient of Determination (R ²)
Kaempferol	0.1 - 20	0.9979
Quercetin	0.1 - 2	0.9992
Isoquercitrin	0.1 - 2	0.9993
Hyperoside	0.1 - 20	0.9996
Rutin	0.1 - 20	0.9998
Myricetin	0.1 - 20	0.9978
Myricetin3-glucoside	0.1 - 20	0.9998
Myricetin3-galactoside	0.1 - 20	0.9999

■ Repeatability Test Results of Tea Leaf Extracts

Seven extracts were prepared from one kind of tea (benifuuki) and repeatability test was performed to confirm validity. Table 4 shows the results.

Table 4 Repeatability Test Results (n=7)

Compound	Repeatability (%RSD)
Kaempferol	- (< LLOQ)
Quercetin	- (< LLOQ)
Isoquercitrin	0.66
Hyperoside	3.33
Rutin	- (< LLOQ)
Myricetin	1.19
Myricetin3-galactoside	2.74
Myricetin3-galactoside	2.79

■ Quantitative Results for Tea Leaves

The extracts of two kinds of tea (Yabukita, Benifuuki) were analyzed to determine the content of flavonols. Fig. 3 shows the chromatograms and Table 5 shows the calculated content of each flavonol in tea leaves.




Fig. 3 Chromatograms of Tea Leaf Extracts

This analytical method was developed in collaboration with the National Agriculture and Food Research Organization (scheduled from April, 2019 to March, 2022) at the Collaborative Research Laboratory for Analysis of Food Functionality in Shimadzu's Healthcare R&D Center. The analytical method and analysis data presented in this report were provided by Mr. Hironori Juichi and Ms. Yayoi Ichiki, researchers at the National Agriculture and Food Research Organization.

Nexera and Shim-pack are trademarks of Shimadzu Corporation in Japan and/or other countries.

Table 5 Flavonol Content in Tea Leaves

Compound	Content (mg/g)	
	Yabukita	Benifuuki
Kaempferol	< LLOQ	< LLOQ
Quercetin	< LLOQ	< LLOQ
Isoquercitrin	0.252	0.131
Hyperoside	0.572	0.359
Rutin	1.975	1.348
Myricetin	n.d. *1	< LLOQ
Myricetin3-galactoside	0.575	0.182
Myricetin3-galactoside	1.006	0.501

*1 Not detected

■ Conclusion

- Using the Nexera series, simultaneous analysis of flavonols was performed.
- The flavonols quantification results show a difference in content depending on the kind of tea leaves.

<References>

- 1) Monobe et al. Quercetin Glycosides-rich Tea Cultivars (Camellia sinensis L.) in Japan. Food Science and Technology Research. 2015, 21 (3), p.333-340.
- 2) Nobuya Shirai. Assay of Flavonol Contents in Tea Leaves and Infusions. Nippon Shokuhin Kagaku Kogaku Kaishi. 2018, 65(7), p. 357-362.



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Shimadzu disclaims any proprietary interest in trademarks and trade names used in this publication other than its own. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.