



GC-MS GCMS-QP[™]2020 NX, GCMS-TQ[™] NX series

Simplified Analysis of Aqueous Short Chain Fatty Acids by GC/MS

Y. Nakagawa

User Benefits

- No derivatization or special consumables required
- ◆ Total run time of about 10 mins including the cooling period

Introduction

Short chain fatty acids are fatty acids with fewer than 6 carbon atoms. Such organic acids are found not only in the environment (e.g. rainwater, atmosphere), but also in the human body. They are produced by the gut microbiota and provide energy for the cells lining the colon. About 95 % of short chain fatty acids in the human body is acetic acid(C2), propionic acid(C3) and butyric acid(C4).

Derivatization is often used in an analysis for short chain fatty acids. In this article, C2-C6 fatty acids in water are analyzed without derivatization or internal standard (Fig. 1).





Materials and Methods

A 10 mM pre-mixed standard solution (MilliporeSigma: CRM46975) was diluted to 4, 10, 40, 100 μ M in water. The prepared calibration solutions are aliquoted into 2 mL vials and filled to the brim in order to suppress the volatile fatty acids escaping into headspace of the 2 mL vial.

Table 1 Instrument Configurations

GC-MS	: GCMS-QP2020 NX
Auto Injector	: AOC-20i Plus
Auto Sampler	: AOC-20s Plus
Analytical Column	: SH-WAX (60 m $ imes$ 0.25 mm l.D., df=0.5 μ m) *1

*1 P/N: 221-75894-60

Table 2 Analytical Conditions				
GC				
Injection Mode	: Split			
Injector Temp.	: 240 °C			
Split Ratio	: 5			
Carrier Gas	: Helium			
Control Mode	: Constant Linear Velocity (34.0 cm/s)			
Column Oven Temp.	: 80 °C (2 min) \rightarrow (40 °C/min) \rightarrow			
	200 °C \rightarrow (25 °C/min) \rightarrow 240 °C (2 min)			
	Total 8.60 mins			
Purge Flow Rate	: 3.0 mL/min			
Sample Injection Volume : 1 µL				
MS				
lon Source Temp.	: 200 °C			
Interface Temp.	: 240 °C			
Measurement Mode	: SIM			
Monitoring lons (m/z)	: Refer to Fig. 3			
Loop Time	: 0.30 seconds			

GC separation is crucial for a group of small compounds with the same functional group. The 10 mM pre-mixed standard was diluted several folds with DI water and the TIC chromatogram was obtained (Fig. 2).



Fig. 2 Total ion chromatogram of a diluted pre-mixed standard solution

Results

Fig. 3 below shows chromatograms obtained from 10 µM solution.



Fig. 3 10 μ M spike with ng/mL conversion in parentheses

Calibration curves (quadratic) were drawn from 4 to 100 μ M/mL. The concentrations in μ g/mL at the smallest and largest calibrator points are tabulated below.

Table 3	Smallest and	largert	calibrator	nointe	in	calibration	CURVAG
I able J	Jinanestanu	largest	cambrator	points		cambration	Cuives

Common d Name	μg/mL conversion			
Compound Name	4 µM	100 µM		
Acetic acid	0.24	6.0		
Propionic acid	0.30	7.4		
lsobutyric acid	0.35	8.8		
Butyric acid	0.35	8.8		
Isovaleric acid	0.41	10.2		
Valeric acid	0.41	10.2		
4-Methyl pentanoic acid	0.46	11.6		
Hexanoic acid	0.46	11.6		

Discussions

Repeatability was measured by 7 consecutive injections of aqueous sample at 10 µM (Table 4).

Table 4 % RSD at 10 µM (n=7)		Table 5 10 μM (n=7) concentrations (bank subtraction used for acetic and propionic acids)				
Compound Name	% RSD		Acetic Acid	Propionic Acid	lsobutyric acid	
Acetic acid	6.0					
Propionic acid	4.4	1	9.57	11.38	10.29	
Isobutyric acid	3.9	2	9.71	11.37	10.17	
Butyric acid	5.3	3	10.56	11.91	9.91	
Isovaleric acid	5.1	4	10.17	11.71	10.10	
Valeric acid	6.9	5	10.82	12.16	10.90	
4-Methyl pentanoic acid	8.5	6	10.61	11.49	9.48	
Hexanoic acid	9.2	7	13.33	13.20	9.81	

It is recommended to use an internal standard such as deuterated acetic acid if lower %RSDs are required. Quantitation results will be biased high at a low concentration region of the calibration curve and may require an internal standard to be injected at a low level. If the concentration of the internal standard is too high, it will not properly reflect the target compounds' behavior at low concentrations.

■ Conclusion

Analysis of aqueous short chain fatty acids is conducted with GCMS-QP2020 NX.

No derivatization or internal standard was used in this experiment and the measured concentration range was $4 \,\mu M$ to 100 µM. Concentration by steam distillation or dilution with water can be possibly explored to bring a target concentration into this measurement range(1, 2).



Fig. 4 Steam distillation

- (1) J.B. Zijlstra, J.Beukema, B.G.Wolthers, B.M.Byrne, A.Groen and J.Dankert, Pretreatment methods prior to gas chromatographic analysis of volatile fatty acids from faecal samples, Volume 78, Issue 2, January 14th 1977, Pages 243-250, ISSN 0009-8981
- (2) D.C.Dyer, A new method of steam distillation for the determination of the volatile fatty acids, including a series of colorimetric qualitative reactions for their identification, J.Biol. Chem. 1917, 28: 445-473

GCMS-QP and GCMS-TQ are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.



For Research Use Only. Not for use in diagnostic procedures. This publication may contain references to products that are not available in your country.

01-00027A-EN First Edition: Mar. 2021 Revision A: Mar. 2023

Please contact us to check the availability of these products in your country

Shimadzu Corporation www.shimadzu.com/an/

are used with trademark symbol "TM" or "@". Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change . without notice.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See http://www.shimadzu.com/about/trademarks/index.html for details.