

Targeted Screening of Coumarins and Furanocoumarins in Essential Oils Utilizing High-Resolution Quadrupole Time-of-Flight Mass Spectrometry

ASMS 2019 MP329

Jennifer Davis¹, Evelyn H Wang¹, Katie Pryor¹,
Jerry Byrne II¹, Helen Hao¹, Christopher Gilles¹
¹ Shimadzu Scientific Instruments, Columbia, MD

Targeted Screening of Coumarins and Furanocoumarins in Essential Oils Utilizing High-Resolution Quadrupole Time-of-Flight Mass Spectrometry

Overview

A 11 min method was developed using high resolution mass spectrometry for targeted screening of coumarins and coumarin derivatives within essential oils.

Introduction

Coumarins and Furanocoumarins are a compound class commonly found in fruit essential oils, such as citrus oil or bergamot oil. Previous research has found these compounds can be linked to phototoxic, mutagenic, or carcinogenic effects. Since this class of compounds is commonly found in essential oils, cosmetics, or bronzing products there is an emerging need to detect their

presence to prevent unnecessary harm to the consumer. Currently GC-MS has shown a limited ability to analyze furanocoumarins due to their relatively polar or heat-labile substituents. During this study a LCMS QToF was used for the screening and identification of coumarins and furanocoumarins in essential oils.




Figure 1 Structural backbone of coumarin and furanocoumarin

Methods

Essential oil samples were prepared by diluting 20 μL of oil into 980 μL of methanol. Ten coumarin standards were obtained in solid form. Each standard was dissolved in methanol and diluted to a final concentration of 500ppb.

Electrospray ionization was used in positive MS scan mode for the initial screening of the prepared samples. A total of thirty-five coumarin and furanocoumarins were

screened for in each essential oil. Primary confirmation was completed based on accurate mass. Secondary confirmation was completed by comparing MSMS spectral data and retention times of the prepared essential oil sample with a neat standard. In the situation where a neat standard was not available secondary confirmation was completed by comparison with metlin's MSMS spectral databases.

Targeted Screening of Coumarins and Furanocoumarins in Essential Oils Utilizing High-Resolution Quadrupole Time-of-Flight Mass Spectrometry

UHPLC conditions (Nexera system)	
Column	: Restek Raptor Biphenyl 100mm×2.1mm, 2.7 μm
Mobile phase A	: Water with 5mM Ammonium formate
B	: Methanol with 5mM Ammonium formate
Flow rate	: 0.4 mL/min
Time program	: B conc. 5%(0-0.5min) - 50%(1min) - 99%(8-9min) – 5%(9.01-11min)
Injection vol.	: 1 uL
Column temperature	: 40°C

MS conditions (LCMS-9030)	
Ionization	: ESI, Positive scan and MSMS mode

Results

A total of 9 essential oils were screened for coumarins and furanocoumarins, representative data is shown for bergamot oil, orange oil, and lemon oil. Utilizing an 11 min gradient method targeted screening of thirty-five coumarins or furanocoumarins was completed for all essential oils. Primary screening of each compound was confirmed within a 4 ppm accuracy of the theoretical mass. The largest presence of coumarins and furanocoumarins was detected in bergamot oil with a total of 22 (Figure 2). Lemon oil and Orange oil were confirmed to have 17 and 10 coumarins, respectively.

After primary screening, secondary confirmation was completed using retention time and fragmentation pattern comparison. Figures 3-5 show examples of each type of secondary confirmation. Ten coumarins were confirmed by retention time and fragmentation pattern comparison with neat standards acquired on the same high resolution mass spectrometer, LCMS-9030. Twenty-five coumarins were confirmed by fragmentation pattern comparison with known databases, such as Metlin or the human metabolome database.

Targeted Screening of Coumarins and Furanocoumarins in Essential Oils Utilizing High-Resolution Quadrupole Time-of-Flight Mass Spectrometry

Primary Confirmation of Coumarins

Table 1. Primary Screening of Essential oils and the ppm difference of the observed mass.

Coumarins	Bergamot Oil		Orange Oil		Lemon Oil	
	Observed Mass	ppm difference	Observed Mass	ppm difference	Observed Mass	ppm difference
4-methyldaphnetin	ND		ND		ND	
Esculetin	179.0340	0.61	ND		179.0340	0.61
6-Hydroxycoumarin	ND		ND		ND	
Isoscooletin	ND		ND		193.0497	0.83
6,7-Dihydroxy-4-methylcoumarin	193.0497	0.83	ND		193.0501	2.90
Daphnetin-7-methylether	193.0497	0.83	ND		ND	
Umbelliferone	ND		ND		ND	
Scopoletin	ND		ND		ND	
5,7-Dihydroxy-4-methylcoumarin	193.0497	0.83	ND		193.0498	1.35
8-Acetyl-6-hydroxy-7-methoxycoumarin	ND		ND		ND	
Fraxidin	223.0604	1.34	223.0603	0.90	223.0604	1.34
Xanthotoxol	203.0346	3.50	203.0339	0.05	ND	
6,7-Dimethylesculetin	ND		207.0654	0.77	207.0654	1.01
Coumarin	147.0443	1.63	ND		147.0443	1.63
8-Acetyl-7-methoxycoumarin	219.0651	-0.41	219.0654	0.96	ND	
Herniarin	177.0550	2.15	177.0546	-0.11	177.0551	2.71
4-methoxycoumarin	177.0552	3.28	ND		177.0546	-0.11
8-acetyl-6,7-dimethoxycoumarin	ND		ND		ND	
3-acetylcoumarin	189.0547	0.42	189.0546	-0.21	189.0547	0.42
7-methylcoumarin	ND		161.0596	-0.81	ND	
Psoralen	187.0390	0.16	ND		187.0390	0.16
Nordalbergin	ND		ND		255.0655	1.22
6-Methoxy-4-methylcoumarin	191.0704	0.68	ND		ND	
7-Methoxy-4-methylcoumarin	191.0701	-0.89	ND		ND	
Xanthotoxin	217.0497	0.74	ND		217.0498	1.20
6-methylcoumarin	ND		161.0597	-0.06	ND	
Dalbergin	269.0813	1.71	ND		ND	
Citropten	207.0660	3.91	ND		207.0655	1.50
Bergapten	217.0501	2.58	217.0496	0.28	217.0496	0.28
Isopimpinellin	247.0607	2.43	ND		247.0606	2.02
7-Ethoxycoumarin	ND		ND		ND	
4-Hydroxycoumarin	163.0392	1.41	163.0390	0.18	163.0391	0.80
4-Ethoxycoumarin	ND		ND		ND	
4-methylumbelliferone	177.0549	1.58	ND		ND	
4-methyl-7-ethoxycoumarin	205.0860	0.39	205.0857	-1.07	ND	

Targeted Screening of Coumarins and Furanocoumarins in Essential Oils Utilizing High-Resolution Quadrupole Time-of-Flight Mass Spectrometry




Figure 2. Representative Extracted Ion Chromatogram for all 22 Coumarins present in Bergamot oil.

Secondary Confirmation by Neat Standard




Figure 3. Overlaid chromatogram for confirmation of Citropten in Bergamot essential oil.
Pink- Citropten neat standard, Black- XIC of 207.0659m/z from Bergamot oil

Targeted Screening of Coumarins and Furanocoumarins in Essential Oils Utilizing High-Resolution Quadrupole Time-of-Flight Mass Spectrometry





Figure 4. MSMS spectral comparison for confirmation of Citropten in Bergamot essential oil.
Top- Citropten neat standard, Bottom- Bergamot oil sample


Targeted Screening of Coumarins and Furanocoumarins in Essential Oils Utilizing High-Resolution Quadrupole Time-of-Flight Mass Spectrometry

Secondary Confirmation by Metlin Spectral Database

Herniarin ESI Q-ToF MSMS spectrum from Metlin, CE= 20V



Herniarin ESI Q-ToF MSMS spectrum from Metlin, CE= 40V



Herniarin ESI Q-ToF MSMS spectrum from Bergamot Oil, CE= 15 to 55V




Figure 5. MSMS spectrum for Herniarin coumarin in Bergamot oil. Top- Herniarin by Metlin database (CE-20V), Middle Herniarin by Metlin database (CE-40V), Bottom- Bergamot oil sample

Targeted Screening of Coumarins and Furanocoumarins in Essential Oils Utilizing High-Resolution Quadrupole Time-of-Flight Mass Spectrometry

Conclusions

Targeted Screening of Coumarins and Furanocoumarins was successfully completed using high resolution mass spectrometry, LCMS-9030. Initial confirmation was completed for Orange, Lemon, Bergamot, Eucalyptus, Lavender, Tangerine, Chamomile, Sandalwood, and Sweet orange oil. Each essential oil contained between 6 to 22 coumarins all confirmed within 4ppm of their respective

accurate masses. Secondary confirmation was completed by either comparing MSMS spectral data from neat standards or by comparing MSMS spectral data from Metlin's database. Further research could be done to quantitate each coumarin if calibration curves were developed using neat standards for all 35 coumarins in the screening panel.

Disclaimer: The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures.

First Edition: July, 2019