


Solutions for

Analysis of Ethylene Oxide and 2-Chloroethanol in Food



Chemical Structure of EtO and 2-CE

Ethylene Oxide (EtO) (Molecular Formula – C_2H_4O)



2-Chloroethanol (2-CE) (Molecular Formula – C₂H₅ClO)

Significance of quantitation of EtO and 2-CE at trace levels in food and other matrices

Even though Ethylene oxide (EtO / EO) is well-known to be a toxic compound with carcinogenic and mutagenic concerns, it has been widely used for fumigation in the food industry because it effectively reduces or eliminates microbiological contamination with bacteria / fungi. Once in contact with food, EtO undergoes various reactions within the matrix, further producing reaction products such as ethylene glycol, 2-Chloroethanol (2-CE) and 2-bromoethanol, which are also toxic in nature. Hence, the use of EtO for food fumigation has been phased out in many countries due to toxicological concerns. In the EU, the use of EtO for the disinfection of foodstuffs, e.g. in storage areas, is not permitted (ECHA, 2020). The EU has also proposed separate maximum residual limits (MRLs) for EtO and its primary metabolite 2-CE in different food and agriculture commodities, ranging from 0.02 to 0.1 mg/kg (Commission Regulation (EU) 2015/868).

The recent recall of food products exported to the EU due to non-compliance with EU regulations has highlighted the importance of quantitation of EtO and 2-CE residues in food.

To ensure the quality and safety of food products, the European Rapid Alert System for Food and Feed (RASFF) prohibits the sales of goods exceeding the MRL values of 0.05 mg/kg (or 50 ppb) for the sum of EtO and 2-CE.

Reference: EURL-SRM - Analytical Observations Report

EU MRLs of EtO and 2-CE for different products are

No.	Product	EU MRLs
1	Teas, cocoa and spices	0.10 mg/kg
2	Nuts, oil fruits and oilseeds	0.05 mg/kg
3	Fruits, vegetables, sugar plants, fungi and pulses	0.02 mg/kg
4	Cereals and products of animal origin	0.02 mg/kg
5	Apicultural products	0.05 mg/kg

Different methods as per EURL-SRM report

Method Reference Components		Internal Standard	Liqu	Mode		
EURL-SRM Method	EtO	2-CE	Optional	QuEChERS-Method (EN 15662)	QuOil-Method (CEN/TS 17062:2019)	Liquid
Tadeo/Bononi Method -		2-CE	-	EtO to 2-CE- Conversion Method		Liquid
- :		2-BE*	-	EtO to 2-BE- Conversion Method		Liquid
Woodrow Method	EtO	-	-	-	-	Headspace
Ayoub Method	EtO	-	-	-	-	SPME

^{*2-}BE: 2-Bromoethanol

Shimadzu Solutions

Shimadzu's solutions for analysis of EtO and 2-CE (Different Techniques / Instruments)

No.	Matrix Covered	Technique	Instrument	Page
1.1		Liquid Injection (EtO and 2-CE)	GCMS-TQ8050 NX with AOC-20i/s	4
1.1		Dynamic Headspace (EtO and 2-CE) - Solvent extraction	GCMS-TQ8050 NX with HS-20 NX	4
1.2	Sesame Seeds	Dynamic Headspace (EtO and 2-CE) - Direct sample in HS vial	GCMS-TQ8040 NX with HS-20 NX	5
1.3	1.3	Liquid Injection (2-CE) -Conversion Method	GCMS-TQ8040 NX with AOC-20i/s	7
		Dynamic Headspace Injection (2-CE) - Conversion Method	GCMS-TQ8040 NX with HS-20 NX	,
2.1		Dynamic Headspace (EtO and 2-CE) - Solvent extraction	GCMS-TQ8040 NX with AOC-30i	9
2.1	1 Chili Powder	Dynamic Headspace (EtO and 2-CE) - Solvent extraction	GCMS-TQ8050 NX with AOC-20i/s	_
2.2		Dynamic Headspace (EtO and 2-CE) - Direct sample in HS vial	GCMS-TQ8050 NX with HS-20 NX	12
3.1	Noodles	Liquid Injection (EtO and 2-CE)	GCMS-TQ8050 NX with AOC-20i/s	15

1.1 Sesame Seeds Analysis



Analysis of EtO and 2-CE in Sesame Seeds

♦ This section deals with various approaches, extraction process and instrument techniques for the analysis of EtO and 2-CE in sesame seeds.

APPROACH	MODE OF ANALYSIS	INJECTOR	INSTRUMENT TECHNIQUE
EURL-SRM Analytical Observation Report (EtO and 2-CE)	Liquid	SPL	AOC-20i/s - GCMS-TQ8050 NX
Headspace - Method 1st (EtO and 2-CE)	Dynamic Headspace	SPL	HS-20 NX - GCMS-TQ8050 NX
Headspace - Method 2 nd (2-CE)	Dynamic Headspace	SPL	HS-20 NX - GCMS-TQ8050 NX
Headspace - Method 3 rd (EtO)	Dynamic Headspace	SPL	HS-20 NX - GCMS-TQ8050 NX



Trace level quantitation of Ethylene Oxide (EtO) and 2-Chloroethanol (2-CE) in sesame seeds by using various GC-MS/MS techniques with their own merits and demerits

► Headspace – Method 1st (EtO and 2 CE): Solvent extraction followed by analysis with the dynamic headspace mode. For simultaneous analysis of EtO and 2 CE in different matrices, this method is used most of the time.

The European Chemical Agency (ECHA) has classified EtO in category 1B, related to carcinogenicity, mutagenicity and reproductive toxicity, and in category 3, related to acute toxicity. The US National Institutes of Health (NIH) classified EtO as "known to be a human carcinogen based on sufficient evidence of carcinogenicity from studies in humans, including epidemiological studies and studies on mechanisms of carcinogens." The US Environmental Protection Agency (EPA) has concluded that EtO is carcinogenic to humans by the inhalation route of exposure. 2-CE is a prominent metabolite of EtO and is an equally hazardous compound. EU MRLs (Maximum Residue Levels as per European Commission) for EtO and 2-CE are different for different commodities. Out of many matrices, the EU MRLs for cereals, pulses and vegetables are the lowest at 0.02 mg/kg. Considering carcinogenicity and no acceptable threshold for exposure, no Acceptable Daily Intake (ADI) was established for EtO and 2-CE. Hence, it is very important to quantitate EtO and 2-CE as low as possible in food matrices.

With reference to the EURL-SRM Analytical Observation Report, Shimadzu has successfully developed and validated methods for trace-level quantification of EtO and 2-CE impurities in sesame seeds using the GCMS-TQ8050 NX with AOC-20i auto injector and AOC-20s auto sampler, and with HS-20 NX dynamic headspace sampler.





Analysis of EtO and 2-CE in Sesame Seeds (Direct Sample In Headspace Vial)

- In this approach, a sesame seed sample was placed directly into a headspace vial and analyzed for its EtO and 2-CE content.
- This method has been successfully validated using the GCMS-TQ8040 NX with HS-20 NX headspace sampler.

APPROACH	MODE OF ANALYSIS	INJECTOR	INSTRUMENT TECHNIQUE
Direct Placement of Sample in Headspace Vial (EtO and 2-CE)	Dynamic Headspace	SPL	HS-20 NX - GCMS-TQ8040 NX

System Description

GC-MS System	GCMS-TQ8040 NX
Headspace Sampler	HS-20 NX (Dynamic Headspace)
Capillary Column	SH-RTX-VMS (60 m X 0.25 mm ID x 1.4 um df)
Ionization Mode	Electron Ionization (EI)
Sample Information	EtO and 2-CE in sesame seeds
Sample Preparation	100 % (Sample in HS vial) + 1 uL Acetonitrile (Diluent)
Spiked Sample Preparation	100 mg Sesame Seeds + 1 uL of 1 ppm Standard
	100 mg Sesame Seeds + 1 uL of 2 ppm Standard
	100 mg Sesame Seeds + 1 uL of 3 ppm Standard

Benefits

- Easier sample preparation with no solvent extraction or any conversion or derivatization process
- · Shorter sample preparation time
- · Very cost effective
- · Outperforms current regulatory limits



Triple Quadrupole Gas Chromatograph Mass Spectrometer

GCMS-TQ8040 NX


- Simplified maintenance allows operators to change parts quickly
- Smart Database Series enables automated method development
- Smart MRM sets method to maximize sensitivity



Chromatography Parameters								
GC-MS System	:	GCMS-TQ804	GCMS-TQ8040 NX					
Liquid Sampler	:	-						
Headspace Sampler	:	HS-20 NX (D	HS-20 NX (Dynamic Headspace)					
Gas Chromatography Parameters								
Capillary Column	:	SH-VMS (60 m X 0.25	mm ID x 1.4 un	n df)				
Injection Mode	:	Split						
Flow Control Mode	:	Column Flow	'					
Carrier Gas	:	Helium						
Column Flow	:	2.00 mL/min						
Linear Velocity	:	36.0 cm/s						
Purge Flow	:	3.0 mL/min						
Split Ratio	:	10						
Diluent	:	Acetonitrile						
Temp. Program		Ramp Rate (°C/min)	Temp. (°C)	Hold Time (min)				
			35.0	5.00				
		30	235.0	0.33				
MS Parameters								
Ionization Mode	:	EI						
Ion Source Temp.	:	250 °C						
Interface Temp.	:	: 240 °C						
CID Gas	:	Argon						
CID Gas pressure	D Gas pressure : 200 kPa							

Headspace Parameters		
Method	:	Dynamic HS
Oven Temp.	:	110 °C
Sample Line Temp.	:	120 °C
Transfer Line Temp.	:	130 °C
Trap Cooling Temp.	:	-20 °C
Trap Desorb Temp.	:	280 °C
Trap Equilib. Temp.	:	-20 °C
Shaking Level	:	5
Multi Inj. Count	:	1
Pressurizing Gas Pressure (kPa)	:	192
Equilibrating Time (min)	:	15
Pressurizing Time (min)	:	0.5
Pressure Equilib. Time (min)	:	0.0
Load Time (min)	:	0.5
Load Equilib. Time (min)	:	0.0
Dry Purge Time (min)	:	0.0
Injection Time (min)	:	10
Needle Flush Time (min)	:	10
GC Cycle Time (min)	:	30
Split Ratio	;	10
Total Flow (mL)	;	25
Trap Tube	:	Tenax Ta

Direct Sample Method - Analysis of EtO and 2-CE (Dynamic Headspace) Summary of results is shown in the table below.



2-CE

Figure 2: Calibration curve, overlay of linearity standards and chromatogram of LOQ solution

Results

Method =>	Dynar	nic Headspace
ivietnod =>	EtO	2-CE
Linearity levels (On column)	10,20 and 30 ppb	10,20 and 30 ppb
Linearity levels (w.r.t sample)	10,20 and 30 ppb	10,20 and 30 ppb
r² (n=3)	0.99696	0.99503
Level 1- conc. (On column)	10 ppb	10 ppb
Level 1- conc. (w.r.t sample)	10 ppb	10 ppb
% RSD (n=3)	2.2	10.1
S/N	26	28
Level 2- conc. (On column)	20 ppb	20 ppb
Level 2- conc. (w.r.t sample)	20 ppb	20 ppb
% RSD (n=3)	2.7	11.7
Level 3- conc. (On column)	30 ppb	30 ppb
Level 3- conc. (w.r.t sample)	30 ppb	30 ppb
% RSD (n=3)	9.7	14.1
Pre-spike 10 ppb-1st (Recovery)	97 %	101 %
Pre-spike 10 ppb-2 nd (Recovery)	105 %	91 %





Analysis of EtO and 2-CE in Sesame Seeds (Conversion)

With reference to the EURL-SRM Analytical Observation Report, the conversion method (EtO=>2-CE) for the trace-level quantification of 2-CE in sesame seeds has been successfully optimized and validated using the GCMS-TQ8040 NX with AOC-20i auto injector and AOC-20s auto sampler, and with HS-20 NX dynamic headspace sampler.

APPROACH	MODE OF ANALYSIS	INJECTOR	INSTRUMENT TECHNIQUE	
EtO to 2-CE Conversion	Liquid	SPL	AOC-20i/s - GCMS-TQ8040 NX	
EtO to 2-CE Conversion	Dynamic Headspace	SPL	HS-20 NX - GCMS-TQ8040 NX	

System Description

GC-MS System	GCMS-TQ8040 NX
Liquid Injector and Sampler	AOC-20i and AOC-20s*
Headspace Sampler	HS-20 NX (Dynamic Headspace)
Capillary Column	SH-Stabilwax DA (60 m X 0.25 mm ID x 0.25 um df)
Ionization Mode	Electron Ionization (EI)
Sample Information	2-CE in sesame seeds (Conversion method)
Sample Preparation	20 % for Liquid / Dynamic Headspace
Reference Method	EURL-SRM Analytical Observation report


^{*} Equivalent data can be obtained using AOC-30i with or without AOC-20s U.

Benefits

- With reference to EURL-SRM report, the conversion method for trace-level quantification of 2-CE has been optimized
 - A) EtO to 2-CE conversion Analyzed by liquid injection
 - B) EtO to 2-CE conversion Analyzed by dynamic headspace injection
- GCMS-TQ8040 NX equipment with AOC-20i auto injector and AOC-20s auto sampler, and with HS-20 NX dynamic headspace sampler is well-suited for the analysis of EtO and 2-CE


Sample Analysis

Conversion of EtO to 2-CE in sesame seeds and analyzed by liquid injection



Chromatography Parameters					Headspace Parameters		
GC-MS System		GCMS-TQ804	10 NX		Method	Ι.	Dynamic HS
Liquid Sampler		AOC-20i and AOC-20s			Oven Temp.		110 °C
· · ·							
Headspace Sampler : HS-20 NX (Dynamic Headspace)			Sample Line Temp.	:	170 °C		
Gas Chromatography Parameters			Transfer Line Temp.	:	180 °C		
Capillary Column	:	SH-Stabilwax (60 m X 0.25	: DA mm ID x 0.25 u	m df)	Trap Cooling Temp.	:	-10 °C
Injection Mode	:	Split			Trap Desorb Temp.	:	280 °C
Flow Control Mode	:	Column Flow	,		Trap Equilib. Temp.	:	25°C
Carrier Gas	:	Helium			Shaking Level	:	5
Column Flow	:	1.00 mL/min			Multi Inj. Count		10
Linear Velocity	:	25.8 cm/s			Pressurizing Gas Pressure (kPa)	:	192
Purge Flow	:	3.0 mL/min			Equilibrating Time (min)	:	15
Split Ratio	:	1 (For liquid i	njection metho	d)	Pressurizing Time (min)	:	0.5
Injector Port Temp	:	230.0 °C			Pressure Equilib. Time (min)	:	0.1
Temp. Program		Ramp Rate (°C/min)	Temp. (°C)	Hold Time (min)	Load Time (min)	:	0.5
			60.0	2.00	Load Equilib. Time (min)	:	0.1
		20	235.0	9.25	Dry Purge Time (min)	:	1
MS Parameters					Injection Time (min)	:	5
Ionization Mode	:	EI			Needle Flush Time (min)	:	15
lon Source Temp.	: 230 °C			GC Cycle Time (min)	:	30	
Interface Temp. : 230 °C			Split Ratio	:	15		
CID Gas : Argon			Total Flow (mL)	:	19		
CID Gas pressure :		200 kPa			Trap Tube	:	Tenax Ta

Conversion Method - EtO =>2-CE (Liquid and Dynamic)






Figure 2: Calibration curve, overlay of linearity standards and chromatogram of LOQ solution for 2-CE (Dynamic)

Results

Summary of results is shown in the table below.

1		
Method =>	Liquid	Dynamic HS
	2-CE	2-CE
Linearity levels (On column)	2,5,10,15 and 20 ppb	2,5,10,15 and 20 ppb
Linearity levels (w.r.t sample)	10,25,50,75 and 100 ppb	10,25,50,75 and 100 ppb
r² (n=3)	0.99919	0.99835
LOQ level conc. (On column)	2 ppb	2 ppb
LOQ level conc. (w.r.t sample)	10 ppb	10 ppb
% RSD (n=6)	7.6	3.8
S/N	16	44
Spiked LOQ level (On column)	2 ppb	2 ppb
Spiked LOQ level (w.r.t sample)	10 ppb	10 ppb
Avg of % recovery (n=3)	101 %	95 %

2.1 Chili Powder Analysis



Analysis of EtO and 2-CE in Chili Powder

 This section deals with various approaches, extraction process and instrument techniques for the analysis of EtO and 2-CE in a chili powder sample.

APPROACH	MODE OF ANALYSIS	INJECTOR	INSTRUMENT TECHNIQUE
EURL-SRM Analytical Observation Report (EtO and 2-CE)	Liquid	SPL	AOC-20i/s - GCMS-TQ8040 NX (AOC-30i - GCMS-TQ8050 NX)*
Headspace - Method 1st (EtO and 2-CE)	Dynamic Headspace	SPL	HS-20 NX - GCMS-TQ8050 NX

^{*}Analysis by GCMS-TQ8050 NX is to show comparison only

System Description

•	
GC-MS System	GCMS-TQ8040 NX
Liquid Injector	AOC-30i
GC-MS System	GCMS-TQ8050 NX
Liquid Injector and Sampler	AOC-20i and AOC-20s
Headspace Sampler	HS-20 NX (Dynamic Headspace)
Capillary Column	RTX-VMS (60 m X 0.45 mm ID x 2.55 um df)
Ionization Mode	Electron Ionization (EI)
Sample Information	EtO and 2-CE in chili powder
Reference Method	EURL-SRM Analytical Observation report

 $^{{\}color{blue} * Chromatographic conditions and sample preparation method are the same as listed in the application note.}\\$

Benefits

- The liquid injection method and developed dynamic headspace method for the analysis of EtO and 2-CE in chili powder matrix have been optimized. Both methods outperform the current regulatory limits
- GCMS-TQ8040 NX / GCMS-TQ8050 NX with auto injector and auto sampler, and with dynamic headspace sampler is a complete tool for the analysis of EtO and 2-CE



Triple Quadrupole Gas Chromatograph Mass Spectrometer

GCMS-TQ8050 NX

- Highly efficient detector and three forms of noise-reduction technologies enable ultra-high sensitivity detection
- Enhanced sensitivity enables ultra-trace level analysis
- Simplified user maintenance allows operators to change parts quickly
- Reduced maintenance frequency and lower long-term operational costs



Solvent Extraction Method for EtO and 2-CE (Liquid)




Figure 1: Calibration curve, chromatogram of LOQ solution and overlay of linearity standards for EtO




Figure 2: Calibration curve, chromatogram of LOQ solution and overlay of linearity standards for 2-CE

Results

Linearity

Summary of calibration standard is shown in Table 1.

Table 1 : Summary for linearity (n=3 for each level)

Maderal	Liquid		
Method =>	EtO	2-CE	
Linearity levels (On column)	5,10,15,20 & 25 ppb		
Linearity levels (w.r.t sample)	10,20,30,40 & 50 ppb		
r ² (n=3)	0.99879 0.99965		

Precision

Summary of precision standard is shown in Table 2.

Table 2: Summary of precision standard (n=6)

Madada	Liquid		
Method =>	EtO	2-CE	
LOQ level conc. (On column)	5 ppb	5 ppb	
LOQ level conc. (w.r.t sample)	10 ppb	10 ppb	
% RSD (n=6)	9.3	3.9	
S/N	15	79	
Highest level conc. (On column)	25 ppb	25 ppb	
Highest level conc. (w.r.t sample)	50 ppb	50 ppb	
% RSD (n=6)	5.6	3.9	

Accuracy

Summary of accuracy is shown in Table 3.

Table 3 : Summary for accuracy (n=3 for each level)

Made at	Liquid		
Method =>	EtO	2-CE	
Spiked LOQ level conc. (On column)	5 ppb	5 ppb	
Spiked LOQ level conc. (w.r.t sample)	10 ppb	10 ppb	
Avg of % recovery (n=3)	98 %	103 %	
Spiked 2 nd level conc. (On column)	10 ppb	10 ppb	
Spiked 2 nd level conc. (w.r.t sample)	20 ppb	20 ppb	
Avg of % recovery (n=3)	96 %	118 %	
Spiked 3 rd level conc. (On column)	15 ppb	15 ppb	
Spiked 3 rd level conc. (w.r.t sample)	30 ppb	30 ppb	
Avg of % recovery (n=3)	99 %	107 %	

Comparison

Summary of comparison for two instruments is shown in Table 4.

Table 4: Summary for comparison

Method => Liquid	GCMSTQ8040	GCMS-TQ8050
iviethod => Liquid	EtO	EtO
LOQ level conc. (On column)	5 ppb	2.5 ppb
LOQ level conc. (w.r.t sample)	10 ppb	5 ppb
% RSD (n=6)	9.3	8.7
S/N	15	32
	2-CE	2-CE
LOQ level conc. (On column)	5 ppb	2.5 ppb
LOQ level conc. (w.r.t sample)	10 ppb	5 ppb
% RSD (n=6)	3.9	5.0
S/N	79	70

Solvent Extraction Method for EtO and 2-CE (Dynamic)




Figure 3: Calibration curve, chromatogram of LOQ solution and overlay of linearity standard for EtO




Figure 4: Calibration curve, chromatogram of LOQ solution and overlay of linearity standard for 2-CE

Results

Linearity

Summary of calibration standard is shown in Table 5.

Table 5 : Summary for linearity (n=3 for each level)

Market .	Liquid		
Method =>	EtO	2-CE	
Linearity levels (On column)	10,20,30,40 & 50 ppb		
Linearity levels (w.r.t sample)	10,20,30,40 & 50 ppb		
r ² (n=3)	0.99877	0.99872	

Precision

Summary of precision standard solutions is shown in Table $\, 6. \,$

Table 6 : Summary for precision (n=6)

Mashad	Liquid		
Method =>	EtO	2-CE	
LOQ level conc. (On column)	10 ppb	10 ppb	
LOQ level conc. (w.r.t sample)	10 ppb	10 ppb	
% RSD (n=6)	2.5	9.7	
S/N	15	29	

Accuracy

Summary of accuracy is shown in Table 7.

Table 7: Summary for accuracy (n=3 for each level)

M. d l	Dynamic Headspace		
Method =>	EtO	2-CE	
Spiked LOQ level conc. (On column)	10 ppb	10 ppb	
Spiked LOQ level conc. (w.r.t sample)	10 ppb	10 ppb	
Avg of % recovery (n=3)	90%	93%	
Spiked 2 nd level conc. (On column)	30 ppb	30 ppb	
Spiked 2 nd level conc. (w.r.t sample)	30 ppb	30 ppb	
Avg of % recovery (n=3)	87%	114%	
Spiked 3 rd level conc. (On column)	50 ppb	50 ppb	
Spiked 3 rd level conc. (w.r.t sample)	50 ppb	50 ppb	
Avg of % recovery (n=3)	71%	105%	



Auto Injector and Auto Sampler for GC and GC-MS

AOC-30i and AOC-20s U

- Achieve consistent results with high reproducibility
- Experience unwavering performance through long-term stability
- Compact design reduces space requirements







Analysis of EtO and 2-CE in Chili Powder (Direct Sample In Headspace Vial)

- ♦ In this approach, a chili powder sample was placed directly into a headspace vial and analyzed for its EtO and 2-CE content.
- This method has been successfully validated using the GCMS-TQ8040 NX with HS-20 NX headspace sampler.

APPROACH	MODE OF ANALYSIS	INJECTOR	INSTRUMENT TECHNIQUE
Place sample in headspace vial (EtO and 2-CE)	Dynamic Headspace	SPL	HS-20 NX - GCMS-TQ8040 NX

System Description

GCMS-TQ8040 NX
HS-20 NX (Dynamic Headspace)
SH-RTX-VMS (60 m X 0.25 mm ID x 1.4 um df)
Electron Ionization (EI)
EtO and 2-CE in chili powder
100 % (As such sample in HS vial) + 1 uL of Acetonitrile (Diluent)
100 mg chili powder + 1 uL of 1 ppm standard
100 mg chili powder + 1 uL of 2 ppm standard
100 mg chili powder + 1 uL of 3 ppm standard
100 mg chili powder + 1 uL of 4 ppm standard
100 mg chili powder + 1 uL of 5 ppm standard

Benefits

- Easier sample preparation with no solvent extraction or any conversion or derivatization process
- · Shorter sample preparation time
- · Highly cost effective
- · Outperforms current regulatory limits

Chromatography Parameters							
GC-MS System	:	GCMS-TQ8040 NX					
Liquid Sampler	:	-					
Headspace Sampler	:	HS-20 NX (Dynamic Headspace)					
Gas Chromatography Parameters							
Capillary Column	:	SH-VMS (60 m X 0.25 mm ID x 1.4 um df)					
Injection Mode	:	Split					
Flow Control Mode	:	Column Flow					
Carrier Gas	:	Helium					
Column Flow	:	2.00 mL/min					
Linear Velocity	:	36.0 cm/s					
Purge Flow	:	3.0 mL/min					
Split Ratio	:	8					
Diluent	:	Acetonitrile					
Temp. Program		Ramp Rate (°C/min) Temp. (°C) Hold Tir					
			35.0	5.00			
		40	235.0	5.00			
MS Parameters	MS Parameters						
Ionization Mode	:	EI					
Ion Source Temp.	:	250 °C					
Interface Temp.	:	240 °C					
CID Gas	:	Argon					
CID Gas pressure	:	200 kPa					

Headspace Parameters		
Method	:	Dynamic HS
Oven Temp.		70 °C
Sample Line Temp.		110 °C
Transfer Line Temp.		120 °C
Trap Cooling Temp.		-20 °C
Trap Desorb Temp.		280 °C
Trap Equilib. Temp.		-20 °C
Shaking Level	:	5
Multi Inj. Count		3
Pressurizing Gas Pressure (kPa)		192
Equilibrating Time (min)		15
Pressurizing Time (min)		0.5
Pressure Equilib. Time (min)		0.1
Load Time (min)		0.5
Load Equilib. Time (min)		0.1
Dry Purge Time (min)		0.0
Injection Time (min)		10
Needle Flush Time (min)		10
GC Cycle Time (min)		28
Split Ratio		8
Total Flow (mL)		21
Trap Tube		Tenax Ta




Gas Chromatograph Mass Spectrometer with Headspace Sampler

GCMS-TQ8040 NX with HS-20 NX Series


- Optimal solution for volatile component analysis
- Upgraded HS-20 NX offers a powerful solution
- User-friendly option for both research and quality control



Direct Sample Method - Analysis of EtO and 2-CE (Dynamic Headspace)



 $\label{thm:continuous} \textbf{Figure 1: Calibration curve, overlay of linearity standards and chromatogram of sample solution } \\$



 $Figure \ 2: Calibration \ curve, \ overlay \ of \ linearity \ standards \ and \ chromatogram \ of \ sample \ solution$

Results

Summary of results is shown in the table below.

	Dynamic Headspace		
Method =>	EtO	2-CE	
Linearity levels (On column)	10,20,30,40 and 50 ppb	10,20,30,40 and 50 ppb	
Linearity levels (w.r.t sample)	10,20,30,40 and 50 ppb	10,20,30,40 and 50 ppb	
r ² (n=3)	0.99908	0.99371	
Level 1 - conc. (On column)	10 ppb	10 ppb	
Level 1 - conc. (w.r.t sample)	10 ppb	10 ppb	
% RSD (n=3)	12.9	3.2	
S/N	42	97	
Sample (Conc. observed)	-	5.8 ppb	
Pre spiked 10 ppb (Conc. observed)	10.8 ppb	17.3 ppb	
Avg of % recovery (n=3)	108 %	109 %	

3.1 Noodles



Analysis of EtO and 2-CE in Ramen (Instant Noodles)

 This section introduces the extraction process and instrumental techniques used for the analysis of EtO and 2-CE in Ramen.

APPROACH	MODE OF ANALYSIS	INJECTOR	INSTRUMENT TECHNIQUE
EURL-SRM Analytical Observation Report (EtO and 2-CE)	Liquid	SPL	AOC-20i and AOC-20s* - GCMS-TQ8050 NX

^{*} Equivalent data can be obtained using AOC-30i with or without AOC-20s U.



Analysis of Ethylene Oxide in Ramen (Instant Noodles) by GC-MS/MS



- GCMS-TQ8050 NX equipment with AOC-20i and AOC-20s provides stable analytical results.
- Excellent reproducibility in actual sample analysis demonstrates an RSD below 5%.
- · Good recovery rates were achieved for EtO (98.9%) and 2-CE (106.3%) in a complex food matrix.

EtO is typically used as a fumigant pesticide to reduce microbial or bacterial contamination. However, it is banned in the European Union (EU) because of its carcinogenic and mutagenic properties. In 2020, Belgium was the first country to raise the alarm about the presence of EtO in imported sesame seed. Since then, EtO has been found in various food additives, including locust bean gum, which is a thickening agent or stabilizer. These ingredients are commonly used in formulations and can be found in products such as flour, cereals, ice cream, chocolate, biscuits, bread, and cheese. As these items are sold by major brands and retailers, it has led to thousands of products being affected and taken off the shelves.

With the surge in demand for EtO detection worldwide, there is an increased need for greater laboratory capacity and relevant analysis methods to ensure food safety. In this Technology Brief, Shimadzu introduces an analysis of EtO with QuEChERS and the GCMS-TQ8050 NX, which is equipped with the AOC-20i and AOC-20s for automation capabilities to meet the requirement for higher throughput and faster turnaround times.



Shimadzu Corporation www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.