

# Application News

## No. A517

### Spectrophotometric Analysis

## Analysis of Discoloration and Coloration Using FTIR and EDX

With just minor discoloration which is now a sufficient reason for complaint, it is no easy task for manufacturers to investigate root cause of discoloration from a wide range of possibilities and pinpoint the correct one. When the quantity of material causing discoloration or coloration is small, extraction is normally performed before analysis. In this article, using the merits of Fourier transform infrared spectrophotometer (FTIR) and energy-dispersive X-ray fluorescence spectrometer (EDX), we identified discoloration materials and coloring compounds by direct sample measurement without pretreatment. Focusing on an analysis of color, this article describes some example analyses that explore causes of discoloration and coloration from organic and inorganic compounds, using an analytical approach that combines FTIR and EDX.

### ■ Analysis of Discolored Paper Using FTIR

Using FTIR, normal paper and discolored paper were analyzed, the results were compared, and the cause of discoloration was explored. Fig. 1 shows images of the paper samples, and Table 1 shows the FTIR analytical conditions used.

(1) Samples

Normal paper and discolored paper



**Fig. 1 Normal Paper (left) and Discolored Paper (right)**

(2) Analysis


Approach: Single reflection ATR

**Table 1 FTIR Analytical Conditions**

Instruments	: IRTracer-100, MIRacle 10 (Diamond prism)
Resolution	: 4.0 cm <sup>-1</sup>
Accumulation	: 40
Apodization	: Happ-Genzel
Detector	: DLATGS

(3) Results

Analysis: Calculation of difference spectrum using data computing function



**Fig. 2 Infrared Spectra and Search Result**

Fig. 2 shows the infrared spectra and difference spectrum (discolored paper - normal paper) for normal paper and discolored paper, and a search result for the difference

spectrum. Although paraffin oil was found after performing a spectrum search, paraffin is normally colorless at room temperature. The discoloration is therefore assumed to be a paraffin oil (e.g., machine oil) containing a discoloration-causing additive that adhered to the paper.

The peak at around 1750 cm<sup>-1</sup> (C=O bond) in the difference spectrum identified by the arrow in Fig. 2 is presumed to be derived from the additive in the paraffin oil.

Because the discolored paper spectrum is dominated by the base paper material, it is difficult to ascertain differences with the normal paper spectrum. Calculating a difference spectrum shows these differences clearly, and simplifies identification of discoloration.

### ■ Analysis of Discolored Resin Product Using EDX

Using EDX, measurements were taken from a normal area and discolored area of a resin product, the results compared, and the cause of discoloration explored. Fig. 3 shows images of the samples, and Table 2 shows the EDX analytical conditions used.

(1) Samples

Resin product with discolored area



**Fig. 3 Resin Product (left), Normal Area (center) and Discolored Area (right)**

(2) Analysis

Approach: <sup>6</sup>C-<sup>92</sup>U qualitative/quantitative analysis


**Table 2 EDX Analytical Conditions**

Instrument	: EDX-8000
X-ray Tube	: Rh target
Voltage / Current	: 50 kV (Al-U) / Auto
Atmosphere	: Vacuum
Measurement Diameter	: 1 mm φ
Integration Time	: 100 sec

(3) Results

Analysis: Calculation of difference profile using blank-corrected function

Detected elements: <sup>13</sup>Al (99.8 %), <sup>26</sup>Fe (0.2 %)



**Fig. 4 Blank-Corrected Qualitative Profile**


Fig. 4 shows the blank-corrected qualitative profile obtained by subtracting the profile of the normal area of resin from the profile of the discolored area. <sup>13</sup>Al and <sup>26</sup>Fe were only detected from the discolored area, with quantitative results showing the discolored area is 99.8 % <sup>13</sup>Al and 0.2 % <sup>26</sup>Fe. Based on these results, the discoloration is assumed to be a fragment of aluminum alloy, and could be from a cutting edge that adhered to the product during the cutting process.

**Analysis of Coloration Using FTIR and EDX**

FTIR and EDX were used to take measurements from white resin and green resin, and identify the coloring compounds in each resin. Fig. 5 shows images of the samples, Table 3-1 shows FTIR analytical conditions, and Table 3-2 shows EDX analytical conditions.

(1) Samples

White resin and green resin



**Fig. 5 White Resin (left) and Green Resin (right)**

(2) Analysis

Approach (FTIR): Single reflection ATR

Approach (EDX): <sup>6</sup>C-<sup>92</sup>U qualitative/quantitative analysis


**table 3-1 FTIR Analytical Conditions**

Instruments	: IRTracer-100, MIRacle 10 (Diamond prism)
Resolution	: 4.0 cm <sup>-1</sup>
Accumulation	: 40
Apodization	: Happ-Genzel
Detector	: DLATGS

**Table 3-2 EDX Analytical Conditions**

Instrument	: EDX-8000
X-ray Tube	: Rh target
Voltage / Current	: 15 kV (C-Sc, S-K)
	: 50 kV (Ti-U, Zn-As, Pb) / Auto
Atmosphere	: Vacuum
Measurement Diameter	: 10 mm φ
Integration Time	: 60 sec

(3) Results of FTIR Qualitative Analysis



**Fig. 6 Infrared Spectra and Search Result**


Fig. 6 shows infrared spectra for white resin and green resin and a spectrum search result. The results showed that polyethylene is the resin component present in both the white resin and green resin. It is presumed that no coloring component was detected because only a very small quantity is present in each sample.

(4) Results of EDX Qualitative/Qualitative Analysis

**Table 4 Quantitative Results**

	Cl	Si	P	Al	S	Cu	C <sub>2</sub> H <sub>4</sub>
White	-	0.010	0.007	0.007	0.003	-	99.97
Green	0.030	0.012	0.007	0.006	0.002	0.002	99.94

Units: % -: Undetected  
Note: Polyethylene (C<sub>2</sub>H<sub>4</sub>) was set as the balance.<sup>1)</sup>




**Fig. 7 Qualitative Profiles**

Table 4 shows the quantitative results and Fig. 7 shows the qualitative profiles. <sup>13</sup>Al, <sup>14</sup>Si, <sup>15</sup>P, and <sup>16</sup>S were detected in both the white resin and green resin, and <sup>17</sup>Cl and <sup>29</sup>Cu were only detected in the green resin. FTIR analysis showed the resins are polyethylene, which is normally transparent, so detected elements are expected to contribute to resin coloration.

First, because <sup>14</sup>Si, <sup>15</sup>P, <sup>13</sup>Al, and <sup>16</sup>S were detected in both resin samples, these elements are presumed to be derived from additives used to improve resin functionality, and not from compounds used to impart color. The filler aluminum silicate (kaolin, Al<sub>2</sub>Si<sub>2</sub>O<sub>5</sub>(OH)<sub>4</sub>) is a likely additive that contains <sup>13</sup>Al and <sup>14</sup>Si. <sup>15</sup>P and <sup>16</sup>S are often found in oxidation inhibitors.

As for <sup>17</sup>Cl and <sup>29</sup>Cu only found in the green resin, green pigments that contain these elements include the organic pigment phthalocyanine green (C<sub>32</sub>Cl<sub>16</sub>CuN<sub>8</sub>), and the inorganic pigment atacamite (Cu<sub>2</sub>(OH)<sub>3</sub>Cl). Ferrous pigments, cadmium pigments, and chrome pigments are also used to impart green coloration, but these materials do not match the elements and composition detected in the sample. EDX qualitative results of 0.03 % <sup>17</sup>Cl and 0.002 % <sup>29</sup>Cu show a higher <sup>17</sup>Cl content than <sup>29</sup>Cu. Phthalocyanine green is probably the best match when the composition of candidate pigments and qualitative results obtained by EDX analysis are compared.

Based on the above findings, the white color in the white resin was assumed to be derived from an additive, and the green color in the green resin from phthalocyanine green. Fig. 8 shows the structure of phthalocyanine green for reference.



Reference: Chemical Risk Information Platform (CHRIP), National Institute of Technology and Evaluation

**Fig. 8 Phthalocyanine Green**

**Conclusion**

Causes of discoloration can be quickly identified with FTIR when discoloration is caused by an organic material, and with EDX when discoloration is caused by an inorganic material. When analyzing a colored material, FTIR was effective in determining its main constituents, and EDX was effective in identifying the pigment.

Combining FTIR and EDX provides a non-destructive and speedy analysis that is effective in the investigation of discoloration and coloration.

References

- 1) Shimadzu Application News No. X255



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.