

# Application Data Sheet

## No. 126

## System Gas Chromatograph

### BID Ultra-Fast Refinery Gas Analyzer Nexis GC-2030 BIDUFRGA

This GC system is for determining the chemical composition of natural gases and similar gaseous mixtures within the composition range shown below. This test method provides data for calculating physical properties of the sample, such as heating value and relative density, or for monitoring the concentrations of one or more of the components in a mixture. A total of 3 valves and 6 columns are used in this GC system. The sample is loaded into three sample loops for determination. Using a pre-column, C6+ components are back-flushed as a single peak. The valve timing then allows the hydrocarbons C1 through to C5 to be separated by an Alumina capillary column and detected by FID. A MS-5A separates H<sub>2</sub>, O<sub>2</sub>, N<sub>2</sub>, CH<sub>4</sub>, CO while CO<sub>2</sub>, C<sub>2</sub>H<sub>4</sub>, C<sub>2</sub>H<sub>6</sub>, C<sub>2</sub>H<sub>2</sub>, H<sub>2</sub>S are separated by Rtx-Q plot column and detected by a BID. The final analysis time is approximately 5.5 minutes. The system includes LabSolutions workstation software and BTU and Specific Gravity calculation software.

#### Analyzer Information

##### System Configuration:

Three valves / three packed columns and three capillary columns with one BID detector and one FID detector

##### Sample Information:

Permanent gas, C<sub>1</sub>-C<sub>6</sub>, H<sub>2</sub>S

##### Concentration Range:

No.	Name of Compound	Concentration Range		Remarks
		Low Conc.	High Conc.	
1	H <sub>2</sub>	0.01%	10.0%	BID
2	O <sub>2</sub>	0.01%	10.0%	BID
3	N <sub>2</sub>	0.01%	10.0%	BID
4	CO	0.01%	10.0%	BID
5	CO <sub>2</sub>	0.01%	10.0%	BID
6	C <sub>2</sub> H <sub>4</sub>	0.01%	10.0%	BID
7	C <sub>2</sub> H <sub>6</sub>	0.01%	10.0%	BID
8	C <sub>2</sub> H <sub>2</sub>	0.01%	10.0%	BID
9	H <sub>2</sub> S	0.01%	30.0%	BID
10	CH <sub>4</sub>	0.01%	80.0%	FID
11	C <sub>3</sub> H <sub>8</sub>	0.001%	5.0%	FID
13	C <sub>3</sub> H <sub>6</sub>	0.001%	5.0%	FID
14	i-C <sub>4</sub> H <sub>10</sub>	0.001%	1.0%	FID
15	n-C <sub>4</sub> H <sub>10</sub>	0.001%	1.0%	FID
16	C <sub>3</sub> H <sub>4</sub>	0.001%	1.0%	FID
17	C <sub>2</sub> H <sub>2</sub>	0.001%	1.0%	FID
18	trans-C <sub>4</sub> H <sub>8</sub>	0.001%	0.5%	FID
19	1-C <sub>4</sub> H <sub>8</sub>	0.001%	0.5%	FID
20	i-C <sub>4</sub> H <sub>8</sub>	0.001%	0.5%	FID
21	cis-C <sub>4</sub> H <sub>8</sub>	0.001%	0.5%	FID
22	i-C <sub>5</sub> H <sub>12</sub>	0.001%	0.5%	FID
23	n-C <sub>5</sub> H <sub>12</sub>	0.001%	0.5%	FID
24	1,3-C <sub>4</sub> H <sub>6</sub>	0.001%	0.5%	FID
25	C <sub>3</sub> H <sub>4</sub>	0.001%	0.5%	FID
26	C <sub>6</sub> plus	0.001%	0.5%	FID

Concentration range may vary depending on the sample. Please contact us for more consultation.

#### System Features

- Versatile software easy GC system operation
- Two channels with FID /BID detectors realizes high-speed analysis
- Linear response, simplifies calibration

Typical Chromatograms




Fig. 1 Chromatogram of BID




Fig. 2 Chromatogram of FID



➤ Please fill out the survey

## Related Products

Some products may be updated to newer models.



➤ **UFRGA Series**  
Nexis GC-2030 BIDUFRGA

## Related Solutions

➤ Price Inquiry

➤ Product Inquiry

➤ Technical Service /  
Support Inquiry

➤ Other Inquiry