

Application Data Sheet

No. 188

System Gas Chromatograph

Trace Oxygenated Hydrocarbons in Liquid Hydrocarbon Streams
 Nexis GC-2030OAS3
 GC-2014OAS3

This method is for determining trace oxygenated hydrocarbons in C4 liquefied petroleum gas (LPG) as described in below compound table. It requires the use of a dedicated gas chromatographic system which is configured with an automatic sampling and backflush technique in multiple columns.

Analyzer Information

System Configuration:

Two valves two SPL injectors / two capillary columns / two FID

Concentration Range:

No.	Name of Compound	Concentration Range	
		Low Conc.	High Conc.
1	Methyl Ether	1ppm	100,000ppm
2	Ethyl Methyl Ether	1ppm	100,000ppm
3	Ethyl Ether	1ppm	100,000ppm
4	Acetaldehyde	1ppm	100,000ppm
5	Methyl Formate	1ppm	100,000ppm
6	tert-Butyl Ethyl Ether	1ppm	100,000ppm
7	tert-Butyl Methyl Ether	1ppm	100,000ppm
8	Isopropyl Ether	1ppm	100,000ppm
9	Propylene Oxide	1ppm	100,000ppm
10	sec-Butyl Methyl Ether	1ppm	100,000ppm
11	Propionaldehyde	1ppm	100,000ppm
12	Butyl Methyl Ether	1ppm	100,000ppm
13	tert-Amyl Methyl Ether	1ppm	100,000ppm
14	n-Propyl Ether	1ppm	100,000ppm
15	Butyl Ethyl Ether	1ppm	100,000ppm
16	Isobutyraldehyde	1ppm	100,000ppm
17	Tetrahydrofuran	1ppm	100,000ppm
18	n-Butyraldehyde	1ppm	100,000ppm
19	Methyl Acetate	1ppm	100,000ppm
20	Tetrahydropyran	1ppm	100,000ppm
21	Trimethylacetaldehyde	1ppm	100,000ppm
22	Methanol	1ppm	100,000ppm
23	Acetone	1ppm	100,000ppm
24	2-Methylbutyraldehyde	1ppm	100,000ppm
25	Isovaleraldehyde	1ppm	100,000ppm
26	Cyclobutanone	1ppm	100,000ppm
27	Methyl Propionate	1ppm	100,000ppm
28	n-Valeraldehyde	1ppm	100,000ppm
29	2-Butanone	1ppm	100,000ppm
30	Ethanol	1ppm	100,000ppm

No.	Name of Compound	Concentration Range	
		Low Conc.	High Conc.
31	2-Ethylbutyraldehyde	1ppm	100,000ppm
32	3,3-Dimethylbutyraldehyde	1ppm	100,000ppm
33	2-Methylvaleraldehyde	1ppm	100,000ppm
34	Methyl Butyrate	1ppm	100,000ppm
35	1,4-Dioxane	1ppm	100,000ppm
36	Hexanal	1ppm	100,000ppm
37	3-Pentanone	1ppm	100,000ppm
38	3,3-Dimethyl-2-butanone	1ppm	100,000ppm
39	2-Pentanone	1ppm	100,000ppm
40	Isopropanol	1ppm	100,000ppm
41	n-Propanol	1ppm	100,000ppm
42	Cyclopropyl Methyl Ketone	1ppm	100,000ppm
43	2-Methyl-3-pentanone	1ppm	100,000ppm
44	3-Methyl-2-pentanone	1ppm	100,000ppm
45	Cyclopentanone	1ppm	100,000ppm
46	4-Methyl-2-pentanone	1ppm	100,000ppm
47	3-Hexanone	1ppm	100,000ppm
48	Isobutanol	1ppm	100,000ppm
49	tert-Butanol	1ppm	100,000ppm
50	sec-Butanol	1ppm	100,000ppm
51	Cyclobutanol	1ppm	100,000ppm
52	2-Hexanone	1ppm	100,000ppm
53	n-Butanol	1ppm	100,000ppm
54	3-Methyl-2-butanol	1ppm	100,000ppm
55	Neopentyl Alcohol	1ppm	100,000ppm
56	3-Pentanol	1ppm	100,000ppm
57	tert-Amyl Alcohol	1ppm	100,000ppm
58	2-Methyl-1-butanol	1ppm	100,000ppm
59	Cyclopentanol	1ppm	100,000ppm
60	2-Pentanol	1ppm	100,000ppm
61	3-Methyl-1-butanol	1ppm	100,000ppm
62	1-Pentanol	1ppm	100,000ppm

Detection limits may vary depending on the sample.
 Please contact us for more consultation.

Methods met:
UOP-960

System Features

- Dual FID channels (One is for detection of target compounds The other one is for checking backflush timing)
- Good repeatability

Typical Chromatograms




Fig. 1 Chromatogram of FID

ID#	Name
1	Methyl ether
2	Ethyl methyl ether
3	Ethyl ether
4	Acetaldehyde
5	Propionaldehyde
6	Isobutyraldehyde
7	n-Butyraldehyde
8	Methanol
9	Acetone
10	2-Methylbutyraldehyde
11	2-Butanone
12	Ethanol
13	Isopropanol + n-Propanol + Cyclopropyl methyl ketone
14	Isobutanol + tert-Butanol + sec-Butanol

Fig. 2 Compound List of The Chromatogram

First Edition: Jnovenber, 2017



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.
The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2017