

Application News

No. C127A

Liquid Chromatography Mass Spectrometry

High-Speed Analysis of Itraconazole in Plasma Using Triple Quadrupole LC/MS/MS (LCMS-8050)

Itraconazole is a triazole antifungal agent that is used widely for the treatment of infections by fungi including the genera *Trichophyton*, *Aspergillus* and *Candida*. This article introduces an example of high-speed analysis of itraconazole and its active metabolite hydroxy itraconazole in plasma using the LCMS-8050 high-sensitivity triple quadrupole mass spectrometer. Although a simple method of sample pretreatment was used that involves only deproteinization, the quantitative results obtained were excellent in terms of accuracy and precision.




Fig. 1 Structural Formula of Itraconazole and Hydroxy Itraconazole

The pretreatment method used is simple and requires no labor-intensive steps such as solid phase extraction, which allows pretreatment to be performed at low cost and in a short period of time.




Fig. 2 Pretreatment Workflow

Sample Pretreatment

The pretreatment workflow is shown in Fig. 2. Known concentrations of itraconazole and hydroxy itraconazole were solvated in a 1:1 water: methanol mixture, and used to prepare standard solutions for addition to plasma. These standard solutions for addition and an acetonitrile solution of an internal standard solution were added to human plasma, the deproteinized supernatant was diluted in an aqueous solution of 1 % acetic acid, and this mixture was used to prepare standard (STD) samples for a calibration curve and QC samples to validate the analysis results. The stable isotope-labeled compound of itraconazole (itraconazole-d3) was used as the internal standard solution.

LC/MS/MS Analytical Conditions

LC/MS/MS analytical conditions are shown in Table 1, and MRM transitions are shown in Table 2. A Shim-pack GIS column was used, and separation was performed in reverse-phase mode. Electrospray ionization (ESI) was used as the method of ionization, and MRM measurements were performed in positive ion mode.

Table 2 MRM Parameters

Compound Name	Polarity	MRM Transition
Itraconazole	+	705.40 > 392.40
Hydroxy itraconazole	+	721.40 > 408.40
Itraconazole-d3	+	710.40 > 397.35

Table 1 Analytical Conditions

Column	: Shim-pack GIS (75 mm L. × 2.1 mm I.D., 3 µm)	Injection Volume	: 3 µL
Mobile Phase	: A) 10 mmol/L Ammonium acetate in water B) Acetonitrile	Probe Voltage	: +4.0 kV (ESI-positive mode)
Time Program	: B. Conc. 65 % (0 - 1.00 min) → 95 % (1.50 - 2.50 min) → 65 % (2.51 - 4.50 min)	DL Temperature	: 300 °C
Analysis Time	: 4.5 min	Block Heater Temperature	: 400 °C
Flowrate	: 0.4 mL/min	Interface Temperature	: 400 °C
Column Temperature	: 40 °C	Nebulizing Gas Flow	: 3 L/min
		Drying Gas Flow	: 5 L/min
		Heating Gas Flow	: 15 L/min

■ Analysis Results

The calibration curves for itraconazole and hydroxy itraconazole in plasma are shown in Fig. 3, and representative chromatograms are shown in Fig. 4. Good linearity was obtained over the 10 to 1000 ng/mL range of sample concentration in plasma, with a correlation coefficient of 0.999 or higher.




Fig. 3 Calibration Curves for Each Compound

The accuracy and repeatability of each point on the calibration curves for itraconazole and hydroxy itraconazole concentration in plasma are shown in Table 3.

Good results were obtained for all points on the calibration curves, with repeatability of 15 % or below and accuracy within 100 ± 10 %.

Table 3 Measurement Results of Standard Samples in Plasma

		Concentration in Plasma (ng/mL)	Accuracy (%)	Concentration Repeatability (%)
Itraconazole	STD1	10	107.3	2.07
	STD2	50	98.9	1.52
	STD3	100	93.7	0.62
	STD4	500	99.1	0.30
	STD5	1000	101.1	0.59
Hydroxy Itraconazole	STD1	10	100.5	11.41
	STD2	50	98.3	4.47
	STD3	100	100.1	0.82
	STD4	500	102.1	1.15
	STD5	1000	99.0	5.55

QC sample measurement results are shown in Table 4. The results obtained were sufficiently precise and accurate, with repeatability of 5 % or lower and concentration accuracy of 80 to 105 %.

Table 4 Measurement Results of QC Samples

		Concentration in Plasma (ng/mL)	Accuracy (%)	Concentration Repeatability (%)
Itraconazole	QC1	25	98.7	0.65
	QC2	250	98.8	0.82
	QC3	750	103.3	0.38
Hydroxy Itraconazole	QC1	25	81.3	3.27
	QC2	250	90.4	4.36
	QC3	750	86.9	3.10




Fig. 4 Representative Chromatograms for Itraconazole and Hydroxy Itraconazole

This Application News was created with the cooperation of the Pharmaceutical Sciences Department of Tohoku University Hospital.

Notes: · The products mentioned in this article have not received approval for use as medical devices based on the Pharmaceutical and Medical Device Act.
· The analytical methods mentioned in this article cannot be used for diagnostic purposes.

Second Edition: Jun. 2016



For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2016

Shimadzu Corporation
www.shimadzu.com/an/

Related Products

Some products may be updated to newer models.



> **LCMS-8050**
High Performance Liquid
Chromatograph Triple Quadr...



> **Shim-pack GIS
Series**
HPLC Column Solutions

Related Solutions

> Price Inquiry

> Product Inquiry

> Technical Service /
Support Inquiry

> Other Inquiry