

Application News

No. **C207A**

Liquid Chromatograph Mass Spectrometry

Analysis of Residual Pesticides (No. 1: in Soybeans) Using Triple Quadrupole LC/MS/MS <LCMSTM-8060>

With a recent increase in the number of regulated pesticides, more effective methods for simultaneous analysis of residual pesticides in food are required.

QuEChERS, which was introduced by the United States Department of Agriculture (USDA) in 2003, is known as a quick and simple pretreatment method and approved as an official method by AOAC and CEN. This method requires no special instruments for extraction of pesticides, but the contaminants that cannot be completely removed by means of purification procedures may affect accurate quantitative analysis. In such cases, sample dilution or review of the purification process is also required.

This article introduces an example of the analysis of 158 pesticides among those specified in the Multi-residue Method I and II for Agricultural Chemicals by LC-MS (Agricultural Products)¹⁾ by measuring these pesticides in the sample solutions pretreated using the QuEChERS method, resulting in good recovery.

M. Kawashima, N. Kato

Sample Pretreatment

The soybean sample was pretreated using the QuEChERS method. The workflow of sample pretreatment is shown in Fig. 1. The concentration of samples extracted was 0.5 g/mL.

PL2005MIX-4, 5, 6, 7, 8, 9 and 10, mixtures of pesticides manufactured by Hayashi Pure Chemical Ind., Ltd., were used as the standard samples. The matrix effect was identified using the matrix standard solution (10 ng/mL pesticide in the solution) made by adding each pesticide to the sample solution pretreated with the QuEChERS method to reach a concentration of 0.02 mg/kg in the soybean extract.




Fig. 1 Pretreatment Workflow

Analytical conditions

The analytical conditions for HPLC and MS are shown in Table 1.

Table 1. Analytical Conditions

| [HPLC conditions] (Nexera TM X2) | |
|---|--|
| Column | : Shim-pack Scepter TM C18-120 (100 mm x 2.1 mm I.D., 3 µm) |
| Mobile phases | : A) 5 mM ammonium formate, 0.02% acetic acid in H ₂ O B) 5 mM ammonium formate, 0.02% acetic acid in MeOH |
| Gradient Program | : B 5% (0-2 min) – B 50% (5 min) – B 97% (13-16 min) – B 5% (16.1-20 min) |
| Flow rate | : 0.3 mL/min |
| Column Temp. | : 40°C |
| Injection volume | : 1 µL |
| [MS conditions] (LCMS-8060) | |
| Ionization | : ESI (Positive and negative mode) |
| Probe Voltage | : +2.0 kV / -1.5 kV |
| Mode | : MRM |
| Nebulizing gas flow | : 3.0 L/min |
| Drying gas flow | : 10.0 L/min |
| Heating gas flow | : 10.0 L/min |
| DL Temp. | : 200°C |
| Heat Block Temp. | : 300°C |
| Interface Temp. | : 200°C |
| Probe position | : +2.0 mm |

MRM Measurement of Matrix Standard Solution

Fig. 2 shows the MRM chromatogram of the matrix standard solution made by adding pesticide standard solution to the soybean extract.




Fig. 2 Example of Peak Detections of 158 Pesticides (Soybean Extract Added to 10 ng/mL Standard Solution)

Recovery

The recovery and peak area repeatabilities (n=6) of the matrix standard solutions for 158 pesticides were determined. The results of determination are shown in Table 2. Details of the recovery are shown in Fig. 3.

The recovery for 156 of 158 pesticides were in the range of 70 to 120%. Even in the test solution containing a high concentration of sample, 98.7% of these pesticides were not significantly affected by the matrix, resulting in good recovery and repeatabilities.




Table 2. Recovery and Peak Area Repeatability of Sample Solutions

| Compound name | Recovery(%) | %RSD | Compound name | Recovery(%) | %RSD | Compound name | Recovery(%) | %RSD |
|----------------------------|-------------|------|------------------------|-------------|-------|-------------------------|-------------|------|
| 1-Naphthaleneacetic Acid | 76.3 | 8.32 | Dymuron | 91.8 | 5.36 | Methoxyfenozide | 92.0 | 6.35 |
| 2,4-D | 97.3 | 4.66 | Epoxiconazole | 90.0 | 1.87 | Metosulam | 102.7 | 7.44 |
| 4-Chlorophenoxyacetic acid | 76.5 | 5.34 | Ethametsulfuron-methyl | 95.3 | 6.52 | Metsulfuron-methyl | 96.5 | 5.54 |
| Abamectin B1a | 93.5 | 1.37 | Ethoxysulfuron | 101.2 | 4.22 | Monolinuron | 96.0 | 2.52 |
| Acibenzolar-S-methyl | 90.0 | 8.84 | Fenamidone | 92.8 | 3.74 | Naproanilide | 91.3 | 5.53 |
| Acifluorfen | 88.7 | 7.50 | Fenhexamid | 92.0 | 6.06 | Naptalam | 95.9 | 9.97 |
| Aldicarb | 92.8 | 4.24 | Fenobucarb | 96.5 | 4.83 | Novaluron | 90.9 | 6.44 |
| Aldoxycarb | 96.5 | 1.15 | Fenoxaprop-ethyl | 88.0 | 2.68 | Oryzalin | 90.5 | 9.95 |
| Anilofos | 93.3 | 3.08 | Fenoxycarb | 92.9 | 3.94 | Oxamyl | 93.6 | 2.71 |
| Aramite | 95.8 | 3.88 | Fenpyroximate E | 93.4 | 3.11 | Oxaziclonemefone | 89.5 | 4.16 |
| Azamethiphos | 93.5 | 4.68 | Fenpyroximate Z | 93.9 | 2.90 | Oxycarboxin | 96.3 | 3.21 |
| Azimsulfuron | 84.9 | 8.46 | Ferimzone(E) | 95.2 | 2.66 | Pencycuron | 95.8 | 3.89 |
| Azinphos-methyl | 95.2 | 3.91 | Ferimzone(Z) | 96.9 | 2.08 | Penoxsulam | 99.9 | 2.92 |
| Azoxystrobin | 93.6 | 6.70 | Flazasulfuron | 97.3 | 4.79 | Pentoxazone | 79.9 | 9.20 |
| Bendiocarb | 99.2 | 2.14 | Florasulam | 97.6 | 7.50 | Phenmedipham | 95.5 | 1.35 |
| Bensulfuron-methyl | 97.9 | 5.57 | Fluazifop | 94.1 | 6.73 | Pirimicarb | 94.6 | 5.96 |
| Benzofenap | 97.4 | 1.42 | Flufenacet | 95.3 | 4.66 | Primisulfuron-methyl | 95.0 | 3.81 |
| Boscalid | 98.0 | 2.53 | Flufenoxuron | 93.1 | 7.23 | Propaquizafop | 93.6 | 3.13 |
| Bromoxynil | 92.7 | 8.96 | Flumetsulam | 101.8 | 6.45 | Propoxycarbazone | 142.5 | 8.39 |
| Butafenacil | 99.0 | 2.63 | Fluridone | 93.6 | 2.22 | Prosulfuron | 99.6 | 6.03 |
| Carbaryl(NAC) | 98.5 | 5.77 | Fluroxypyr | 91.2 | 8.65 | Pyraclostrobin | 96.5 | 4.62 |
| Carbofuran | 93.5 | 5.35 | Fomesafen | 103.4 | 4.83 | Pyrazolynate | 93.7 | 2.96 |
| Carpropamid | 94.4 | 3.72 | Foramsulfuron | 115.1 | 8.65 | Pyrazosulfuron-ethyl | 96.8 | 8.99 |
| Chloridazon | 92.1 | 2.82 | Forchlorfenuron | 92.0 | 5.85 | Pyrifthalid | 95.7 | 4.52 |
| Chlorimuron-ethyl | 101.8 | 7.76 | Furametpyr | 97.5 | 3.02 | Quizalofop-ethyl | 81.4 | 1.31 |
| Chloroxuron | 95.5 | 5.95 | Furathiocarb | 93.5 | 1.76 | Silafluofen | 84.9 | 6.25 |
| Chlorsulfuron | 96.9 | 5.36 | Gibberellic acid | 63.5 | 10.61 | Simeconazole | 95.1 | 2.63 |
| Chromafenozide | 95.1 | 1.65 | Halosulfuron-methyl | 80.2 | 7.45 | SpinosynA | 100.9 | 6.04 |
| Cinosulfuron | 98.4 | 9.43 | Haloxypop | 82.7 | 7.70 | SpinosynD | 105.9 | 4.00 |
| Clodinafop acid | 91.9 | 8.00 | Haloxypop | 85.0 | 8.78 | Sulfentrazone | 86.4 | 7.29 |
| Clofentezine | 84.7 | 4.06 | Hexaflumuron | 96.4 | 7.41 | Sulfosulfuron | 97.4 | 7.18 |
| Clomeprop | 87.6 | 3.77 | Hexythiazox | 93.8 | 3.31 | Tebufenozide | 95.3 | 5.07 |
| Cloprop | 97.8 | 9.36 | Imazalil | 106.8 | 3.52 | Tebuthiuron | 91.6 | 4.51 |
| Cloquintocet-mexyl | 97.8 | 3.64 | Imazaquin | 95.5 | 4.25 | Teflubenzuron | 87.9 | 7.57 |
| Cloransulam-methyl | 101.9 | 6.16 | Imazosulfuron | 94.0 | 5.67 | Tetrachlorvinphos | 94.2 | 3.80 |
| Clothianidin | 85.9 | 5.42 | Imidacloprid | 89.9 | 1.12 | Thiabendazole | 94.0 | 3.67 |
| Cumyluron | 98.5 | 2.16 | Indanofan | 94.3 | 3.07 | Thiacloprid | 94.3 | 1.85 |
| Cyazofamid | 95.7 | 1.39 | Indoxacarb | 99.9 | 4.16 | Thiamethoxam | 96.0 | 1.72 |
| Cyclanilide | 96.8 | 4.10 | Iodosulfuron-methyl | 93.0 | 7.59 | Thidiazuron | 82.8 | 7.17 |
| Cycloate | 94.9 | 3.31 | Ioxynil | 98.8 | 7.08 | Thifensulfuron-methyl | 96.6 | 6.43 |
| Cycloprothrin | 72.6 | 5.13 | Iprovalicarb | 95.6 | 3.46 | Thiodicarb | 95.8 | 2.97 |
| Cyclosulfamuron | 96.8 | 5.74 | Isoxaflutole | 92.8 | 6.43 | Tralkoxydim 1 | 104.0 | 5.25 |
| Cyflufenamid | 91.9 | 1.72 | Lactofen | 90.5 | 2.10 | Tralkoxydim 2 | 93.9 | 4.25 |
| Cyprodinil | 94.6 | 3.10 | Linuron | 95.4 | 3.54 | Triasulfuron | 96.8 | 3.77 |
| Diallate | 94.1 | 4.05 | Lufenuron | 93.2 | 4.48 | Tribenuron-methyl | 94.1 | 7.77 |
| Dichlorprop | 97.5 | 9.08 | MCPA | 96.1 | 4.16 | Triclopyr | 94.5 | 7.21 |
| Diclomezine | 100.7 | 8.89 | MCPB | 86.6 | 2.15 | Tridemorph 1 | 97.5 | 4.28 |
| Diclosulam | 95.9 | 2.23 | Mecoprop+Mecoprop-P | 85.2 | 2.79 | Tridemorph 2 | 96.3 | 2.18 |
| Diflubenzuron | 87.4 | 3.15 | Mepanipyrim | 94.5 | 3.97 | Trifloxysulfuron | 96.3 | 7.75 |
| Dimethirimol | 94.7 | 3.20 | Mesosulfuron-methyl | 95.1 | 3.50 | Triflurumuron | 92.9 | 3.70 |
| Dimethomorph(E) | 98.1 | 2.86 | Methabenzthiazuron | 96.6 | 2.11 | Triflurosulfuron-methyl | 99.5 | 5.49 |
| Dimethomorph(Z) | 98.1 | 2.86 | Methiocarb | 95.0 | 4.16 | Triticonazole | 94.2 | 2.68 |
| Diuron | 96.6 | 2.34 | Methomyl | 97.8 | 1.44 | | | |

1) Ministry of Health, Labour and Welfare: Testing Method of Agricultural Chemical Residues in Food, Feed Additives or Components of Animal Pharmaceuticals (PFSB/DFS Notification No. 1129002)P

LCMS, Nexera and Shim-pack Scepter are trademarks of Shimadzu Corporation in Japan and/or other countries.



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Shimadzu disclaims any proprietary interest in trademarks and trade names used in this publication other than its own. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

First Edition: Jan. 2020
Second Edition: Jun. 2020

Related Products

Some products may be updated to newer models.



> LCMS-8060

High Performance Liquid
Chromatograph Triple Quadr...

Related Solutions

> Price Inquiry

> Product Inquiry

> Technical Service /
Support Inquiry

> Other Inquiry