

# Application Data Sheet

## No.49

### System Gas Chromatograph

## Permanent Gas Analysis System Nexis GC-2030PGAS1 GC-2014PGAS1

This method is for determination of the permanent gases in various pure gases. Sample is introduced into two sample loops of two 10-port valves. One sample loop is for determination of He and H<sub>2</sub>, the other is for determination of Ar, O<sub>2</sub> and N<sub>2</sub>. Both 10 ports have Porapak-N backflush columns to remove CO, CO and hydrocarbons. Nitrogen gas is used for He and H<sub>2</sub> analysis, and Helium is used for Ar, O<sub>2</sub> and N<sub>2</sub> analysis. Two molecular sieve capillary columns carry out the separations. The analysis time is approximately 9 minutes. The system includes LabSolutions GC workstation software.

#### Analyzer Information

##### System Configuration:

Two valves / capillary column with Two TCD detectors

##### Sample Information:

He, H<sub>2</sub>, Ar, O<sub>2</sub>, N<sub>2</sub>

##### Concentration Range:

| No. | Name of Compound | Concentration Range |            |
|-----|------------------|---------------------|------------|
|     |                  | Low Conc.           | High Conc. |
| 1   | He               | 0.005%              | 10%        |
| 2   | H <sub>2</sub>   | 0.005%              | 10%        |
| 3   | Ar               | 0.005%              | 10%        |
| 4   | O <sub>2</sub>   | 0.005%              | 20%        |
| 5   | N <sub>2</sub>   | 0.005%              | 50%        |

Detection limits may vary depending on the sample. Please contact us for more consultation.

#### Typical Chromatograms

- 15 minutes analysis for hydrocarbons analysis can be carried out
- Single FID channel with split/splitless injector
- Liquid sample is measured through internal sample loop in the liquid sampling device

#### System Features




Fig. 1 Chromatogram of TCD




Fig. 2 Chromatogram of TCD