

Analysis of Liquid Fertilizer by ICPE™-9820

■ Introduction

Fertilizers contain many elements which are useful in crop growth in addition to nitrogen, potassium, and phosphorus, which are the principal ingredients. As an analysis method for fertilizers, Japan's Ministry of Agriculture, Forestry and Fisheries (MAFF) established Official Methods of Analysis of Fertilizers in 1992. However, because much time had passed since that standard was established, Food and Agricultural Materials Inspection Center (FAMIC) issued Testing Methods for Fertilizers, which introduced new analysis conditions and analysis methods suited to the times. Testing Methods for Fertilizers includes inductively coupled plasma (ICP) atomic emission spectrometry as an analysis method for water-soluble metallic ingredients ^{(1), (2)}.

This article introduces an example of analysis of the water-soluble ingredients in liquid fertilizers using a Shimadzu ICPE-9820 simultaneous ICP atomic emission spectrometer.

Y. Matsui

■ Analysis Samples

Four commercially-available liquid fertilizers were used in the analysis.

■ Sample Pretreatment

Based on Testing Methods for Fertilizers, 1 g of the sample was weighed, introduced into a 100 mL volumetric flask and diluted to the marked line with pure water, after which the solution was filtered with filter paper. A set amount of the filtrate was sampled and introduced into a 100 mL volumetric flask, and the measurement solution was prepared by adding hydrochloric acid (HCl) and yttrium (Y) as an internal standard element and diluting the solution to the marked line with pure water.

■ Instrument and Measurement Conditions

A Shimadzu ICPE-9820 simultaneous ICP atomic emission spectrometer was used in the measurements.

Table 1 shows the measurement conditions.

Table 1 Measurement Conditions

| | |
|-------------------------|----------------------------|
| Instrument | : ICPE-9820 |
| High frequency output | : 1.2 (kW) |
| Plasma gas flow rate | : 10 (L/min) |
| Auxiliary gas flow rate | : 0.6 (L/min) |
| Carrier gas flow rate | : 0.7 (L/min) |
| Nebulizer | : Nebulizer 10 |
| Chamber | : Cyclone chamber |
| Torch | : Mini-torch |
| View directions | : Radial (RD) / Axial (AX) |

■ Analysis

Table 2 shows the view direction and measured wavelength for each element. The measured wavelengths were based on Testing Methods for Fertilizers.

Table 2 View Direction and Measured Wavelength for Each Element

| Element | View direction | Measured wavelength (nm) |
|---------|----------------|--------------------------|
| B | Axial | 249.773 |
| Ca | Radial | 393.366 |
| Co | Axial | 228.616 |
| Cu | Axial | 327.396 |
| Fe | Axial | 259.940 |
| K | Radial | 766.490 |
| Mg | Radial | 279.553 |
| Mn | Axial | 257.610 |
| Mo | Axial | 202.030 |
| P | Radial | 178.287 |
| Zn | Axial | 213.856 |

■ Study of Conditions

Next, we studied the conditions for determining whether ionization interference occurs or not and conducting a simultaneous analysis of the wavelengths prescribed in Testing Methods for Fertilizers. When analyzing samples with high concentrations of alkali metal elements, ionization interference due to self-absorption requires attention, and it is known that ionization interference causes deflection of the calibration curve in axial view observation. Fig. 1 shows the calibration curves for radial (RD) view and axial (AX) view observation of K. Obvious deflection of the calibration curve can be seen in the result for AX observation. In RD view observation, ionization interference can be suppressed and the calibration curve displays satisfactory linearity. RD view observation also enables analysis at the prescribed wavelength in the case of Ca, which has good sensitivity. Thus, combined use of AX and RD view observation makes it possible to conduct simultaneous analysis of elements over the full range of concentrations from low to high without changing the dilution factor depending on elements.

■ Analysis Results

Table 3 shows the quantitation results for the respective samples. Table 4 shows the results of the dilution tests and spike and recovery tests of the samples. To check for interference, dilutions tests were conducted for samples with high concentrations of the contents, and spike and recovery tests were conducted for samples with low concentrations.

Satisfactory results were obtained for all samples, as the results of the dilution test and spike and recovery test were within 90% to 110% in all cases.




Fig. 1 Calibration Curves of K by Radial (RD) View and Axial (AX) View Observation

Conclusion

Eleven water-soluble elements in four liquid fertilizers were measured in accordance with Testing Methods for Fertilizers (2019). Use of the Shimadzu ICPE-9820 enabled simultaneous analysis over the full range of concentrations from low to high by combined use of high sensitivity axial (AX) view observation and radial (RD) view observation for analysis of high concentration ingredients. In addition, dilution tests and spike and recovery tests demonstrated that simultaneous analysis of low concentration ingredients is also possible without interference caused by high concentration ingredients.

<References>

- (1) Testing Methods for Fertilizers (2019) (Food and Agricultural Materials Inspection Center (Incorporated Administrative Agency))
- (2) Keisuke Aoyama, Simultaneous Determination of Water-Soluble Principal Ingredients (W-P₂O₅, W-K₂O, W-MgO, W-MnO and W-B₂O₃) in Liquid Fertilizer using Inductively Coupled Plasma-Optical Emission Spectrometry (ICP-OES) (FAMIC, Research Report of Fertilizer, No. 8 (2015))

Table 3 Analysis Results of Samples

Unit: mg/L

| Element | Liquid fertilizer A | Liquid fertilizer B | Liquid fertilizer C | Liquid fertilizer D | Limit of Quantitation 10 σ |
|---------|---------------------|---------------------|---------------------|---------------------|-----------------------------------|
| B | 0.180 | 0.022 | 0.004 | 0.061 | 0.002 |
| Ca | < | 0.22 | 0.04 | 0.06 | 0.02 |
| Co | < | < | < | < | 0.002 |
| Cu | 0.059 | < | < | < | 0.003 |
| Fe | 0.141 | 0.485 | 0.007 | 0.015 | 0.001 |
| K | 46 | 47 | 69 | 48 | 1 |
| Mg | 0.073 | 0.390 | 0.253 | 0.654 | 0.001 |
| Mn | 0.054 | 0.013 | 0.002 | 0.051 | 0.0003 |
| Mo | 0.062 | 0.008 | 0.005 | 0.007 | 0.004 |
| P | 54 | 48 | 31 | 48 | 1 |
| Zn | 0.062 | 0.005 | 0.004 | 0.008 | 0.001 |

< : Below LOQ

Table 4 Results of Spike and Recovery Tests and Dilution Tests of Samples

Unit: %

| Element | Liquid fertilizer A | Liquid fertilizer B | Liquid fertilizer C | Liquid fertilizer D |
|---------|---------------------|---------------------|---------------------|---------------------|
| B | 106 | 100 | 99 | 99 |
| Ca | 100 | 100 | 99 | 100 |
| Co | 104 | 106 | 102 | 102 |
| Cu | 104 | 102 | 99 | 100 |
| Fe | 108 | 98* | 104 | 103 |
| K | 102* | 101* | 99* | 102* |
| Mg | 99 | 100 | 100 | 100 |
| Mn | 107 | 105 | 103 | 104 |
| Mo | 104 | 101 | 100 | 98 |
| P | 99* | 95* | 99* | 99* |
| Zn | 106 | 104 | 102 | 103 |

* : Result of dilution test

ICPE is a trademark of Shimadzu Corporation in Japan and/or other countries.

First Edition: Feb. 2020



For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Shimadzu disclaims any proprietary interest in trademarks and trade names used in this publication other than its own. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Shimadzu Corporation

www.shimadzu.com/an/