

Application Data Sheet

No.83

GC-MS

Gas Chromatograph Mass Spectrometer

Simultaneous Analysis of 66 Pesticides by GC-MS Using Hydrogen Carrier Gas

Helium gas is typically used as a carrier gas for GC-MS. However, in recent years, the insufficient supply of helium gas and its increase in price have caused problems, prompting global efforts to address these problems. The use of hydrogen gas as an alternative carrier gas holds promise for addressing the issue.

This data sheet presents the application of hydrogen carrier gas to the simultaneous analysis of a standard sample of 66 pesticides using the GCMS-QP2010 Ultra, which features a Dual Inlet Turbo Molecular Pump.

Experimental

Standard samples of pesticides (from Wako Pure Chemical Industries, Ltd.) were diluted with methylene chloride to prepare 0.005 and 0.5 mg/L standard solution mixtures. Table 1 shows the analytical conditions.

Table 1 Analytical Conditions

GC-MS:		GCMS-QP2010 Ultra												
Column:		Rtx®-5MS (Length 30 m, 0.25 mm I.D., df = 0.25 µm) (Restek Corporation, P/N: 12623)												
Glass Liner:		Custom SkyLiner, Splitless (Restek Corporation, P/N: 567366)												
[GC]														
Injection Temp.:		250 °C								[MS]				
Column Oven Temp.:		80 °C (2 min) → (20 °C /min) → 180 °C → (5 °C /min) → 280 °C (3 min)								Interface Temp.:		250 °C		
Injection Mode:		Splitless (High Pressure Injection 250 kPa, 2.3 min)								Ion Source Temp.:		230 °C		
Flow Control Mode:		Linear Velocity (65 cm/sec)								[Scan Conditions]				
Injection Volume:		2 µL								Mass Range:		m/z 70 - 360		
Sampling Time:		2 min								Scan Event Time:		0.5 sec		
										[SIM Conditions]				
										Monitoring m/z:		See below.		
										SIM Event Time:		0.3 sec		
SIM Monitoring m/z														
Compound Name			m/z		Compound Name			m/z		Compound Name			m/z	
Dichlorvos			185.0 109.0		Chlorothalonil			266.0 264.0		Pendimethalin			252.0 281.0	
Dichlobenil			171.0 173.0		Iprobenfos			204.0 91.0		Methyl dymron			107.0 119.0	
Etridiazole			213.0 211.0		Bromobutide			120.0 119.0		Isofenphos			213.0 185.0	
Chloroneb			193.0 191.0		Terbucarb			220.0 205.0		Captan			79.0 117.0	
Isoproc carb			136.0 121.0		Simetryn			213.0 170.0		Dimepiperate			145.0 119.0	
Molinate			126.0 98.0		Tolclofos-methyl			265.0 125.0		Procymidone			96.0 283.0	
Fenobucarb			150.0 121.0		Alachlor			188.0 160.0		Methidathion			145.0 85.0	
Trifluralin			306.0 290.0		Metalaxyl			206.0 160.0		alpha-Endosulfan			241.0 195.0	
Benfluralin			292.0 276.0		Dithiopyr			354.0 306.0		Butamifos			286.0 200.0	
Pencycuron			125.0 180.0		Fenitrothion			277.0 260.0		Naproramide			128.0 72.0	
Dimethoate			125.0 87.0		Esprocarb			91.0 222.0		Flutolanil			173.0 281.0	
Simazine			201.0 186.0		Malathion			173.0 93.0		Isoprothiolane			189.0 118.0	
Atrazine			215.0 200.0		Thiobencarb			100.0 72.0		Pretilachlor			238.0 262.0	
Propyzamide			175.0 173.0		Fenthion			278.0 125.0		Buprofezin			105.0 175.0	
Pyroquilon			130.0 173.0		Chlorpyrifos			314.0 197.0		Isoxathion			177.0 105.0	
Diazinon			304.0 179.0		Phthalide			243.0 241.0		beta-Endosulfan			195.0 241.0	
Disulfoton			89.0 97.0		Dimethametryn			212.0 255.0		Mepronil			119.0 269.0	

Analytical Results

The mass spectral patterns of some compounds are known to change when hydrogen carrier gas is used. To evaluate this, the 0.5 mg/L standard solution was measured in scan mode and checked the resulting mass spectra. Fig. 1 shows the mass spectra for dichlorvos. No change was evident in the mass spectral pattern for the target components.

In addition, it has been reported that using hydrogen carrier gas results in lower sensitivity than when helium is used. To test this, the 0.005 mg/L standard solution was measured in SIM mode using both gases and then compared the S/N ratios. For all the components, the S/N ratio for hydrogen carrier gas was lower, primarily due to the baseline of the chromatogram being higher. Fig. 2 shows the SIM mass chromatograms for dichlorvos, for which the S/N ratio was particularly reduced.

Table 2 shows the repeatability results (n = 5) at a concentration of 0.005 mg/L. For all components, %RSD was less than 10 %.




Fig. 1 Mass Spectra for 0.5 mg/L Dichlorvos
(Top: hydrogen carrier gas, Bottom: helium carrier gas)




Fig. 2 SIM Mass Chromatograms for 0.005 mg/L Dichlorvos
(Top: hydrogen carrier gas, Bottom: helium carrier gas)

Table 2 Repeatability (n = 5, 0.005 mg/L)

Compound Name	%RSD	Compound Name	%RSD	Compound Name	%RSD	Compound Name	%RSD
Dichlorvos	3.96	Chlorothalonil	1.67	Pendimethalin	4.70	Chlorimethifen	4.18
Dichlobenil	2.19	Iprobenfos	7.26	Methyl dymron	4.93	Edifenphos	7.76
Etridiazole	2.47	Bromobutide	3.79	Isofenphos	8.79	Propiconazole-1	7.72
Chloroneb	2.61	Terbucarb	5.08	Captan	6.86	Endosulfan	3.28
Isoprocab	4.12	Simetryn	5.60	Dimepiperate	7.67	Propiconazole-2	7.08
Molinat	1.04	Tolclofos-methyl	2.59	Procymidone	6.65	Thenylchlor	3.85
Fenobucarb	3.34	Alachlor	5.80	Methidathion	5.84	Pyributicarb	5.60
Trifluralin	1.44	Metalaxyl	8.22	alpha-Endosulfan	9.82	Iprodione	4.20
Benfluralin	2.04	Dithiopyr	1.67	Butamifos	8.81	Pyridaphenthion	6.36
Pencycuron	3.25	Fenitrothion	5.68	Naproramide	9.34	EPN	7.40
Dimethoate	3.68	Esprocarb	7.13	Flutolanil	3.51	Piperophos	8.33
Simazine	6.79	Malathion	7.01	Isoprothiolane	9.47	Bifenox	9.05
Atrazine	5.35	Thiobencarb	2.40	Pretilachlor	7.38	Pyriproxyfen	1.71
Propyzamide	4.82	Fenthion	4.63	Buprofezin	9.01	Mefenacet	6.06
Pyroquilon	3.94	Chlorpyrifos	6.57	Isoxathion	3.54	Cafenstrole	7.55
Diazinon	7.06	Phthalide	2.41	beta-Endosulfan	2.97		
Disulfoton	7.54	Dimethametryn	2.80	Mepronil	5.13		

Summary

Using hydrogen carrier gas for the simultaneous analysis of 66 pesticides showed no change in the mass spectral pattern. Since the flow rate is high compared to when helium carrier gas is used, however, the baseline was raised, resulting in lower sensitivity. Nevertheless, the results demonstrated that the GCMS-QP2010 Ultra, featuring a Dual Inlet Turbo Molecular Pump, provides satisfactory sensitivity that allows application of hydrogen gas. However, be sure to confirm the required sensitivity and quantitative performance before replacing helium with hydrogen carrier gas.

Hydrogen gas is combustible, so it is necessary to handle it carefully. For safety measures, please refer to our web site: <http://www.shimadzu.com/an/gc/support/faq/bombe/bombe1.html>

First Edition: July, 2013



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2013