

Application Notes

DNA/ RNA Microchip Electrophoresis System (MultiNA)

H1N1 virus detection using Seeplex® FluA ACE subtyping kit

Introduction

H1N1 is currently the most striking issue in the public mind. Unlike seasonal influenza, H1N1 virus as the subtype of influenza A was classified by World Health Organization (WHO) as pandemic influenza as it has caused sustained community level outbreaks in multiple parts of the world. To overcome this issue, it is necessary to detect the exact subtypes of influenza A virus rapidly and reliably for the appropriate medicine to be given at early stage. Here we report on the application of MCE-202 (MultiNA) in combination with the Seeplex® FluA ACE subtyping kit to detect five different influenza virus A using DPO™ (Dual Priming Oligonucleotide) based multiplex PCR method (Seeplex® PCR).

Materials

DNA-1000 Reagent Kit from Shimadzu Corporation, Seeplex® FluA ACE Subtyping kit reagent from Seegene Corporation, SYBR® Gold nucleic acid gel stain from Invitrogen Corporation, DNA 1000 ladder (phiX174DNA-HaeIII digest) from Promega and TE buffer from Nacalai Tesque

Method



Specimen Collection (e.g.: Nasopharyngeal swab)

Nucleic acid isolation and Reverse transcription



DPO based multiplex PCR amplification



Electrophoresis of PCR samples on MultiNA

Results & Discussions




Fig. 1 Gel images and electropherograms of phiX174DNA-HaeIII digest (Promega), conc. 10ng/μL

Figure 1 showed reproducibility of MultiNA separation using phiX174DNA-HaeIII digest as ladder (multiple injections) with DNA-1000 kit. The internal upper and lower markers virtually eliminate chip to chip variation.




Fig. 2 Gel images and electropherograms of multiplex PCR products of five different influenza A virus using Seeplex® FluA ACE subtyping kit. DNA plasmid was used as the internal control. phiX174DNA-HaeIII digest (Promega) was used as the ladder.

Figure 2 showed that using the Seeplex/MultiNA, a panel of four different influenza virus A and one human common influenza virus were clearly detected and identified according to their subtyping, for example: 527 bp from Human FluA-H3, 440 bp from Human FluA-H1, 313 bp from Avian FluA-H5, 257 bp from Swine FluA-H1 and 204 bp from Universal Influenza A.

Authors: Djohan Kesuma, Chua Ai Ming and Zhan Zhaoqi