

Application Data Sheet

No. 8

GCMS

Gas Chromatograph Mass Spectrometer

Analysis of Brominated Flame Retardants Using a Pyrolysis GC-MS System

This datasheet introduces a sample measurement of Decabromodiphenyl ether (Deca BDE), a brominated flame retardant, using the EGA/PY-3030D Multi-Shot Pyrolyzer and the GCMS-QP2010 Ultra.

Analysis Conditions and Results

A 0.5 mg of shaved polystyrene (NMIJ CRM 8108-a) containing the brominated flame retardant was used as the measurement sample. The FASST (Fast Automated Scan/SIM Type) measurement mode, which is capable of simultaneous Scan and SIM measurements, was used. The analysis conditions are shown in Table 1.

Table 1: Analysis Conditions

Pyrolyzer	: EGA/PY-3030D Multi-Shot Pyrolyzer		
GC-MS	: GCMS-QP2010 Ultra		
Column	: Ultra ALLOY-PBDE (15 mL. × 0.25 mmI.D., 0.05 μm)		
[Pyrolyzer]			
Pyrolysis furnace temperature	: 200 °C -> (20 °C/min) -> 340 °C(3 min)		
Interface temperature	: Manual (320 °C)		
[GC]			
Vaporization chamber temperature	: 320 °C		
Column oven temperature	: 80 °C (1 min) -> (40°C/min) -> 200 °C -> (15°C/min) -> 330°C (5 min)		
Injection mode	: Split		
Carrier gas	: Helium		
Control mode	: Pressure (100kPa)		
Purge flow rate	: 3.0 ml/min		
Split ratio	: 20		
		[MS]	
		Interface temperature	: 320 °C
		Ion source temperature	: 230 °C
		Solvent elution time	: 1.5 min
		Tuning mode	: Normal
		Measurement mode	: FASST (simultaneous Scan/SIM measurements)
		Scan mass range	: <i>m/z</i> 35-1000
		Scan event time	: 0.3 sec
		SIM monitoring <i>m/z</i>	: 799.30, 801.30, 719.40, 721.40
		SIM event time	: 0.2 sec
		SIM micro-scan width	: 0.5 <i>u</i>

