

Application News

No. A622

Spectrophotometric Analysis

Measurement of Color, Bitterness, and Diacetyl Content of Beer Using UV-Vis Spectrophotometer

Among alcoholic beverages, beer is extremely popular and enjoys a strong market, and diverse beers with distinctive individual characters, exemplified by craft beers, have sold well in recent years. Although sensory evaluation is an essential technique, evaluations using various types of analytical instruments are also effective for ensuring the uniform quality of these beers.

The UV-Vis spectrophotometer is one useful instrument for beer analysis. The standard international methods for sensory evaluation of beer are the methods of the American Society of Brewing Chemists (ASBC) and the European Brewery Convention (EBC), and in Japan, the Brewery Convention of Japan (BCOJ), centering on major beer producers, has established analytical methods for beer.

Among the 12 inspection items in the ASBC Methods of Analysis, color (Beer-10) and bitterness (Beer-23) are inspected by UV-Vis spectrophotometry. UV-Vis spectrophotometry also makes it possible to analyze the diacetyl, iron, polyphenol, and other components contained in beer.

This article introduces measurements of the color, bitterness, and diacetyl content of commercially-available beers using a Shimadzu UV-1900i UV-Vis spectrophotometer.

K. Kawahara, K. Maruyama

■ Measurement Samples (Color and Bitterness)

Four types of commercially-available beer were prepared, and their color and bitterness were measured. Fig. 1 shows the appearance of the measurement samples, and Table 1 shows the details of each sample.



Fig. 1 Appearance of Measurement Samples

Table 1 Types and Features of Measurement Samples

	Type	Features
A	Belgian white	Color is extremely light, and has a spicy, acidic taste. Has a clean, refreshing taste.
B	Pilsner	World's most widely-consumed type. The main beers of major beer producers are this type. Yellow-gold color.
C	Indian pale ale (IPA)	Has a strong aroma and bitter taste due to the large amount of hops used. Alcohol content is also high.
D	Pale ale	Gold to copper color. Rich aroma is a distinctive feature.

■ Color

Beer has a wide range of colors, from the light-colored beers called white ale to black beers such as stout. The units used to express the color intensity of beer are the SRM (Standard Reference Method) system specified by the ASBC and the EBC (European Brewery Convention) system of the EBC. In these standards, the values are calculated using the absorbance of a single wavelength. Fig. 2 shows a reference diagram of the colors in these two systems.




Fig. 2 Reference Diagram of Colors in SRM and EBC Systems

The colors of the four beers shown in Fig. 1 were measured using a Shimadzu UV-1900i. The measurement method referred to ASBC Methods of Analysis, Beer Methods, Beer-10A⁽¹⁾ and the BCOJ Beer Analysis Methods, 8.15 Bitterness Unit⁽²⁾. Fig. 3 shows the flowchart of the measurement, and Fig. 4 shows the measurement instrument. Although absorbance value at 430 nm is used in color calculations, absorbance value at 700 nm is also measured to check turbidity.




Fig. 3 Flowchart of Color Measurement



Fig. 4 UV-1900i UV-Vis Spectrophotometer

Table 2 shows the absorbance values of the four samples at 430 nm and 700 nm and the results of the measurements of SRM and EBC calculated by Eq. 1 and Eq. 2, respectively. Due to the observable turbidity (cloudiness) of beer D, this beer was measured after filtration. From the appearance in Fig. 1, the color becomes darker in the order of A, B, C, D. The SRM and EBC values also increased in the same order.

Table 2 Results of Color Measurements of Samples

Sample	Absorbance value (430 nm)	Absorbance value (700 nm)	SRM value	EBC value
A	0.310	0.008	3.9	7.7
B	0.482	0.008	6.1	12.1
C	0.725	0.013	9.2	18.1
D	1.078 (after filtration)	0.051 (after filtration)	13.7	27.0

In addition to the single-point (single wavelength) measurement methods introduced here, the tristimulus method is also used to measure the color of beer. In the tristimulus measurement method, after measuring the transmission spectrum in the visible light region, L^* , a^* , and b^* are obtained by the CIELAB method. This method makes it possible to acquire color-related information that cannot be obtained by the single-point method.

■ Bitterness

Bitterness is one of the important components which determine the flavor of beer, and originates from the hops used to give beer its aroma. Alpha (α) acids are eluted from the hop flowers when the hops are boiled, and are then converted to iso- α acids, which are the bitterness components of beer.

In measurements of bitterness, the International Bitterness Units (IBU) of the sample is obtained by solvent elution of the bitterness components in the beer, followed by measurement of the absorbance value at 275 nm.

Here, the bitterness of the four types of beer shown in Fig. 1 was measured using the UV-1900i. The measurement method referred to the ASBC Methods of Analysis, Beer Methods, Beer-23A⁽³⁾. Fig. 5 shows the flowchart of the measurement.




Fig. 5 Flowchart of Bitterness Measurement

Table 3 shows the results of the bitterness measurements of the samples. Beer C had the highest value, and beer A had the lowest value. Beer C is an Indian pale ale (IPA), which has a strong aroma and bitterness due to the large amount of hops used. The measurement results also showed that this beer had the highest bitterness value in comparison with the other samples.

Table 3 Results of Bitterness Measurements of Samples

Sample	Absorbance value (275 nm)	Bitterness unit (BU)
A	0.253	12.5
B	0.661	33.0
C	1.146	57.5
D	0.838	42.0

■ Diacetyl Level

The unpleasant flavor that spoils the inherent good taste of beer is called "off-flavor" (4). This unpleasant off-flavor is known to occur when beer has high contents of diacetyl and 2,3-pentanedione (both compounds are included under the general name vicinal diketone: VDK). If the VDK content exceeds a certain level, the beer will have a buttery flavor called "butterscotch." However, VDKs can be reduced and deodorized by dehydrogenase in the presence of yeast. Therefore, in the beer production process, a maturation process in which VDKs are consumed by yeast is necessary after the fermentation process is completed. Thus, checking the VDK concentration after the fermentation process is completed is extremely important for the flavor of beer (5).

VDK evaluations are mainly conducted by sensory evaluation, but quantitative analysis of diacetyl, which is one VDK, is also possible with a UV-Vis spectrophotometer if pretreatment is conducted.

In this experiment, five types of beer which had completed the fermentation process were prepared, and changes in the concentration of diacetyl with time were measured using the UV-1900i. The measurement method referred to the ASBC Methods of Analysis, Beer Methods, Beer-25C (6). Fig. 6 shows the flowchart of the measurement.




Fig. 6 Flowchart of Diacetyl Measurement

Table 4 shows the measurement results of the standard samples, and Fig. 7 shows the calibration curve. Table 5 shows the measurement results for diacetyl in each sample. It was found that the concentration of diacetyl decreases with time after the fermentation process is completed, and the lowest concentration is reached before bottling. From this, it can be said that providing an adequate maturation period after completion of the fermentation process is necessary in order to reduce the diacetyl level in beer. (Beer is generally matured for 7 to 10 days at approximately the same temperature as the fermentation temperature⁽⁵⁾.)

Table 4 Measurement Results of Standard Samples

Sample concentration (mg/L)	Absorbance value (530 nm)
0.5	0.004
1.0	0.007
1.5	0.013
3.0	0.026
4.0	0.032




Fig. 7 Calibration Curve of Standard Samples

Table 5 Measurement Results of Diacetyl in Samples

Time after completion of main fermentation	Absorbance value (530 nm)	Diacetyl concentration (mg/L)
4 th day	0.023	2.81
5 th day	0.019	2.30
7 th day	0.015	1.90
Before bottling	0.007	0.74

Conclusion

The color, bitterness, and diacetyl concentration of beer samples were obtained by using a UV-1900i UV-Vis spectrophotometer. These values are extremely important indexes for maintaining the quality of beer. Use of a UV-Vis spectrophotometer enables quantitative evaluation of various beer indexes.

<Acknowledgement>

We wish to express our thanks to President Narihiro Suzuki and Mr. Motoki Sasaki of Nikenjyamochi Kadoya Honten (Ise Kadoya Brewery) for their cooperation in conducting these measurements.

<References>

- (1) American Society of Brewing Chemists, St. Paul, ASBC Methods of Analysis, Beer-10A, (1958, rev. 1975)
- (2) Brewers Association of Japan, Methods of Analysis of BCOJ (Revised Edition), Brewing Society of Japan, 8.15 Bitterness (2013)
- (3) American Society of Brewing Chemists, St. Paul, ASBC Methods of Analysis, Beer-23A, (1968, rev. 1975)
- (4) Narihiro Suzuki, "Beer Guy! Creating the World's Best Beer with Wild Yeast," Shinchosha (2019)
- (5) Toru Kishimoto, "The recent findings in the off-flavors of beer," J. Japan Association on Odor Environment, Vol. 44, No. 1, p. 13-20 (2013)
- (6) American Society of Brewing Chemists, St. Paul, ASBC Methods of Analysis, Beer-25C, (1964, rev. 1987)