

Simultaneous Analysis of Chlorate, Bromate, Perchlorate and Chlorite in Tap Water Using a Triple Quadrupole LC/MS/MS

Chlorate is derived from impurities in sodium hypochlorite, which is used as a disinfectant in the process of water purification or is generated in the same process. Bromate is also known to be derived from impurities or generated in the water purification process, especially when ozone disinfection is used.

Chlorate and bromate, which may have an adverse effect on human health, have been added to the water quality standards of Japan for tap water. The analytical method has been designated in "Method Determined by the Minister of Health, Labour and Welfare on the Basis of the Ordinance Provisions Relating to Water Quality Standards (Notification No. 261 issued by the Ministry of Health, Labour and Welfare of Japan (MHLW) in 2003)".

In addition to bromate, specified as an item contained in annex table 18-2 "Liquid Chromatography – Mass Spectrometry", chlorate, which had been specified as a component to be contained in Annex table 13 of "Simultaneous Analysis by Ion Chromatography", was added to the above annex table in response to the 2020 amendment, meaning that these components can now be analyzed simultaneously (Notification No. 95 issued by the MHLW in March 2020).

Using LCMS™-8060, this article introduces an example of simultaneous analysis of four components, consisting of perchlorate (as a tentative item under consideration) and chlorite (as a complementary item with a target value), in addition to chlorate and bromate for which the water quality control standards are already established, based on the method specified in "the annex table 18-2".

H. Niwa, H. Horiike




Fig. 1 MRM Chromatograms of (a) Chlorate, (b) Bromate, (c) Perchlorate and (d) Chlorite

Table 1 Analytical Conditions

Column	: Acclaim™ Trinity™ P1 (100 mm × 3.0 mm I.D., 3 μm)		
Mobile phases	: A 200 mmol/L ammonium acetate - water containing 0.5 % acetic acid B Acetonitrile		
Time program	: B.conc 95 % (0 min) → 90 % (3-10 min) → 5 % (10.5-15.5 min) → 95 % (16-22 min)		
Flow rate	: 0.4 mL/min		
Column temperature	: 35 °C		
Injection volume	: 1 μL		
Ionization	: ESI - Negative		
Probe voltage	: -1.0 kV		
DL temperature	: 150 °C		
Block heater temperature	: 400 °C		
Interface temperature	: 400 °C		
Nebulizing gas flow	: 2 L/min		
Drying gas flow	: 5 L/min		
Heating gas flow	: 20 L/min		
MRM transition	Chlorate ion	<i>m/z</i> 82.95 > 51.00 (Quantifier ion)	82.95 > 66.95 (Qualifier ion)
	Bromate ion	<i>m/z</i> 126.90 > 110.95 (Quantifier ion)	128.90 > 113.00 (Qualifier ion)
	Chlorite ion	<i>m/z</i> 66.95 > 50.95 (Quantifier ion)	66.95 > 34.95 (Qualifier ion)
	Perchlorate ion	<i>m/z</i> 98.95 > 66.95 (Quantifier ion)	98.95 > 82.95 (Qualifier ion)

■ MRM Chromatograms of Chlorate, Bromate, Perchlorate and Chlorite (Analysis of Standard Mixture)

Fig. 1 shows the MRM chromatograms obtained from a standard mixture of chlorate, bromate, perchlorate and chlorite at concentrations of 0.005 mg/L (standard value: 0.6 mg/L), 0.0005 mg/L (standard value: 0.01 mg/L), 0.0005 mg/L (target value: 0.025 mg/L) and 0.03 mg/L (target value: 0.6 mg/L), respectively. The chromatograms indicate that analysis of these components is possible at concentrations of 1/10 or lower of these standard or target values.

The analytical conditions are shown in Table 1.



LCMS™-8060

Linearity of Calibration Curves for Chlorate, Bromate, Perchlorate and Chlorite (Analysis of Standard Mixture)

Fig. 2 shows calibration curves created by an external standard method for chlorate, bromate, perchlorate and chlorite in the concentration ranges of 0.005 - 0.05 mg/L (4 points), 0.0005 - 0.005 mg/L (4 points), 0.0005 - 0.005 mg/L (4 points) and 0.03 - 0.6 mg/L (5 points), respectively. Good linearity was obtained with correlation coefficients R of > 0.999 for all calibration curves.

Analysis of Tap Water

A spike-and-recovery test was performed using tap water (from Kanagawa Prefecture). Collected tap water, which was dechlorinated by adding ethylenediamine, was spiked with chlorate, bromate, perchlorate and chlorite to prepare standard solutions with volumes of 10 µg/L, 1 µg/L, 1 µg/L and 60 µg/L, respectively. Fig. 3 shows MRM chromatograms of blank tap water and spiked test solutions obtained from analysis of these samples. Additionally, Table 2 shows the quantitative results and the recoveries and peak area reproducibilities (%RSDs) for the samples.

Conclusion

The use of LCMS in accordance with the method specified in "annex table 18-2" enabled a precise analysis of chlorate and bromate that have been added to the water quality standards, and perchlorate (a tentative item under consideration) and chlorite (a complementary item for water quality control with a target value).




Fig. 2 Calibration Curves of (a) Chlorate, (b) Bromate, (c) Perchlorate and (d) Chlorite




Fig. 3 MRM Chromatograms of Tap Water (Blank) and Spiked Samples

Table 2 Spike-and-Recovery Test Results Using Tap Water

Component	Spiked conc. (µg/L)	Quantitative conc. of tap water (blank) (µg/L)	Quantitative conc. of spiked samples (N=5, µg/L)	Recovery (N=5)	%RSD (N=5)
Chlorate	10	29.16	41.19	120 %	2.12
Bromate	1	0.10	1.06	96 %	6.32
Perchlorate	1	0.10	1.18	108 %	3.77
Chlorite	60	N.D.	63.06	105 %	2.85

LCMS is a trademark of Shimadzu Corporation in Japan and/or other countries. ACCLAIM and Trinity are either trademarks or registered trademarks of DIONEX CORPORATION in the United States and/or other countries. Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. These trademarks may not be accompanied by the symbol "TM" or "®" in this publication.

First Edition: Jul. 2020



For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Shimadzu disclaims any proprietary interest in trademarks and trade names used in this publication other than its own. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Shimadzu Corporation

www.shimadzu.com/an/