

## Application News

# No. M299

### Gas Chromatography Mass Spectrometry

## Evaluation of Food Deterioration by Multivariate Analysis

During the development of food products, physiochemical tests, microbial tests, and sensory evaluations are conducted, and the shelf life of the food is set based on their objective indexes. In markets characterized by a large number of products and short product cycles, like at present, safety and quality are adequately ensured during the shelf life of products, but it is difficult to grasp the distinctive features of deterioration of each individual product. In addition, consumers also have various understandings of shelf life after the product is opened. For these reasons, methods for predicting deterioration over time are demanded.

Deterioration of food is caused by microorganisms, water, oxygen, temperature, and light, resulting in decomposition of proteins, carbohydrates, and fats, into low molecular compounds. For example, it is known that proteins decompose into peptides and amino acids and then form amines and organic acids.

The purpose of this article was to track the condition of deterioration of milk after the product was opened. Amino acids, organic acids, sugars, and other hydrophilic components were identified using the Smart Metabolites Database™, and changes over time in those components were analyzed with the Traverse MS multivariate analysis software program.

K. Kawakita, Y. Sakamoto

### ■ Sample Preparation and Analysis Conditions

Commercially-available milk (paper carton) was prepared as the measurement sample. In order to understand changes due to elapsed time, samples were taken on day 0 (date of opening) and 3, 7, 14, and 21 days after opening. The milk was stored in a refrigerator except when removed for sampling.


For measurements, 100 µL of milk was measured out in a 1.5 mL microtube, and 10 µL of aqueous Ribitol (0.2 mg/mL) was added as an internal standard. Next, 500 µL of a mixed solvent of water:methanol:chloroform = 1:2.5:1 was added as an extractant for the hydrophilic components. After shaking for 30 min at 37 °C, the sample was fractionated by centrifugation (4 °C, 3,000 g, 10 min), and 450 µL of the supernatant (water/methanol phase) was collected. 400 µL of ultrapure water was added to the supernatant, centrifugation was conducted again for 5 min, and 500 µL of the supernatant was taken. After vaporizing the methanol with a centrifugal evaporator, the sample was sufficiently exsiccated by freeze-drying. 200 µL of a methoxyamine-pyridine solution (20 mg/mL) was added to the residue, and the sample was shaken for 90 min at 30 °C. Following this, 100 µL of N-Methyl-N-(trimethylsilyl) trifluoroacetamide (MSTFA) was added, the solution was shaken for 30 min at 37 °C. The analysis sample was measured by the Smart Metabolites Database, which contains MRM transitions for 475 metabolites. Table 1 shows the measurement conditions.

**Table 1 Measurement Conditions**

|                          |  |
|--------------------------|--|
| GC-MS                    | : GCMS-TQ™8040 NX                                    |
| Autoinjector             | : AOC™-20i + s                                       |
| Column                   | : BPX-5<br>(Length 30 m, 0.25 mm I.D., df = 0.25 µm) |
| [GC]                     |  |
| Vaporizing chamber temp. | : 250 °C   |
| Column oven temp.        | : 60 °C (2 min) → (15 °C/min) → 330 °C (3 min)       |
| Injection mode           | : Split  |
| Split ratio              | : 30   |
| Carrier gas              | : He   |
| Carrier gas control mode | : 39.0 cm/s (Constant linear velocity)               |
| Injection volume         | : 1 µL   |
| [MS]                     |  |
| Ion source temp.         | : 200 °C   |
| Interface temp.          | : 280 °C   |
| Data sampling mode       | : MRM  |
| Loop time                | : 0.25 s   |

### ■ Results

Fig. 1 shows the TICC (total ion current chromatogram) of the milk immediately after opening. A total of 106 components (94 compounds), including lactose and other sugars, amino acids and saturated fatty acids such as palmitic acid, were detected. The detected compounds are listed in Table 2.



**Fig. 1 Analysis Result of Milk Immediately After Opening**

Table 2 Compounds Detected from Milk

|                                   |                        |  |                                   |                            |                         |                      |
|-----------------------------------|------------------------|--|-----------------------------------|----------------------------|-------------------------|----------------------|
| Sugars/sugar alcohols/sugar acids | 2-Deoxy-glucose        | Arabinose                              | Arabitol                          | Erythrulose                | Fructose                | Fucose               |
|                                   | Galactose              | Galacturonic acid                      | Gluconic acid                     | Glucose                    | Glucuronic acid         | Glyceric acid        |
|                                   | Glycerol               | Inositol                               | Lactitol                          | Lactose                    | Maltose                 | Mannose              |
|                                   | N-Acetylmannosamine    | Ribitol (internal standard)            | Ribonic acid                      | Ribose                     | Ribulose                | Sucrose              |
|                                   | Threonic acid          | Xylitol                                | Xylose                            | Xylulose                   |                         |                      |
| Sugar phosphates                  | 3-Phosphoglyceric acid | Dihydroxyacetone phosphate             | Fructose 6-phosphate              | Glucose 6-phosphate        | Glycerol 2-phosphate    | Glycerol 3-phosphate |
|                                   | Mannose 6-phosphate    | Ribose 5-phosphate                     | Ribulose 5-phosphate              | Sedoheptulose 7-phosphate  |                         |                      |
| Amino acids                       | 2-Aminopimelic acid    | 3-Aminoglutaric acid (β-Glutamic acid) | 3-Aminopropanoic acid (β-Alanine) | 4-Aminobutyric acid (GABA) | 4-Hydroxyproline        | 5-Oxoproline         |
|                                   | Alanine                | Aspartic acid                          | Creatinine                        | Cysteine                   | Dimethylglycine         | Glutamic acid        |
|                                   | Glycine                | Hypotaurine                            | Isoleucine                        | Leucine                    | Lysine                  | Phenylalanine        |
|                                   | Proline                | Serine                                 | Threonine                         | Tyrosine                   | Ureidosuccinic acid     | Valine               |
| Saturated fatty acids             | Caproic acid           | Decanoic acid                          | Lauric acid                       | Myristic acid              | Nonanoic acid           | Octanoic acid        |
|                                   | Palmitic acid          | Stearic acid                           |                                   |                            |                         |                      |
| Organic acids                     | 2-Ketoglutaric acid    | 2-Hydroxyglutaric acid                 | 3-Hydroxybutyric acid             | 3-Hydroxyisovaleric acid   | 3-Hydroxypropionic acid | Aconitic acid        |
|                                   | Benzoic acid           | Citric acid                            | Fumaric acid                      | Glycolic acid              | Glyoxylic acid          | Lactic acid          |
|                                   | Malic acid             | Pyruvic acid                           | Succinic acid                     |                            |                         |                      |
| Others                            | 2-Aminoethanol         | Allantoin                              | Dopamine                          | O-Phosphoethanolamine      | Orotic acid             | Pantothenic acid     |
|                                   | Phosphoric acid        | Uracil                                 | Urea                              | Xanthosine monophosphate   |                         |                      |

Multivariate Analysis

Traverse MS (Reifycs Inc.) multivariate analysis software is an integrated analysis program for realizing high speed multivariate analysis using MRM data acquired by GC-MS/MS. Traverse MS is the real solution as it allows you to carry out peak integration and identification quickly, accurately, and easily by using a novel peak recognition algorithm developed specifically for MRM analysis data and a highly-operable graphical user interface. In order to evaluate variations in metabolism in vast amounts of MRM data sets, Traverse MS integrates statistical analysis tools including graphical display of peak areas of sample groups, principal component analysis, hierarchical clustering analysis, and metabolic pathway analysis. The qgd files (MRM data) acquired by GCMSolution™ can also be read directly into Traverse MS.

An analysis by Traverse MS was conducted targeting the above-mentioned 106 components detected from the milk. The analysis was carried out with 3 samples from each sampling day after opening, and the intensities of the compounds detected in each sample were normalized by using Ribitol as an internal standard. Fig. 2 shows the result of PCA (Principal Component Analysis) score plot by Traverse MS. Focusing on the first principal component (horizontal axis, PC1), the results are arranged in the order of day 0, day 3, day 7, day 14, and day 21 from the positive (right) side to the negative (left) side, and the results for days 0 and 3 are distributed on the positive side, while those for days 14 and 21 are on the negative side. It was suggested that the first principal component expresses changes over time.






Fig. 2 Distribution of Samples by PCA Score Plot



**Fig. 4 Results of Analysis of Variance (ANOVA)**

**0: Time of Opening, 1: 3 Days After Opening, 2: 7 Days After Opening, 3: 14 Days After Opening, 4: 21 Days After Opening**  
**(Only some components with  $p < 0.01$  are shown. The bars graphs show the class average, and the error bars indicate the standard deviation in the class.)**

## Conclusion

A simultaneous analysis of the hydrophilic components in food was possible by using the Shimadzu Smart Metabolites Database, and components that increased or decreased during deterioration could be identified by conducting widely targeted multivariate analysis using the analysis results.

This suggested the possibility that some organic acids, amino acids, and sugars can be used as markers which provide indexes of the freshness or deterioration of foods. Use in combination with objective evaluations for product development is expected to provide a more detailed understanding of the features of food products.

Smart Metabolites Database, GCMS-TQ, and GCMSsolution are trademarks of Shimadzu Corporation in Japan and/or other countries. AOC is a trademark of Shimadzu Corporation.

First Edition: Sep. 2020



**For Research Use Only. Not for use in diagnostic procedures.**

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Shimadzu disclaims any proprietary interest in trademarks and trade names used in this publication other than its own. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2020

Shimadzu Corporation

[www.shimadzu.com/an/](http://www.shimadzu.com/an/)

## Related Products

Some products may be updated to newer models.



> GCMS-TQ™8040

NX

Triple Quadrupole Gas Chromatograph  
Mass Spectrometry

## Related Solutions

> Food Metabolomics

> Price Inquiry

> Product Inquiry

> Technical Service /  
Support Inquiry

> Other Inquiry