

Application News

No. X273

X-Ray Analysis

ICH Q3D Elemental Impurities Analysis of Tablets by EDX - Verification Based on USP <233> ELEMENTAL IMPURITIES-PROCEDURES -

In the United States, the United States Pharmacopoeia General Test Chapters USP <232> and <233> have been applied to new drug products since January 1, 2018. ICP-MS and ICP-AES are recommended as the analysis procedures in the chapters. However, if the validation requirements are met, the alternative procedure can be substituted for the recommended analysis procedures^{(1) (2) (3)}.

Therefore, the appropriateness of Energy Dispersive X-ray Fluorescence Spectrometer was verified referring to "Limit Procedures" in USP <233>^{*1}.

The used instrument was EDX-7000 and the test was conducted by evaluating of elemental impurities in oral drug (tablets). As the target concentration (allowable concentration), the value obtained by dividing 30% of the PDE^{*2} value by the maximum daily dose was set.

*1 Since the measurement sample in EDX is "Powder", the "Solution" part described in the test method was replaced with "powder"

*2 PDE: Permitted Daily Exposure

T. Nakao, H. Nakamura

Elements

The measurement elements conformed to the method described in ICH Q3D. The seven elements of "Class 1 and Class 2A" are essential for risk assessment. The intentionally added element was assumed to be Pd in "Class 2B".

- Class1 : As, Hg, Pb, Cd
- Class2A : V, Co, Ni
- Class2B : Pd

Samples

Standardization : Cellulose powder spiked with standard solution for ICP

Test sample : as shown in Table 1

Table 1 Sample details under test

Tablets	Oral solid dosage form (Ethical drug)
Active ingredient content	30 mg / 0.5 g
Dosage form	Plain tablet
Daily amount of drug product	0.5 g
Main Component	Hydroxypropylcellulose

Target Concentration

The target concentration was 30% of the PDE values divided by daily amount of drug product. These values for each element are shown in Table 2.

Table 2 PDE Values and Target Concentrations of Oral Preparation

Elements	Unit	As	Hg	Pb,Cd	Co	V,Pd	Ni
PDE Values of Oral Preparation	μg/day	15	30	5	50	100	200
Target Concentrations	μg/g	9	18	3	30	60	120

Sample Preparation

Test samples were crushed into powders. Standardization, Standard powders, and Spiked sample powders 1 and 2 were prepared by spiked standard solutions for ICP. These details are shown in Table 3.

Table 3 Samples Used for Verification

Sample Name	Detail	Spiked Concentration	Number Produced (n)
Standardization	Cellulose Powder	0J, 0.5J, 1.5J	Each 1
Standard Powder	Cellulose Powder	1.0J	3
Spiked Sample Powder 1	Test sample	Target Concentration	3, 6
Spiked Sample Powder 2	Test sample	80% of Target Concentration	3
Unspiked Sample Powder	Test sample	Unspiked	3, 6

Sample Pretreatment

Each samples were introduced into a sample container lined with a polypropylene film. These samples are shown in Fig. 1.



Fig. 1 Standard Powder (Left), Spiked Sample Powder1 (Right)

Calibration Curve

The calibration curve was created with 3 points of 0J, 0.5J and 1.5J. Fig. 2 shows the calibration curve and correlation coefficient (R). The correlation coefficient R is 0.999 or more, and good linearity is obtained.





Fig. 2 Calibration Curve and Correlation Coefficient (R)

Table 4 Validation Requirements, Acceptance Criteria, and Results

Validation Requirements	Procedure	Acceptance Criteria	Results	Judgment
Detectability	(1) Standard powder: n = 3 measured 3 times each Spiked sample powder 1: n = 3 measured 3 times each	The average value of Spiked sample powder1 should be within $\pm 15\%$ of the average value of Standard powder	【Table 5 Results (1)】	Pass
	(2) Spiked sample powder 2: n = 3 measured 3 times each	Average value of Spiked sample powder 2 < Average value of Standard powder	【Table 5 Results (2)】	Pass
Specificity	- Comparison with Unspiked samples - Spectrum confirmation - Matrix components and coexisting elements: removal and correction of overlapping effects by them	Specific detection for matrix components and coexisting elements (Meeting detection sensitivity requirements)	【Fig. 3】 - Applying overlap correction for [As] for [Pb], and [Co] for [Fe] - Compared with Unspiked sample, each element spectral peak of added sample 1 is clear	Pass
Repeatability	Measure 6 of Spiked sample powder 1	Relative standard deviation (RSD) $\leq 20\%$	【Table 6】	Pass

Table 5 Detectability

Elements	As	Hg	Pb	Cd	V	Co	Ni	Pd	Judgment
Spiked Concentration	9.0	18.0	3.0	3.0	60	30	120	60	
(A) Standard Powder 1.0J	9.0	18.2	3.0	2.9	58.5	29.9	120.5	59.1	Pass
(B) Unspiked Sample Powder	<0.47	<0.26	<0.63	<1.15	<2.25	<1.17	<0.52	<0.63	
(C) Spiked Sample Powder 1 (Target Concentration)	9.3	17.8	3.1	2.9	56.7	29.3	117.9	59.5	
Results(1) [(C)/(A) -1] x100 [%]	+3.6	-2.2	+2.6	0.0	-3.1	-2.0	-2.2	+0.7	
Spiked Concentration	7.2	14.4	2.4	2.4	48	24	96	48	Pass
(D) Spiked Sample Powder 2 (80% of Target Concentration)	7.5	13.9	2.4	2.1	45.9	23.3	95.1	48.2	
Results(2) Relationship of (D)<(A)	<	<	<	<	<	<	<	<	

**Fig. 3 Specificity****Table 6 Repeatability**

Elements	As	Hg	Pb	Cd	V	Co	Ni	Pd	Judgment
Average of quantitative values by continuous repeated measurement	9.3	17.5	3.0	2.9	56.1	28.9	119.0	60.4	Pass
Standard Deviation	0.09	0.09	0.06	0.19	0.83	0.22	0.76	0.32	
RSD [%]	0.9	0.5	1.9	6.6	1.5	0.8	0.6	0.5	

Validation Results

Validation requirements, procedure, acceptance criteria, and results are shown in Table 4. Each result's details are shown in Table 5 and 6, and Fig. 3.

Conclusion

EDX-7000 was verified as the alternative procedure referring to "Limit Procedures" in USP <233>. It can be applied to the management of formulations and drug substances with similar compositions.

Along with the recommended ICP-MS and ICP-AES, cost reduction can be expected by using EDX according to the type and dosage of the drug product.

Table 7 Measurement Condition

Instrument	: EDX-7000, Sample Turret (option)
Elements	: As, Hg, Pb, Cd, V, Co, Ni, Pd
Analysis group	: Quantitative
Detector	: SDD
X-ray tube	: Rh target
Tube voltage	: 50 [kV]
Tube current	: Auto [μA]
Collimator	: 10 [mmφ]
Primary filter	: #1 (Cd, Pd), #2(V), #4 (As, Hg, Pb, Co, Ni)
Atmosphere	: Air
Integral time	: 1,800 [s] x 3 (#1, #2, #4)
Dead time	: Max.30 [%]

References

- (1) <232> Elemental Impurities-Limits
- (2) <233> Elemental Impurities-Procedures
- (3) Application News No.X271
ICH Q3D Elemental Impurities Analysis of Drug Substances by EDX
- (4) ICH HARMONISED GUIDELINE, GUIDELINE FOR ELEMENTAL IMPURITIES Q3D(R1) (Final version Adopted on 22 March 2019)

First Edition: Apr. 2020



For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Shimadzu disclaims any proprietary interest in trademarks and trade names used in this publication other than its own. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2020

Shimadzu Corporation

www.shimadzu.com/an/

Related Products

Some products may be updated to newer models.



> **EDX-8100**
Energy Dispersive X-ray Fluorescence Spectrometer

Related Solutions

> Food Dietary Restrictions

> Small Molecule Pharmaceutical

> Life Science

> Price Inquiry

> Product Inquiry

> Technical Service / Support Inquiry

> Other Inquiry