

Application News

No. G297

Gas Chromatography

Comparison of Separation by MICROPACKED-ST Columns

Barrier discharge ionization detectors (BID) are capable of detecting almost all compounds, excluding helium and neon, with a higher sensitivity compared to general detectors such as thermal conductivity detectors (TCD) and flame-ionization detectors (FID). When analyzing inorganic gases together with small hydrocarbons, simultaneous analysis with high sensitivity is possible by using a BID as the detector and the MICROPACKED-ST as the column. Until now we have introduced a number of example analyses using MICROPACKED-ST columns and now recently, 1.0 m and 3.0 m long columns have become available, enhancing the applicability of the columns.

This article introduces examples of separation using MICROPACKED-ST columns of different lengths: 1.0 m, 2.0 m, and 3.0 m.

S. Uchiyama

Column Information and Analysis Conditions

A MICROPACKED-ST column is a stainless steel tube with an inner diameter of 1.0 mm and is packed with SHINCARBON ST. Therefore, carrier gas can be controlled by inputting information of a capillary column with an approximated flow resistance.

Table 1 shows the common analysis conditions and Tables 2 to 4 show the analysis conditions for MICROPACKED-ST columns 1.0 m, 2.0 m, and 3.0 m long respectively. Although BIDs are easily influenced by changes in the carrier gas flow rate, the Nexis GC-2030 is capable of constant flow rate control, making it possible to maintain a constant flow rate with no need to create a pressure program.

Table 1 Common Analysis Conditions

Model	: Nexis GC-2030
Detector	: BID-2030
Inj. Mode	: Split 1:4
Inj. Temp.	: 150 °C
Carrier Gas	: 7 mL/min
Purge flow	: 3 mL/min
Det. Temp.	: 280 °C
Discharge Gas	: 50 mL/min (He)
Inj. Volume	: 1.0 mL (Using MGS-2030) 50 µL (Using a gas-tight syringe)

Table 2 Analysis Conditions for the 1.0 m MICROPACKED-ST

Column	: MICROPACKED-ST 1.0 m × 1.0 mm I.D. (Input 125 m × 0.50 mm I.D. and df = 15 µm for flow rate calculation)
Carrier Gas *	: He, 152.8 kPa (2.0 min) - 21.7 kPa/min - 242.2 kPa (0 min) - 13.5 kPa/min - 269.3 kPa (2.87 min) (Pressure Mode)
Column Temp.	: 35 °C (2.0 min) - 40 °C/min - 200 °C (0 min) - 25 °C/min - 250 °C (2.87 min)

Table 3 Analysis Conditions for the 2.0 m MICROPACKED-ST

Column	: MICROPACKED-ST 2.0 m × 1.0 mm I.D. (Input 250 m × 0.50 mm I.D. and df = 10 µm for flow rate calculation)
Carrier Gas *	: He, 226.8 kPa (2.5 min) - 15.2 kPa/min - 390.1 kPa (9.95 min) - 11.2 kPa/min - 405.1 kPa (5.42 min) (Pressure Mode)
Column Temp.	: 35 °C (2.5 min) - 20 °C/min - 250 °C (0 min) - 15 °C/min - 270 °C (5.42 min)

Table 4 Analysis Conditions for the 3.0 m MICROPACKED-ST

Column	: MICROPACKED-ST 3.0 m × 1.0 mm I.D. (Input 250 m × 0.50 mm I.D. and df = 30 µm for flow rate calculation)
Carrier Gas *	: He, 282.0 kPa (4.5 min) - 9.1 kPa/min - 450 kPa (9.98 min) (Pressure Mode)
Column Temp.	: 35 °C (4.5 min) - 10 °C/min - 270 °C (10 min)

* With GC-2010 Plus, approximated to a constant flow rate using a pressure program

Checking Separation of Oxygen and Nitrogen

Fig. 1 shows the chromatograms obtained by analyses after injecting 50 µL of indoor air using a gas-tight syringe. With the 1.0 m and 2.0 m columns, separation was incomplete with a resolution of under 1.5. On the other hand, the 3.0 m column achieved complete separation of oxygen and nitrogen with a resolution of over 1.5.




Fig. 1 Comparison of Separation of Oxygen and Nitrogen

Standard Gas Analysis of Inorganic Gases and Small Hydrocarbons of Low and High Concentrations

Low-concentration standard gas (1.0 mL) of inorganic gases and small hydrocarbons was analyzed by using the MGS-2030 manual gas sampler (P/N: 221-78990-41) and SPLITTER-INJ (P/N: 221-78280-41). The obtained chromatograms are indicated in Fig. 2. The analysis time with the 1.0 m column is almost half the analysis time with the 2.0 m column, showing that reduction of analysis time is possible. In addition, the rise in the baseline is small compared to that with the 2.0 m column. This contributes to noise suppression during heating and therefore enables high-sensitivity analysis of propane and propylene. On the other hand, with the 3.0 m column, propylene of 5 ppm could not be detected due to the characteristics of the packing material used for MICROPACKED-ST columns; propylene easily adsorbs to the packing material. We then analyzed high-concentration standard gas of inorganic gases and small hydrocarbons by injecting 50 µL using a gas-tight syringe. The obtained chromatograms are indicated in Fig. 3. The chromatograms indicate that by reducing the injection volume, analysis of high-concentration samples is possible.

The following are guidelines for selecting each column.

- 1.0 m : Where separation of oxygen and nitrogen is not necessary
For high-speed analysis, quantification of trace amounts of propylene
- 2.0 m : General analysis of inorganic gases and small hydrocarbons
- 3.0 m : Where separation of oxygen and nitrogen is necessary
Measurement of impurities in high-concentration substances

Analysis of inorganic gases and small hydrocarbons can be optimized by selecting a column length appropriate for the purpose of analysis.




Fig. 2 Comparison of Chromatograms of Low-concentration Standard Gas of Inorganic Gases and Small Hydrocarbons




Fig. 3 Comparison of Chromatograms of High-concentration Standard Gas of Inorganic Gases and Small Hydrocarbons

First Edition: Nov. 2017



For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Shimadzu disclaims any proprietary interest in trademarks and trade names used in this publication other than its own. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Shimadzu Corporation

www.shimadzu.com/an/

© Shimadzu Corporation, 2017

Related Products

Some products may be updated to newer models.



> Nexis™ GC-2030
Gas Chromatograph

Related Solutions

> Environment

> Price Inquiry

> Product Inquiry

> Technical Service /
Support Inquiry

> Other Inquiry