

Disulfide Bond Characterization of Monoclonal Antibody (mAb) Using Q-TOF Mass Spectrometer

Yonghai Lu and Zhaoqi Zhan

Application Development & Support Centre, Shimadzu (Asia Pacific), Singapore

Introduction

Monoclonal antibody (mAb) is emerging as the fastest-growing category of biotherapeutics with a wide range of therapeutic and diagnostic applications. Higher-order structure of mAb plays a critical role in the efficacy and safety. For example, the number of disulfide bonds and their positions are critical quality attributes (CQAs) for mAb, because incorrect disulfide linkage formation can cause a loss of biological activity or even can elicit an immune response from the host. Herein, we report a LCMS-based method to precisely characterize disulfide bonds in mAb biosimilar by a comparative analysis of non-reduced and reduced conditions. The method uses ProteaseMAX™ surfactant to denature the protein, and trypsin to digest with/without reduction and alkylation. The peptides were gradient-eluted and analyzed using a Shimadzu LCMS™-9030 Q-TOF mass spectrometer for MS scan and MS/MS analysis.

Experimental

A. Reduced condition

mAb solution: 5 mg/mL of bevacizumab biosimilar in 50 mmol/L Tris-HCl (pH 8.0) buffer.

Sample solution: the sample was prepared by taking 20 µL of mAb solution and diluted 5 times with 50 mM ammonium bicarbonate (ABC) solution.

Denaturation: 10 µL ProteaseMAX™(0.5%, w/w) was added and incubated at 60 °C to denature the mAb.

Reduction and Alkylation: 10 µL Dithiothreitol (DTT, 0.2 M) was added to reduce disulfide bonds. Alkylation was done by adding 30 µL iodoacetamide (IAM, 0.2 M).

Dilution: dilute the sample with 50 mM ABC solution to the final volume of 478 µL.

Digestion: 20 µL sequencing grade trypsin was added for protein digestion at 37°C for overnight.

Stop: 2 µL trifluoroacetic acid (TFA) was added to stop trypsin activity.

Analysis: the analytical conditions on LCMS-9030 (Q-TOF) are shown in Table 1.

B. Non-Reduced condition

The procedure of non-reduced condition is the same as it used in reduced condition except without the step of reduction and alkylation.

Table 1. Analytical conditions on LCMS-9030 (Q-TOF)

Column	: Shim-pack™ GISS-HP, 3 µm, 150 × 3.0
Mobile phase	mm : (A) 0.1% FA + 0.01% TFA in water (B) 0.1% FA + 0.01% TFA in acetonitrile
Flow rate	: 0.5 mL/min
Gradient program	: B Conc. 0% (0-2 min) → 15% (10 min) → 35% (23 min) → 45% (30 min) → 75% (35-40 min) → 0% (40.1-45 min).
Column temp.	: 40°C
Injection volume	: 20 µL
Interface	: Heated ESI (positive mode)
MS Mode	: MS scan
Interface voltage	: 4.5 kV
TOF mass range	: 100 – 2000 (m/z)
Heat block temp.	: 400°C
DL temp.	: 250°C
Interface temp.	: 300°C
Nebulizing gas	: N ₂ , 3 L/min
Drying gas	: N ₂ , 10 L/min
Heating gas	: Zero air, 10L/min




Figure 1. Disulfide bond linkage structure of bevacizumab




Figure 2. Comparative total ion chromatograms (TIC) of disulfide bond peptides of bevacizumab biosimilar under non-reduced and reduced conditions. The peak # refers to tables 3 and 4.




Figure 3. MS (a) and MS/MS (b) spectra of the non-reduced disulfide bond peptide (S-S-2) of bevacizumab biosimilar.




Table 3. Detection of non-reduced disulphide bond peptides of bevacizumab biosimilar

Peak No.	RT (min)	Peptide [AA numbers]	Peptide m/z	Adduct Ion
S-S-1	24.7	Light chain, [19-23-42] VTITCSASQDISNYLNWYQQKPGK FSGSGSGTDFLTITSLQPEDFATYYCQQYSTVPWTFGQGTK [62-88-103], Light chain	1469.6837	[M+5H]5+
S-S-2	20.4	SGTASVVCLLNFFYPRL [127-134-142], Light chain VYACEVTHQGLSPVTK [191-194-207], Light chain	889.9428	[M+4H]4+
S-S-3	19.6	LSCAASGYTFTNYGMNWVR [20-22-38], Heavy chain AEDTAVYYCAK [88-96-98], Heavy chain	1124.4967	[M+3H]3+
S-S-4	24.3	Heavy chain, [140-150-153] STSGGTAALGCLVK DYFPEPVTVSWNSGALTSGVHTFPAVLQSSGLYSLSSVTVPSSSLGTQTYICNV NHKPSNTK [154-206-216], Heavy chain	1584.3884	[M+5H]5+
S-S-5	15.2	TPEVTCVVVDVSHEDPEVK [262-267-280], Heavy chain CK [327-328], Heavy chain	777.0393	[M+3H]3+
S-S-6	17.8	NQVSLTCLVK [367-373-376], Heavy chain WQQGNVFCFSVMHEALHNHYTQK [423-431-445], Heavy chain	962.2111	[M+4H]4+
S-S-7	2.5	GEC [212-214], Light chain SCDK [225-226-228], Heavy chain	379.1272	[M+2H]2+
S-S-8	23.7	THTCPPCPAPELLGGPSVFLFPPKPK [229-232-235-254], Heavy chain THTCPPCPAPELLGGPSVFLFPPKPK [229-232-235-254], Heavy chain	1364.7007	[M+4H]4+

Table 4. Detection of reduced disulphide bond (i.e., Cys-containing) peptides of bevacizumab biosimilar

Peak No.	RT (min)	Peptide [AA numbers]	Peptide m/z	Adduct Ion
S1	20.9	R.VTITCSASQDISNYLNWYQQKPGK.A [19, 42] Light chain	934.4560	[M+3H] ³⁺
S2	25.3	R.FSGSGSGTDFLTITSLQPEDFATYYCQYSTVPWTFGQGTK.V [62, 103] Light chain	1554.0400	[M+3H] ³⁺
S3	23.3	K.SGTASVVCLLNNFYPR.E [127, 142] Light chain	899.4510	[M+2H] ²⁺
S4	13.6	K.VYACEVTHQGLSSPVTK.S [191, 207] Light chain	938.4664	[M+2H] ²⁺
S5	2.5	R.GEC.- [212, 214] Light chain	365.1119	[M+H] ⁺
S6	20.7	R.LSCAASGYTFTNYGMNWVR.Q [20, 38] Heavy chain	1099.4914	[M+2H] ²⁺
S7	12.5	R.AEDTAVYYCAK.Y [88, 98] Heavy chain	645.7863	[M+2H] ²⁺
S8	15.7	K.STSGGTAAALGCLVK.D [140, 153] Heavy chain	661.3433	[M+2H] ²⁺
S9	25.0	K.DYFPEPVTVSWNSGALTSGVHTFPAVLQSSGLYSLSSVTVPSSTLTQTYICNV NHKPSNTK.V [154, 216] Heavy chain	1679.0810	[M+4H] ⁴⁺
S10	2.1	K.SCDK.T [225, 228] Heavy chain	509.2016	[M+H] ⁺
S11	22.0	K.THTCPPCPAPELLGGPSVFLFPPKPK.D [229, 254] Heavy chain	948.8226	[M+3H] ³⁺
S12	17.3	R.TPEVTCVVVDVSHEDPEVK.F [262, 280] Heavy chain	713.6804	[M+3H] ³⁺
S13	2.0	K.CK.V [327, 328] Heavy chain	307.1429	[M+H] ⁺
S14	16.8	K.NQVSLTCLVK.G [367, 376] Heavy chain	581.3172	[M+2H] ²⁺
S15	17.2	R.WQQGNVFCFSVMHEALHNHYTQK.S [423, 445] Heavy chain	934.4253	[M+3H] ³⁺

□ Results and Discussion

A. Disulfide bond linkages in bevacizumab

Typical structure of bevacizumab is shown in **Figure 1**. It contains a total of 16 disulfide bonds, 12 of which are intra-chain linkages (four in the light chain and eight in the heavy chain), whereas the other 4 disulfide bonds are inter-chain linkages (two links light chain and heavy chain, and two links the two heavy chain in the hinge region). By employing reduction reagent DTT, 15 Cys-containing tryptic peptides could be produced from the mAb (reduced disulfide bonds). Without reduction by DTT, 8 intact disulfide bond linked peptides could be produced (non-reduced disulfide bonds).

B. Detection of disulfide bond peptides

The S-S bond of bevacizumab biosimilar were cleaved in reduced sample, while the intact S-S linked peptides remained in the non-reduced sample. To identify these disulfide bond peptides, both reduced and non-reduced samples were prepared and analyzed by LCMS-9030 (Q-TOF) under the same analytical conditions.

Figure 2 shows the total ion chromatograms (TIC) of both non-reduced and reduced samples. The mirror plot clearly indicates the major different peaks between the two samples. Compared to the reduced sample, there were several large peptide molecules with longer elution times measured in the non-reduced sample. By comparing with theoretical masses of tryptic disulfide bond peptides using the Skyline s/w, all the 8 intact S-S linked peptides (S-S-1/8) and the 15 Cys-containing peptides (S1/15) from bevacizumab biosimilar were correctly measured by LCMS-9030 (Q-TOF) with <3 ppm mass error. The peptide sequences and their accurate masses are shown in **Tables 3 and 4**.

C. MS/MS sequencing of disulfide bond peptides

De novo sequencing analysis of all the disulfide bond peptides listed in Tables 3 and 4 were performed. Here we show the results of non-reduced S-S-2 and reduced S-3 as representative examples only. **Figure 3** shows *de novo* sequencing of non-reduced disulfide bond peptide (S-S-2) from light chain of bevacizumab biosimilar through a MS/MS spectrum with the fragment ion annotation. The peptide sequence and accurate *m/z* values of multiple-charged ions (2+, 3+, 4+) are displayed in the top panel. MS/MS spectrum of the 4+ ion (*m/z* 889.9445) is shown in the bottom panel, with the explanations (y ion) of matching fragments based on the peptide sequence. In addition to the non-reduced disulfide bond peptides, we also conduct *de novo* sequencing for reduced disulfide bond peptides. **Figure 4** shows the MS/MS spectra of the reduced peptide (S3) from the S-S-2 in the reduced sample.

□ Conclusions

A straightforward LCMS-based method for accurate disulfide bond peptide characterization of mAb biosimilar was established on LCMS-9030 (Q-TOF). The MS/MS spectra with fragmentation data provide high confidence results on sequencing analysis. The demonstrated performance for bevacizumab biosimilar in detection and *de novo* sequencing of non-reduced and reduced disulfide bonds signifies its practicability for the structural characterization of mAb biosimilars.

□ Acknowledgement

The authors sincerely thank Dr Atsuhiko Toyama from Marketing Innovation Centre, Shimadzu for his valuable suggestions and discussion during this work.

LCMS and Shim-pack are trademarks of Shimadzu Corporation in Japan and/or other countries. ProteaseMAX is a trademark of Promega Corporation.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".



SHIMADZU (Asia Pacific) Pte. Ltd
79 Science Park Drive, #02-01/08 Cintech IV, Singapore 118264,
www.shimadzu.com.sg; Tel: +65-6778 6280 Fax: +65-6778 2050

For Research Use Only. Not for use in diagnostic procedures.

Contents and/or instrumentations in this application may not be available in some countries under each regulation.
Please contact the local representatives in details.

Copyright © 2019 SHIMADZU (Asia Pacific) Pte. Ltd.
All rights reserved. No part of this document may be reproduced in any form or by any means without permission in writing from SHIMADZU (Asia Pacific) Pte. Ltd.

Related Products

Some products may be updated to newer models.



> **LCMS-9030**

Quadrupole Time-of-Flight Liquid Chromatograph Mas...

HPLC Column Solutions - Reversed phase and Normal phase -

Shim-pack

> **Shim-pack GISS Series**

HPLC Column Solutions

Related Solutions

> **Biopharmaceutical**

> **Small Molecule Pharmaceutical**

> **Life Science**

> **Price Inquiry**

> **Product Inquiry**

> **Technical Service / Support Inquiry**

> **Other Inquiry**