

# Application News

## No. P106

### Electron Probe Microanalyzer

## Trace Element Mapping of Plagioclase

### ■ Introduction

The minerals contained in rocks record physical and chemical information from the time when the rocks were formed. Shimadzu EPMA is useful for investigating the precise chemical composition and 2-dimensional distribution of the elements in these rock-forming minerals. Here, element mapping of plagioclase in lava from Mt. Unzen volcano in Japan was carried out using a Shimadzu EPMA™ electron probe microanalyzer (EPMA-8050G, hereinafter, EPMA). As a result, it was possible to discover the microstructure in the rim part of the plagioclase.

The thin-slice specimens used here were provided by Mr. Hideto Yoshida (Technical Director, Technical Division, School of Science) Department of Earth and Planetary Science, Graduate School of Science, The University of Tokyo.

R. Ogawa, T. Ono

### ■ Element Mapping of Plagioclase

Plagioclase, the principal component mineral of volcanic rock, is a solid solution consisting of anorthite ( $\text{CaAl}_2\text{Si}_2\text{O}_8$ ) and albite ( $\text{NaAlSi}_3\text{O}_8$ ). In the cooling process of high-temperature magma, first, crystals with a large anorthite component form, and then, as cooling proceeds, the composition changes to one with a large albite component. Areas where the chemical composition changes discontinuously appear as stripes. This is referred to as "zonal texture."

Fig. 1 shows the result of a wide-area EPMA mapping analysis of the principal elements of plagioclase, O, Na, Al, Si, and Ca, and the trace elements Mg, K, and Fe. Concentric zoning can be seen, and as a general tendency, the core area has a large anorthite content, while the rim area has an albite-rich composition. However, the outermost rim has an anorthite-rich composition, which is termed "reverse zoning," and a reaction rim can be observed (Fig. 2). These facts suggest that the temperature, pressure, and abundance of  $\text{H}_2\text{O}$  changed during the crystal cooling process, and interaction with other magma also occurred.




Fig. 1 Wide-Area Element Mapping of Plagioclase

## ■ Mapping of Rim Area of Plagioclase

In the rim area of phenocrysts in rocks, various phenomena occur as a result of contact with other magma following once cooling. To elucidate those changes, the rim area of the plagioclase was mapped under higher magnification. As a result, it was found that counts of Na decreases in the rim area, while counts of Al, Ca, and trace elements Mg, K, and Fe increase (Fig. 2).

Further enlargement of the rim area revealed that the contents of the principal elements of plagioclase, namely, Na, Al, and Ca are small in parts where the K content is large, and in contrast, particles of Mg and Fe had collected in those parts (Fig. 3).




Fig. 2 Element Mapping of Rim Area of Plagioclase




Fig. 3 Element Mapping of Rim Area of Plagioclase (Enlarged Views)

## ■ Trace Element Mapping of Rim Area of Plagioclase

High magnification mapping revealed that Mg- and Fe-containing inclusions that concentrated in parts of the rim with high K contents existed at a number of scattered points as microscopic particles with a maximum size of several 100 nm and a minimum size of several 10 nm or less (Fig. 4).

As this experiment shows, it was possible to obtain valuable information for understanding the history of interaction between phenocrysts and magma by high resolution mapping with the Shimadzu EPMA.




Fig. 4 Mapping of Inclusions Consisting of Mg and Fe

EPMA is a trademark of Shimadzu Corporation in Japan and/or other countries.

First Edition: Feb. 2020



For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Shimadzu disclaims any proprietary interest in trademarks and trade names used in this publication other than its own. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2020

Shimadzu Corporation

[www.shimadzu.com/an/](http://www.shimadzu.com/an/)

# Related Products

Some products may be updated to newer models.



> EPMA-8050G

Electron Probe Microanalyzer

# Related Solutions

- > Geology, Mineralogy

> Chemicals

> Environment

> Soil

> Price Inquiry

> Product Inquiry

> Technical Service / Support Inquiry

> Other Inquiry