

## Application News

# No. C147A

### nSMOL™ Antibody BA Kit

## LCMS Bioanalysis of Antibody Drugs Using Fab-Selective Proteolysis nSMOL - Part 3 - Nivolumab analysis -

### ■ nSMOL™ Antibody BA Kit Features

nSMOL is Shimadzu's completely new and breakthrough technology that enables selective proteolysis of the Fab region of monoclonal antibodies. This technique facilitates method development independent of the variety of antibody drug and achieves a paradigm shift in the bioanalysis of antibody drugs.

Furthermore, this is the only method with respect to antibody drugs that has fulfilled the criteria of "Guideline on Bioanalytical Method Validation in Pharmaceutical Development" for low MW drug compounds issued by the Japanese Ministry of Health, Labour and Welfare. Shimadzu also offers optimization methods and protocols, and nSMOL can be applied to clinical research at various institutions.

### ■ Method Validation for Nivolumab Bioanalysis

Cancer cells have been found to evade immune surveillance mechanism through the expression of immunosuppressive ligands, and avoid cytotoxicity from immune cells.

Nivolumab was developed by Dr. Honjo et al. as a breakthrough medicine that activate immune cells by blocking PD-1 mediated inhibitory signals\*. Innovative drugs that apply these immunological mechanisms are named as immune checkpoint inhibitors, and many drug discovery for this field now continue to progress around the world.

These medicines are used in a cancer chemotherapy which act on advanced and complex immunological mechanisms. Therefore, it is important to progress integrative clinical trials in order to develop more efficient treatments by using many clinical indexes and biomarkers.

Shimadzu has applied the nSMOL and performed analytical validation of Nivolumab for the pharmacokinetic monitoring into early clinical implementations.

### ■ Quantitation Peptides of Nivolumab

Peptide	MRM transition	Purpose
P <sub>14</sub> R	512.1>292.3 (b3+)	For quantitation (IS)
	512.1>389.3 (b4+)	For structure confirmation
	512.1>660.4 (b6+)	For structure confirmation
ASGITFSNSG MHWVR	550.8>661.5 (y11++)	For quantitation
	550.8>746.4 (y13++)	For structure confirmation
	550.8>785.4 (y6 +)	For structure confirmation

* Quantitation range in human plasma	: 0.15 to 300 µg/ml
Averaged accuracy	: 100.4 %

### ■ Analysis Conditions for Nivolumab Using the nSMOL


#### <Sample Processing Protocol>

With the nSMOL, the same sample processing protocol can be applied to all antibody drugs. For details, refer to Shimadzu Application News (Trastuzumab analysis).

#### <LCMS Analysis Conditions>

[LC] NexeraX2 System	
Column	: Shim-pack GLSS C18 (50 mm × 2.1 mm)
Column oven	: 50 °C
Solvent A	: 0.1 % formic acid/water
Solvent B	: 0.1 % formic acid/acetonitrile
Gradient	: 1 %B (1.5 min)/1-40 %B (3 min)/95 %B (1 min)/1 %B (1 min)
Flow rate	: 0.4 mL/min
Injection	: 10 µL
[MS] LCMS-8050, 8060	
Ionization	: ESI Positive
DL	: 250 °C
Heat Block	: 400 °C
Interface	: 300 °C
Nebulizer gas	: 3 L/min
Drying gas	: 10 L/min
Heating gas	: 10 L/min


### ■ Structure Configuration of Nivolumab Candidate Signature Peptides Identified in nSMOL Reactions



**Fig. 1 Structure Configuration of Nivolumab Tryptic Peptides**

Detected peptides are indicated in red (heavy chain) and green (light chain). Fv-selective proteolysis has been progressing by nSMOL.

## ■ MRM Chromatograms



**Fig. 2** MRM Chromatograms of ASGITFSNSGMHWVR (Blue), and P<sub>14</sub>R Internal Standard (Black) (in Human Plasma)

## ■ Full Validation Results for Nivolumab

## &lt;Precision and accuracy&gt;

Set Concentration [µg/ml]	Data Average (N = 15)	Accuracy (%)	CV (%)
2.93	2.97	101	7.51
200	202	101	6.75

## &lt;Freeze-thaw test&gt;


Set Concentration [µg/ml]	Data Average (N = 5)	Accuracy (%)	Temperature (°C)
2.93	2.73	95.6	-20
200	183	96.1	-20

## &lt;Long-term stability test&gt;

Set Concentration [µg/ml]	Data Average (N = 5)	Accuracy (%)	Temperature (°C)
2.93	3.03	104	-20
200	213	107	-20

## &lt;Processed sample stability for 48 h&gt;

Set Concentration [µg/ml]	Data Average (N = 5)	Accuracy (%)	Temperature (°C)
2.93	3.08	105	5
200	195	97.6	5



**Fig. 3** Nivolumab Calibration Curve

## ■ Observations, Conclusions, and References

Although eight candidate signature peptides including CDRs were obtained using nSMOL, only the peptide ASGITFSNSGMHWVR indicated a positive correlation to drug concentration. This indicates that the sequence homology of fully human antibodies and endogenous IgGs is extremely similar.

In order to set suitable bioanalysis conditions, peptide candidates with structural specificity must be strictly selected. By utilizing Fv-selective reactions, Shimadzu nSMOL greatly facilitates the development of assay methods. The lower limit of quantitation is 0.15 µg/ml and the same assay method can be used from preclinical to clinical trials.

## &lt;References&gt;

- \* Ishida Y, Agata Y, Shibahara K, and Honjo T., *EMBO J*, 1992, 11(11):3887  
 Iwamoto N et al. *Analyst*, 2014, DOI:10.1039/c3an02104a  
 Iwamoto N et al., *J. Chromatogr. B*, 2016, DOI:10.1016/j.jchromb.2016.04.038

## &lt;Chief Scientists&gt;

Noriko Iwamoto, Ph.D. and Takashi Shimada, Ph.D., Technology Research Laboratory

Notes: The product described in this document has not been approved or certified as a medical device under the Pharmaceutical and Medical Device Act of Japan.

It cannot be used for the purpose of medical examination, treatment or related procedures.

First Edition: Mar. 2017



**For Research Use Only. Not for use in diagnostic procedures.**

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2017

Shimadzu Corporation

www.shimadzu.com/an/