

Application News

No. i217

Dynamic Ultra Micro Hardness Tester

Hardness Test for Zirconia Particles

Introduction

Fracture strength in a compression test is generally used in strength evaluations of fine particles such as various types of industrial raw materials. Recently, however, there have also been requests for evaluation techniques of the hardness of particles.

As a test method, specimens for hardness evaluations are prepared by embedding the fine particles which are the object of measurement in resin and then polishing the specimen surface. This article introduces an example in which the hardness of zirconia powder was measured using a Shimadzu DUH™-211S dynamic ultra micro hardness tester.

Zirconia powder is now used in a diverse range of products, beginning with materials for dental treatment, taking advantage of its characteristic of high hardness.

Testing Apparatus and Specimens

The testing apparatus used in this experiment was a Shimadzu DUH-211S dynamic ultra micro hardness tester (Fig. 1 shows the appearance of a representative system). Two types of zirconia powder with different particle diameters (100 μm and 30 μm) were used in the evaluation. Samples were prepared by fixing the particles in an embedding resin and polishing the specimen surface. Table 1 shows the details of the samples used in the test.

Table 1 Test Samples

1) Sample name	Zirconia powder	
2) Sample symbol	A	B
3) Particle diameter (μm)	100	30
4) Sample dimensions (mm)	φ25 × 10 t (embedding resin)	

Test Conditions

Table 2 shows the test conditions (loading jig and loading conditions). Loading was applied with a diamond triangular pyramid indenter having a tip angle of 115 degrees. Fig. 2 shows an image of the loading method.

Table 2 Test Conditions

1) Testing machine	DUH-211S dynamic ultra micro hardness tester (see Fig. 1)
2) Indenter	The diamond triangular pyramid indenter with tip angle 115 degrees
3) Type of test	Load-unload test
4) Test force (mN)	49
5) Loading rate (mN/s)	2.685
6) Holding time (s)	5




Fig. 2 Loading Image



Fig. 1 Overview of DUH-211S Dynamic Ultra Micro Hardness Tester

Test Results

Fig. 3 shows the results of the hardness measurements by a load-unload test conducted under the conditions shown on the previous section (average of the measured values at multiple points for both samples A and B). Fig. 4 shows the test force-indentation depth curves for representative measurement points.

Hardness test results (average values)					
Sample name	Sample symbol	Fmax [mN]	hmax [μm]	HMT115 [N/mm ²]	HV*
Zirconia	A	49.37	0.543	5355.1	1155.4
	B	49.40	0.742	2914.6	864.5

Note 1) The symbols shown in the above table have the following meanings.
 Fmax: Maximum test force
 hmax: Maximum indentation depth
 HMT115: Martens hardness obtained with a 115 degrees triangular pyramid indenter
 $HMT115 = F_{max} / (26.43 \times h_{max}^2)$
 HMT115 is value which includes machine compliance correction (Cf correction) by the rigidity of the testing machine and indenter area function correction (Ap correction) for the influence of the roundness of the indenter tip.
 HV*: Converted Vickers hardness (reference value)

From these results, it can be understood that sample A has a higher Martens hardness (HMT115) and a higher converted Vickers hardness (HV*) than sample B.




Fig. 3 Test Results (Measured Values)






Fig. 4 Test Results (Test Force-Indentation Depth Curves)

Fig. 5 shows enlarged images of the indentations of the two samples after unloading, as photographed with a microscope. If a hardness test is conducted with a sharp-tipped indenter under a condition in which fine particles are loosely distributed and not fixed, accurate measurement will be impossible due to scattering of the particles. However, as described in this example, it is possible to prepare samples

with a polished fine and flat surface, even when testing fine particles, if the particles are embedded in resin and polished. Because the DUH-211S dynamic ultra micro hardness tester is capable of measuring hardness by forming indentations with extremely low force while aiming at the center of the particles, it is an effective tool for testing fine particles, as shown in this experiment.



Sample A (particle diameter: 100 μm)



Sample B (particle diameter: 30 μm)

Fig. 5 Results of Microscopic Observation of Test Specimens

DUH is a trademark of Shimadzu Corporation in Japan and/or other countries.

First Edition: Feb. 2021



For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Shimadzu disclaims any proprietary interest in trademarks and trade names used in this publication other than its own. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Shimadzu Corporation

www.shimadzu.com/an/