

Application News

Gas Chromatograph Mass Spectrometer GCMS-TQ™8050 NX

A Sample Prep-Free Analysis of Triglycerides and Fatty Acids with “Smart IS+” and “SMCI+” - Authentication Studies -

Chua Chun Kiang¹, Jackie¹, Riki Kitano²

¹Shimadzu (Asia Pacific) Pte Ltd., ²Shimadzu Corporation Japan

User Benefits

- ◆ “Smart IS+” and “SMCI+” enable derivatization-free procedures for qualitative GC/MS analysis of fatty acids and triglycerides
- ◆ Easily switch between both electron ionization and positive chemical ionization modes with “Smart IS+” setup and “SMCI+” with methanol reagent gas produces characteristic mass spectra for saturated fatty acids
- ◆ “Smart IS+” and “SMCI+” support quick preliminary authentication study of vegetable oil, e.g., coconut oil and sunflower oil

■ Introduction

Triglycerides are the primary constituents of body fat in animals as well as vegetable fat. Triglycerides are tri-esters that constitute glycerol and three fatty acids. The analysis of triglycerides and fatty acids is crucial to refining the current physiological understanding of animals and plants.

The direct analysis of free fatty acids (FFAs) and triglycerides (TGs) with gas chromatography/mass spectrometry (GC/MS) technique requires specialized columns. Otherwise, multi-step transesterification or derivatization procedures are necessary to enable GC/MS analysis. The usage of a direct probe as a sample inlet in GC/MS provides an alternative technique that eliminates the need for tedious sample preparation. This benefit paves the way for a quick and direct method to achieve preliminary detection and identification of FFAs and TGs.

This article demonstrates the application of a direct probe in conjunction with a Smart EI/CI ion source (Smart IS) or solvent mediated chemical ionization (SMCI) unit for the analysis of FFAs and TGs. The mass spectra of FFAs and TGs generated by Smart IS and SMCI units will be evaluated, followed by examples on their usage for preliminary authentication study of vegetable oils.

■ Measurement Conditions and Samples

Analytical Setup.

The analytical results in this report were generated using a Direct Sample Inlet (DI) probe in conjunction with a Smart IS or SMCI unit. The combination of DI with Smart IS or SMCI unit is hence known as Smart IS and “SMCI+” in this article, respectively (Fig. 1).

The DI probe is designed to be able to fit a miniature sample vial at its tip. The sample vial is thereafter placed close to the ion source and subsequently heated up according to a temperature program. The chemicals in the sample vial are hence volatilized and ionized in the ion source.

Smart IS is a 2-in-1 ion source that enables both electron ionization (EI) and positive chemical ionization (PCI) modes. PCI is achieved with the usage of isobutane gas as a reagent gas. Due to the simplicity of switching between two different ionization modes with Smart IS, the PCI mode attained with Smart IS is referred to as quick chemical ionization (QCI).

On the other hand, the SMCI unit enables PCI mode with conventional PCI ion source and methanol as the reagent gas. Usage of methanol allows safe (i.e., it eliminates the use of flammable and toxic reagent gases such as methane, isobutane,



Fig. 1 Polymode Ionization setup inclusive of “Smart IS+” and “SMCI+”.

Name	Abbreviation	Formula	MW
Fatty Acids			
Decanoic Acid	C10:0	C ₁₀ H ₂₀ O ₂	172
Lauric Acid	C12:0	C ₁₂ H ₂₄ O ₂	200
Myristic Acid	C14:0	C ₁₄ H ₂₈ O ₂	228
Palmitic Acid	C16:0	C ₁₆ H ₃₂ O ₂	256
Stearic Acid	C18:0	C ₁₈ H ₃₆ O ₂	284
Oleic Acid	C18:1	C ₁₈ H ₃₄ O ₂	282
Triglycerides			
Trilaurin	LLL	C ₃₉ H ₇₄ O ₆	639
Trimyristin	MMM	C ₄₅ H ₈₆ O ₆	723
Tripalmitin	PPP	C ₅₁ H ₉₈ O ₆	807
Tristearin	SSS	C ₅₇ H ₁₁₀ O ₆	891
Triolein	OOO	C ₅₇ H ₁₀₄ O ₆	885

and ammonia) and convenient adoption of PCI mode in routine GC/MS analysis.

A total of 6 free fatty acids (i.e., decanoic acid C10:0, lauric acid C12:0, myristic acid C14:0, palmitic acid C16:0, stearic acid C18:0, and oleic acid C18:1) and 5 triglycerides (i.e., trilaurin LLL, trimyristin MMM, tripalmitin PPP, tristearin SSS, and triolein OOO) were analyzed in this study.




Fig. 2. EI, QCI, and SMCI mass spectra of FFAs collected with "Smart IS+" and "SMCI+". SI: Library similarity index.




Fig. 3. EIT profiles of QCI mode for a mixture of FFAs.

Experimental Condition.

Standard solutions of FFAs and TGs were prepared to a concentration of 5000 ppm in ethanol and acetone, respectively. 1 μ L of the standard solution was introduced into the DI sample vial for analysis. A mixture sample of FFAs and TGs were prepared by introducing 1 μ L of each standard solution into the DI sample vial. The samples were left to dry before analysis.

The DI probe was heated at 20 $^{\circ}$ C/min to 100 $^{\circ}$ C, then 40 $^{\circ}$ C/min to 450 $^{\circ}$ C and held for 7 min. The ion source temperature was set to 230 $^{\circ}$ C. Ionization mode used included EI, QCI (isobutane), and SMCI (methanol). Scan mode was performed in the range of m/z 50-900 with a scan speed of 3333.

Cold-pressed coconut oil and regular sunflower oil purchased from a local supermarket were used for the analysis.

Results and Discussion

“Smart IS+” and “SMCI+” Mass Spectra of FFAs.

Using “Smart IS+” and “SMCI+”, the FFAs were first analyzed with three different ionization modes, EI, QCI, and SMCI, to establish the mass spectra. The consolidated mass spectra of the FFAs collected with the three ionization modes are shown in Fig. 2.


The EI mass spectra of the FFAs were successfully matched to the NIST mass spectral library with high similarity index scores of >85. Despite that, the intensities of molecular ion peaks for C10:0 and C18:1 were rather low. The base peaks in the mass spectra were small fragmentation products.

In the QCI mass spectra of saturated FFAs, the protonated molecule, $[M+H]^+$ ion, was detected as the base peak. The mass spectrum of unsaturated FFA C18:1 indicated the $[M-OH]^+$ ion as the base peak at m/z 265. These characteristic ions could potentially be exploited to enable the rapid identification of FFAs.

Subsequent analysis with SMCI mode revealed the presence of three characteristic ions in the mass spectra of saturated FFAs. Peaks corresponding to $[M-OH]^+$, $[M+H]^+$, and $[M-OH+OMe+H]^+$ ions were found in the mass spectra. The presence of $[M-OH+OMe+H]^+$ ion peak suggested a possible reaction between the FFAs and methanol reagent gas, resulting in the formation of fatty acid methyl ester (FAME).

In contrast, the esterification product was not observed for C18:1. All in all, the characteristic ions found in QCI and SMCI mass spectra could also serve as confirmatory signals for the rapid identification of FFAs. Furthermore, the concurrent presence of $[M-OH]^+$, $[M+H]^+$, and $[M-OH+OMe+H]^+$ ions in SMCI mass spectra could further serve as an indicator for the presence of saturated FFAs.

A mixture of FFAs was subsequently analyzed to evaluate their elution profile. The extracted ion thermogram (EIT) profiles of $[M+H]^+$ ion for the FFAs collected with QCI mode are shown in Fig. 3. The FFAs eluted within 225 $^{\circ}$ C with a visible resolution, starting from shorter carbon chain species. The order of elution of the saturated FFAs was from C10:0, C12:0, C14:0, C16:0 to C18:0. C18:1 eluted slightly earlier than C18:0, possibly due to lower volatilization temperature as a result of unsaturation in the carbon chain.



(continued on next page)




Fig. 4. EI, QCI, and SMCI mass spectra of TGs collected with "Smart IS+" and "SMCI+". The formula of FA refers to $C_{12}H_{24}O_2$ for LLL, $C_{14}H_{28}O_2$ for MMM, $C_{16}H_{32}O_2$ for PPP, $C_{18}H_{36}O_2$ for SSS, and $C_{18}H_{34}O_2$ for OOO.

"Smart IS+" and "SMCI+" Mass Spectra of TGs.

The consolidated mass spectra of the TGs collected with the three ionization modes are shown in Fig. 4. The EI mass spectra of the TGs contained distinct peaks corresponding $[M-FA]^+$ ion peak resulting from the α -cleavage of a fatty acid chain from the triglyceride. In addition, each triglyceride has its characteristic ions, such as m/z 183 and 257 for LLL, m/z 211 and 285 for MMM, m/z 239 and 313 for PPP, m/z 267 and 341 for SSS, and m/z 265 and 339 for OOO. These characteristic peaks were also observed in the QCI and SMCI mass spectra. The $[M-FA]^+$ ion peaks in QCI mass spectra were detected as base peak, which was also the case for SMCI mass spectra, except for SSS.

A mixture of TGs was then analyzed to evaluate their elution profile. The extracted ion thermogram profiles of $[M-FA]^+$ ion for the TGs collected with QCI mode are shown in Fig. 5. Like the elution order of FFAs, the TGs eluted with a visible resolution, starting from the shorter carbon chain species within the temperature range of 150 to 350 °C. The order of elution was LLL, MMM, PPP, OOO, and SSS. OOO eluted at an almost similar temperature as SSS.

Authentication of Vegetable Oils.

By using the EIT of the unique $[M-FA]^+$ ion of each TGs on GC/MS data of vegetable oil, it is possible to determine the authenticity of the oil. Fig. 6 shows the TIT profiles of coconut and sunflower oil typically used for cooking. Coconut oil is known to contain a high amount of saturated fatty acids while sunflower oil is known for its high content of unsaturated fatty acids. The fatty acids are usually present in the form of triglycerides. The most prevalent saturated fatty acid in coconut oil is C12:0, which correlated to LLL. On the other hand, the most prevalent unsaturated fatty acid in sunflower oil is C18:1, which correlated to OOO.

As a result of analyzing coconut oil with SMCI mode (Fig. 6, top), distinctive $[M-FA]^+$ ion peaks of LLL and MMM were observed while only a slight amount of PPP was present. As expected, unsaturated triglycerides, i.e., OOO, was barely observed. A closer look at the mass spectrum at the peak apex of the TIT profile, within the range of m/z 300-600, indicated the $[M-FA]^+$ ions of LLL and MMM. The remaining peaks, i.e., m/z 383, 411, 467, and 523 could arise from isomers of triglycerides.




Fig. 5. EIT profiles of SMCI mode for a mixture of TGs.



Coconut Oil

Rich in saturated FAs



Sunflower Oil

Rich in unsaturated FAs




Fig. 6. TIT and EIT profiles of SMCI mode for (top) coconut oil and (bottom) sunflower oil. Inset shows the mass spectrum at peak apex of the corresponding TIT profile.

Subsequent analysis of sunflower oil with SMCI mode (Fig. 6 bottom) indicated the distinctive presence of $[M-FA]^+$ ion for OOO and a trace amount of PPP. LLL and MMM species were not observed. Similarly, a closer look at the mass spectrum at the peak apex of the TIT profile revealed the presence of peaks such as m/z 599, 601, and 604, which could arise from isomers of triglycerides.

Conclusion

The "Smart IS+" and newly introduced "SMCI+" enable a direct and quick qualitative analysis of fatty acids and triglycerides, which conventionally require tedious derivatization steps prior to GC/MS analysis. The "Smart IS+" setup delivers convenience in switching between electron ionization and positive chemical ionization mode of analysis. On the other hand, the "SMCI+" setup delivers utmost convenience and safety to carry out positive chemical ionization since it utilizes methanol, which is a common laboratory solvent, as the reagent gas. The characteristic ion in the mass spectra of triglycerides further enabled quick preliminary authentication study of coconut oil and sunflower oil, supported by the "Smart IS+" and "SMCI+" setups.

GCMS-TQ is a trademark of Shimadzu Corporation in Japan and/or other countries.



Shimadzu Corporation
www.shimadzu.com/an/

SHIMADZU (Asia Pacific) Pte. Ltd,
www.shimadzu.com.sg

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

04-JMST-202-EN

First Edition: Mar. 2021

Related Products

Some products may be updated to newer models.



> GCMS-TQ™8050
NX
Triple Quadrupole Gas Chromatograph
Mass Spectrometry



> DI-2010
Direct Sample Inlet Unit for GCMS-TQ
Series, GCMS-...



> GC/MS CI , CI/NCI
and SMC
Measurement

Related Solutions

> Small Molecule
Pharmaceutical

> Food and Beverages

> Price Inquiry

> Product Inquiry

> Technical Service /
Support Inquiry

> Other Inquiry