

## Application News

SPM-Nanoa™ Scanning Probe Microscope (Atomic Force Microscope)

### The Integrated High Performance Optical Microscope Realizes Smooth Location on SPM

K. Kuroda, S. Moriguchi

#### User Benefits

- ◆ The high performance optical microscope integrated into the SPM enables precise target location.
- ◆ Accurate target capturing eliminates the need for rework in locating the field of view, realizing more efficient operations.
- ◆ Patterned Si substrates can be easily observed without sample preparation.

#### ■ Introduction

Scanning Probe Microscope [SPM (AFM)] is a microscope that enables easy, high-resolution 3D observation at the nano meter scale in air. The range of application is extensive, from hard materials such as semiconductor substrates and metals to soft materials such as biological samples. As the miniaturization of electronic devices and bio-devices progresses, further miniaturization is required for the parts comprising these devices. For example, in the semiconductor industry, the line width of circuits used to be 10 μm when integrated circuits were first developed, but now it is 14 nm. Due to such micronization of parts, the selection of locations to observe with microscopes is now required to be even more precise compared to the past. Since SPM has a high resolution and high magnification, it is necessary to strictly specify the observed location. Although an observed location for SPM is usually selected using an optical microscope, precise positioning has not conventionally been possible, due to the lack of resolution on the optical microscope. However, SPM-Nanoa integrates a high-performance optical microscope, thereby facilitating precise location of the target on the micrometer.

#### ■ SPM-Nanoa

SPM is a microscope that observes and measures the 3D topography and local physical properties of a sample at high magnification by scanning the sample surface with a tiny probe (cantilever). Equipped with an advanced high-sensitivity detection system and automatic observation function as standard, SPM-Nanoa is a new SPM which makes it easier, more detailed, and faster to "observe" what you want to observe. It powerfully assists you in diverse tasks, from observation of the topography of microscopic areas to the measurement of physical properties. Fig.1 shows the external view of SPM-Nanoa. The three features of SPM-Nanoa are as follows:

- ① Automatic observation:  
Adjusts the optical axis of the laser, sets the conditions during observation, and automates image processing.
- ② Extensive functionality:  
Sharp and clear image both on optical microscope and SPM.
- ③ Saves time:  
Various support functionality achieves fast observation.

In this article, among ② high performance features, we introduce a case where the observation location was selected on the micrometer using the high performance optical microscope, and the microstructure was smoothly observed by SPM.




Fig. 1 SPM-Nanoa™ Scanning Probe Microscope

#### ■ Observation Sample and Conditions

In this article, a patterned Si substrate was observed. Fig. 2 shows the image observed by the optical microscope. An image of a probe is drawn in the center of the substrate, and a portion of the image was observed with SPM-Nanoa. It was possible to make a smooth transition from wide range observation with the integrated high-performance optical microscope to high-resolution observation with SPM, without losing sight of the target.




Fig. 2 Patterned Si Substrate

Table 1 shows the observation conditions. We used the optical microscope to zoom into the image from 600 μm × 600 μm to 200 μm × 200 μm, and then used SPM to further zoom in to 20 μm × 20 μm and 1 μm × 1 μm.

Table 1 Observation Conditions

Instrument	: SPM-Nanoa Scanning Probe Microscope
Scanner	: Wide range scanner (125 μm)
Observation mode	: Contact mode
Observation field	: Optical microscope 600 μm × 600 μm
	200 μm × 200 μm
	SPM 20 μm × 20 μm
	1 μm × 1 μm

## ■ Observation Results

Fig. 3 shows the observation results. Figs. 3(a) and (b) show the images observed by the optical microscope. From these observed images, it can be seen that the probe design appearing in the overall substrate image shown in Fig. 2 was composed of various types of pictograms.

Figs. 3(c) and (d) show the surface topography images observed by SPM, and Fig. 3(e) shows the cross-sectional profile. Figs. 3(c) and (e) show that the step height forming a pictogram is about 200 nm. From Figs. 3(d) and (e), it can also be seen that these pictograms are made up of etched holes with a diameter of 218.3 nm and a depth of about 24.2 nm.

## ■ Conclusion

With its integrated high performance optical microscope, the SPM-Nanoa enables precise location of targets on the order of micrometers. This feature realizes efficient high-resolution observation in SPM.

Sample provided by NanoSensors




Fig. 3 Surface Topography Images of Patterned Si Substrate Observed by Optical Microscope and SPM

SPM-Nanoa is a trademark of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.



Shimadzu Corporation  
Analytical & Measuring Instruments Division  
Global Application Development Center

[www.shimadzu.com/an/](http://www.shimadzu.com/an/)

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

01-00133-EN

First Edition: Sep. 2021

➤ Please fill out the survey

## Related Products

Some products may be updated to newer models.



### ➤ SPM-Nano

Scanning Probe Microscope/Atomic  
Force Microscope

## Related Solutions

### ➤ Electronics

### ➤ Price Inquiry

### ➤ Product Inquiry

### ➤ Technical Service / Support Inquiry

### ➤ Other Inquiry