

Application News

Fatigue and Endurance Testing Machine EMT-1kNV-50/30

Improvement of Test Control of Fatigue Testing Machine Subject to Horizontal Load using a Gas Bearing

S. Bressan

User Benefits

- ◆ Improvement of test control when the testing machine is subject to horizontal load.
- ◆ Tests according to the ASTM F384 standard can be performed with accurate control.

■ Introduction

In fatigue testing, samples are cyclically loaded with different patterns, directions, frequencies, and magnitudes to simulate the in-use conditions of mechanical components. In some cases, a load perpendicular to the actuator's movement direction (horizontal force) is applied, depending on the test. In such conditions, frictional forces occur in the actuator's internal bearing, and achieving smooth test control becomes difficult.

A newly developed gas bearing has been employed instead of the traditional standard bearing, resulting in smoother test control. In this article, an electromagnetic fatigue testing machine was loaded with a horizontal load, and the output waves obtained using a standard bearing and a gas bearing were compared to evaluate the control accuracy.

■ Testing Equipment

Tests were conducted with an electromagnetic-type fatigue testing machine. These testing machines are equipped with an electromagnetic actuator having an extremely high frequency response, allowing for highly precise fatigue tests (the controller was a Controller 4830). The testing machine was equipped with a standard (oil-less) and a gas bearing. In gas bearings, pressurized air creates a thin gas layer between the actuator and the bearing, reducing the friction caused by lateral forces and improving the test control accuracy.

■ Experimental Procedure

Tests were conducted using the conditions shown in Table 1 to evaluate the effect of the gas bearing on the control accuracy. As shown in Fig. 3, a 100 N horizontal force was applied to the actuator with a weight and pulley. The tests were stroke-controlled with a sinusoidal load path. The details of the tests are shown in Table 1, and the test results are shown in Fig. 4. Although the maximum and minimum set values were applied correctly, the actuator return movement was delayed when using the standard bearing. In the case of the gas bearing loaded with pressurized air at 0.5 MPa, the movement of the actuator appears smoother. When the air pressure was increased to 0.65 MPa, the stroke signal assumed a shape close to an ideal sine wave.



Fig. 1 Fatigue Testing Machine




Fig. 2 Gas Bearing

Table 1 Testing Conditions for the Gas Bearing Evaluation

Waveform	Sine wave
Control	Stroke
Max. Stroke	1 mm
Min. Stroke	-1 mm
Test Frequency	1 Hz
Horizontal Load	100 N
Bearing	Standard Bearing Gas bearing (0.5 MPa, 0.65 MPa)

■ ASTM F384 Compression Bending Test

Tests were conducted using jigs for orthopedic fracture fixation devices according to the ASTM F384 standard. A photograph of the apparatus employed is shown in Fig. 5. In the ASTM F384 test, the specimen is tilted, and the actuator is subjected to both vertical and horizontal components of the force. A test with a maximum force of 50 N and a stress ratio of 0.1 was conducted. The details of the tests are shown in Table 2. The waveform obtained using the gas bearing loaded with 0.65 MPa is the closest to an ideal sinusoidal waveform, similar to the previous test data shown in Fig. 4.




Fig. 3 Application of Horizontal Load to the Testing Machine




Fig. 4 Test Results

■ Conclusions

These tests compared the output wave (stroke and force) of a testing machine subject to horizontal load when using a standard and an air bearing. In the case of the standard and the gas bearing loaded with pressurized air at 0.5 MPa, the actuator movement was not smooth due to the friction caused by the lateral force. However, when the air pressure was raised to 0.65 MPa, the output wave resulted in an ideal sinewave. Therefore, when conducting tests where horizontal forces occur (such as the ASTM F384), using a gas bearing allows accurately controlled tests to be conducted.




Fig. 5 ASTM F384 Bending Compression Test

Table 2 Test Conditions for ASTM F384 Bending Compression Test

Waveform	Sine wave
Control	Force
Maximum Load	50 N
Stress Ratio	0.1
Test Frequency	1 Hz
Bearing	Standard bearing Gas bearing (0.65 MPa)




Fig. 6 Test Results



Shimadzu Corporation
Analytical & Measuring Instruments Division
Global Application Development Center

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

01-00292-EN

First Edition: Mar. 2022

➤ Please fill out the survey

Related Products

Some products may be updated to newer models.



➤ EMT Series

Electromagnetic Force Fatigue and
Endurance Testin...

Related Solutions

➤ Automotive

➤ Price Inquiry

➤ Product Inquiry

➤ Technical Service /
Support Inquiry

➤ Other Inquiry