

## Using the DPiMS-8060 Mass Spectrometer to Analyze Drugs in Plasma (1) -A Quantitative Analysis of Everolimus-

Kana Waki, Hidekazu Saiki

### User Benefits

- ◆ The concentration dependence of drugs in plasma is easily analyzed with this method.
- ◆ The analytical results obtained are unaffected by column condition or degradation.

### Introduction

Measuring the concentration of drugs in blood or plasma is a type of analysis frequently performed in the context of research and clinical sampling. Accordingly, techniques are needed for obtaining results quickly and easily. Currently, LC/MS systems are typically used to measure drug concentrations in samples. However, because of the columns used in LC/MS systems, matrix components originating from metabolic products, such as proteins in blood and plasma, must be carefully removed. If the samples are not pretreated in this way, they can cause degradation or changes to the condition of columns, which can affect the analytical results. Furthermore, inadequate pretreatment of samples can lead to instrument contamination, increasing the frequency of maintenance required.

This article introduces an analysis method using the LCMS-8045 system equipped with the DPiMS-8060 probe electrospray ionization (PESI) unit. The concentration of everolimus in commercial plasma is then measured with only simple deproteinization pretreatment.

### Principles behind the PESI method and Applicable samples

When an LC/MS system is used to analyze drugs in plasma, the deproteinization process must be performed very carefully to avoid degradation of flow lines including injector and the column. Therefore, the process from pretreatment to analysis takes approximately 30 minutes per sample, although the time can vary depending on the mobile phase and analytical method used. The PESI method used here is a direct ionization technique that does not involve the use of a column. This eliminates any impact on the analytical results due to column degradation. Furthermore, only a few dozen pL of solution are analyzed per analysis, so this method minimizes instrument contamination and the impact of matrix effects that could interfere with ionization. Given its advantages, the PESI method is ideal for simple measurements such as measuring drug concentrations in samples.




Fig. 1 DPiMS™-8060 System and the Principles behind the PESI Method




Fig. 2 Sample Preparation Process




Fig. 3 Everolimus Structure

## Measuring the Everolimus Standard Sample

One of the molecular target drug everolimus was used in this experiment (Fig. 3). Although plasma contains a large number of matrix components that interfere with ionization, it was confirmed that the concentration in plasma could be quantified using the PESI method under these conditions.

Standard samples of everolimus dissolved in a 50 % ethanol-water solution (v/v) were analyzed in Q3 Scan mode with the instrument parameter settings in Tables 1 and 2. The results showed that everolimus was detected primarily as three types of adduct ions (Fig. 4 A). A product ion scan was then performed. The Na adduct ion ( $m/z$  980.8), which had the highest intensity of the three types of ions detected by Q3, was set as the precursor ion. The analysis results showed that specific fragment ions were generated when a CE value of -51.0 was specified (Fig. 4 B). Of these, the fragment ion with the highest signal intensity ( $m/z$  389.3) was used for measuring everolimus in plasma.

## Measuring Everolimus in Plasma

Samples were prepared as per the sequence in Fig. 2. Plasma samples with various final concentrations of everolimus were prepared. Next, equal amounts of ethanol were added, and the solutions were mixed in a vortex mixer. The mixed solutions were centrifuged, and then collected supernatant for use as analytical samples. In terms of the analytical conditions, the values configured are indicated in Tables 1 and 2. The calibration curve prepared from the MRM mode sample analysis results is shown in Fig. 5.

The calibration curve for everolimus in plasma was good linearity over the concentration ranges from 3 to 50 ng/mL, with  $R^2 = 0.9972$  coefficient of determination.

Table 1 PESI Actuation Conditions

Ionization Position	: -37 mm
Ionization Stop Time	: 200 msec
Sampling Position	: -45.0 mm
Sampling Stop Time	: 50 msec
Probe Speed	: 250 mm/s
Probe Acceleration	: 0.63 G

Table 2 Mass Spectrometer Analytical Conditions

DL Temperature	: 250 °C
Heat Block Temperature	: 50 °C
Interface Voltage	: + 2.45 kV (ESI-positive mode)
Sample Acquisition Stop time	: 50 msec
Scan Range	: $m/z$ 10 – 1200 (Q3 Scan)
	: $m/z$ 10 – 500 (Product ion scan)
Scan Speed	: 5000 u/sec





Fig. 4 Everolimus Standard Sample Measurement Results (5 ppm)



Everolimus Concentration (ng/mL)	Mean Peak Area	SD	%RSD
3	1712	162.3	9.5
5	3635	467.2	12.8
10	5857	943.2	16.1
20	10756	1158.6	10.8
30	20487	2898.4	14.1
50	34640	4256.4	12.3
100	66089	13650.1	20.6

Fig. 5 Results for the Measurement of Everolimus in Plasma

A: Calibration Curve for Everolimus Concentration in Plasma (n = 3);

B: %RSD Value Calculated from Mean Area and Standard Deviation (SD) Values for Each Concentration

DPiMS is a trademark of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.



Shimadzu Corporation

www.shimadzu.com/an/

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu.

See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.