

# Application News

High Performance Liquid Chromatograph Mass Spectrometer LCMS-9050

# Rapid Screening Analysis of On-Line Purified Residual Pesticides in Crop Extract Using Accurate Mass Information

Junna Nakazono and Tetsuo lida

## **User Benefits**

- Enables analysis using positive/negative ion polarity switching by the LCMS-9050.
- Enables comprehensive analysis of compounds using the accurate mass and retention time information.
- This method using an online purification column enables quick and easy sample preparation.

# Introduction

The pesticide residue limits in food are strictly regulated by each region or country in order to protect the people's health. In recent years, the number of pesticides regulated around the world is on the rise.

The QuEChERS method published by the USDA (United States Department of Agriculture) in 2003, which has been adopted as an official method by organizations such as the Association of Analytical Communities (AOAC) and Committee of European Normalization (CEN), enables efficient extraction of pesticides without the need for special equipment. However, there are problems with reproducibility being dependent on operators and the time required for the work.

In this article, an example of comprehensive analysis of residual pesticides in spinach using a Revive In-Line Sample Preparation (ILSP) column<sup>1)</sup> (RESTEK Co.) and the quadrupole time-of-flight mass spectrometer LCMS-9050 (Fig. 1) is introduced. The Revive ILSP column, which enables on-line cleanup of compounds derived from the matrix, makes the sample preparation quicker and easier. It is also expected that this method can contribute to the reduction of waste and costs associated with pretreatment.



Fig. 1 Exterior of LCMS-9050

# Sample Preparation

Commercially available spinach and pesticides standard solution 74, 75 (Kanto Chemical Co., Inc.) were used for this analysis. The detailed preparation processes are shown in Fig. 2. 1.0 g of spinach that was frozen with dry ice and crushed was put in a 15 mL tube, and 3 mL of 1 % acetic acid in acetonitrile was added, then the tube was shaken. Subsequently, the tube was centrifuged after waiting for 2 minutes, and the supernatant was collected for sampling. Compared with the QuEChERS method, this method enables quicker and easier sample preparation. For example, it is expected that the time needed for sample preparation can be reduced to approximately one third by using this method when fourteen samples of spinach are pretreated<sup>1</sup>).

The recovery rate for losses in the preparation process and matrix effects were evaluated by adding a fixed concentration of pesticides standard solution to the spinach before extraction.



Fig. 2 Workflow for Sample Preparation (Left: Revive ILSP, Right: QuEChERS)

# Analytical Conditions

For the analysis of pesticides, the method included in the LC/MS/MS Method Package Residual Pesticides Ver. 3 was applied to the LCMS-9050. The HPLC and MS conditions are shown in Table 1, and the outline of this system and the LC time program are shown in Figs. 3 and 4 respectively.

Hydrophobic compounds from the matrix in injected samples are separated from pesticides by the Revive ILSP column, and the pesticides are eluted to the analytical column (upper of Fig. 3). On the other hand, the compounds from the matrix are backflushed by a pump for washing and eluted to a drain, after switching a 6-port valve when the targeted pesticides have been eluted from the Revive ILSP column (bottom of Fig. 3 and Fig. 4). In this way, the Revive ILSP column can be reused.

Tab	le 1 Analytical Conditions
UHPLC (Nexera <sup>™</sup> X3 system	n)
Analytical Column:	Shim-pack <sup>™</sup> Velox Biphenyl (100 mm L <sup>×</sup> 2.1 mm l.D., 2.7 μm) P/N: 227-32015-03
In-Line Sample: Preparation Column:	Revive ILSP Pesticides Single 5 $ imes$ 2.1 mm Cartridge (RESTEK)
Mobile Phase A:	2 mM Ammonium formate-0.002 % Formic acid-Water
Mobile Phase B, C (Wash):	2 mM Ammonium formate-0.002 % Formic acid-Methanol
Gradient Program:	B conc. 3 % (0 min)-10 % (1 min)-55 % (3 min)- 100 % (10.5-12 min)-3 % (12.01-15 min)
Flowrate (A & B):	0.4 mL/min
Flowrate (C):	0 mL/min (0-5.49 min) – 1 mL/min (5.5-7.5 min) – 0.4 mL/min (7.51-11 min) - 0 mL/min (11.01-15
Switching Valve Position:	nin) load (0-6.49 min) – wash (6.5-10.5 min) – load (10.51-15 min)
Injection Volume:	2 μL (Co-injection 40 μL Water)
MS (LCMS-9050)	
lonization:	ESI (Positive, Negative)
TOF-MS:	<i>m/z</i> 50-800
Nebulizing Gas Flow:	2.0 L/min
Drying Gas Flow:	10.0 L/min
Heating Gas Flow:	10.0 L/min
DL Temp.:	150 °C
Block Heater Temp.:	300 °C
Interface Temp.:	200 °C
Probe Position:	+2 mm



# Analysis of Pesticides Standard Solution by the LCMS-9050

Table 2 shows the mass error of 79 compounds in pesticides standard solution diluted to 50 ppb. All the compounds were detected within  $\pm 1$  mDa mass error. Theoretical *m*/*z* values of pesticides were calculated using LabSolutions Insight Explore<sup>TM</sup>.

Fig. 5 shows the total ion current chromatogram (TICC) of the pesticides standard solution. Also, Fig. 6 shows the extracted ion chromatogram (XIC) of 79 compounds in pesticides standard solution and blank solvent.

Table 2 List of Pesticide Compounds

Compound	Molecular Formula	Selected	Theoretical	Mass Error	Retention Time
A anto as in sid			<u>m/z</u>	(mDa)	(min)
Acetamiprid	$C_{10}H_{11}CIN_4$	[/VI+H] <sup>+</sup>	223.0745	-0.2	4.81
Dimethoate	$C_5H_{12}NO_3PS_2$		230.0069	0.0	4.18
Bromacil	$C_9H_{13}BrN_2O_2$	[M-H]	259.0087	0.4	4.80
Propoxur	$C_{11}H_{15}NO_3$	[M+H] <sup>+</sup>	210.1125	-0.2	5.01
Isouron	$C_{10}H_{17}N_3O_2$	[M+H] <sup>+</sup>	212.1394	-0.1	4.97
Fluometuron	$C_{10}H_{11}F_3N_2O$	[M+H] <sup>+</sup>	233.0896	-0.1	4.81
Pyraclonil	$C_{15}H_{15}CIN_6$	[M+H] <sup>+</sup>	315.1120	-0.1	7.07
Metalaxyl/Metalaxyl-M*	$C_{15}H_{21}NO_4$	[M+H] <sup>+</sup>	280.1544	-0.1	6.51
Methidathion	$C_6H_{11}N_2O_4PS_3$	[M+H] <sup>+</sup>	302.9692	-0.2	7.17
Flumioxazin	C <sub>19</sub> H <sub>15</sub> FN <sub>2</sub> O <sub>4</sub>	[M+H] <sup>+</sup>	355.1089	0.0	8.50
Chlorbufam	C <sub>11</sub> H <sub>10</sub> CINO <sub>2</sub>	[M+H] <sup>+</sup>	224.0473	-0.5	5.91
Ethiprole	C <sub>13</sub> H <sub>9</sub> Cl <sub>2</sub> F <sub>3</sub> N <sub>4</sub> OS	[M+H] <sup>+</sup>	396.9899	0.3	5.95
Paclobutrazol	C <sub>15</sub> H <sub>20</sub> CIN <sub>3</sub> O	[M+H] <sup>+</sup>	294.1368	-0.1	6.18
Barban	$C_{11}H_9CI_2NO_2$	[M+H] <sup>+</sup>	258.0083	-0.3	6.68
Benthiavalicarb-isopropyl	C <sub>18</sub> H <sub>24</sub> FN <sub>3</sub> O <sub>3</sub> S	[M+H] <sup>+</sup>	382.1595	-0.1	6.74
Tiadinil	C <sub>11</sub> H <sub>10</sub> CIN <sub>3</sub> OS	[M-H]⁻	266.0160	0.0	6.55
Triadimenol	C14H18CIN3O2	[M+H] <sup>+</sup>	296.1161	-0.4	6.30
Triflumizole Metabolite	C <sub>12</sub> H <sub>14</sub> CIF <sub>3</sub> N <sub>2</sub> O	[M+H] <sup>+</sup>	295.0820	-0.2	4.81
Prometryn	$C_{10}H_{19}N_5S$	[M+H] <sup>+</sup>	242.1434	-0.2	6.66
Tetraconazole	C <sub>13</sub> H <sub>11</sub> Cl <sub>2</sub> F <sub>4</sub> N <sub>3</sub> O	[M+H] <sup>+</sup>	372.0288	-0.2	6.90
Flusilazole	C <sub>16</sub> H <sub>15</sub> F <sub>2</sub> N <sub>3</sub> Si	[M+H] <sup>+</sup>	316.1076	0.0	7.65
Bensulide	C14H24NO4PS3	[M+H] <sup>+</sup>	398.0678	-0.1	7.99
Flubendiamide	C <sub>23</sub> H <sub>22</sub> F <sub>7</sub> IN <sub>2</sub> O <sub>4</sub> S	[M-H] <sup>-</sup>	681.0160	0.1	7.27
Kresoxim-methyl	C18H19NO4	[M+H] <sup>+</sup>	314.1387	-0.1	8.30
Pvrazoxvfen	C20H16Cl2N2O3	[M+H] <sup>+</sup>	403.0611	0.1	8.90
Famoxadone	C <sub>22</sub> H <sub>18</sub> N <sub>2</sub> O <sub>4</sub>	[M+NH₄] <sup>+</sup>	392,1605	0.1	8.60
Phoxim	C12H15N2O3PS	[M+H] <sup>+</sup>	299.0614	-0.2	8.47
Trichlamide	C12H14Cl2NO2	[M-H]	338.0123	0.6	7.40
Metconazole	C17H22CIN3O	[M+H] <sup>+</sup>	320,1524	-0.1	7.56
Pyraclofos	C14H10CIN2O2PS	[M+H] <sup>+</sup>	361 0537	0.1	8 44

Compound	Molecular Formula	Selected	Theoretical m/z	Mass Error (mDa)	Retention Time
Bitertanol	C20H23N3O2	[M+H] <sup>+</sup>	338.1863	-0.2	7.94
Pyrazophos	C14H20N3O5PS	[M+H]+	374.0934	0.0	9.38
Diflufenican	C19H11F5N2O2	[M+H] <sup>+</sup>	395.0814	0.0	8.26
Pentoxazone	C17H17CIFNO4	[M+NH₄] <sup>+</sup>	371.1168	0.1	9.24
Tolfenpyrad		[M+H] <sup>+</sup>	384.1473	0.0	9.08
Pyributicarb	$C_{19}H_{22}C_{19}N_2O_2S$	[M+H] <sup>+</sup>	331 1475	0.1	9.47
Chlorpyrifos		[M+H] <sup>+</sup>	349 9336	-0.2	9.25
Ftoxazole		[M+H] <sup>+</sup>	360 1770	0.1	9 39
Cvenopyrafen	C24H21N2O2	[M+H] <sup>+</sup>	394 2489	0.1	945
Spirodiclofen	C21H24Cl2O4	[M+H] <sup>+</sup>	411 1125	-0.1	9.62
3-Hydroxycarbofuran	C12H15NO4	[M+H] <sup>+</sup>	238 1074	-0.1	4.06
Cymoxanil	C7H10N4O2	[M+H] <sup>+</sup>	199.0826	-0.4	4 35
Phosphamidon		[M+H] <sup>+</sup>	300 0762	-0.1	5.04
Terbacil		[M-H] <sup>-</sup>	215 0593	-0.2	4 86
XMC (3.5-xylyl methylcarbamate)	CioHioNOo	[M+H]+	180 1019	-0.2	5 25
Flutriafol	$C_{10}H_{13}NO_2$	[M+H]+	302 1100	0.2	5.25
Fensulfothion		[M+H]+	302.1100	0.1	5.70
Triforine (isomer-1)		[M+H]+	432 9321	0.0	5.68
Triforine (isomer-2)	$C_{10}H_14C_6N_4O_2$	[M+H]+	432.0321	0.2	5.00
Diethofencarb		[M+H]+	768 1544	-0.3	6.20
Fludiovonil	C14H21NO4	[NI-H]-	200.1344	-0.5	6.00
Mandipropamid		[M+H]+	412 1310	0.2	7.53
Pyriminobac-methyl (E)	$C_{23} H_{22} C_{11} O_4$	[M+H]+	362 13/7	-0.1	7.55
Malathian			221 0424	-0.1	7.00
Malatinon Promobutido dibromo	$C_{10}\Pi_{19}U_{6}PS_{2}$		221.0424	-0.1	7.54
Fluenicelide			254.1052	0.0	6.10
Triadimatan			204 1004	0.1	6.70
			294.1004	-0.1	0.79
Pramphop-methyl			312 0059	-0.1	7.55
Gerfentrezene ethul			312.0958	0.0	0.07
Carlentrazone-etnyi			429.0702	0.1	8.00 7.16
Dimethaneuryn	$C_{11}\Pi_{21}N_5S$		250.1591	-0.2	7.10
Pentniopyrad	$C_{16}H_{20}F_3N_3OS$	[M+H]	360.1352	0.0	6.88
l'ebuconazole	$C_{16}H_{22}CIN_{3}O$	[M+H] <sup>+</sup>	308.1524	0.0	7.24
Benalaxyi	$C_{20}H_{23}NO_3$	[M+H] <sup>+</sup>	326.1751	0.0	8.40
	$C_{15}H_{14}Cl_2N_2O_3$	[MI+H] <sup>+</sup>	341.0454	-0.2	8.19
ISOXATNION	C <sub>13</sub> H <sub>16</sub> NO <sub>4</sub> PS	[M+H] <sup>+</sup>	314.0611	0.0	8.56
Prochloraz	$C_{15}H_{16}CI_{3}N_{3}O_{2}$	[M+H]	3/6.0381	0.0	8.70
Pirimiphos-methyl	$C_{11}H_{20}N_3O_3PS$	[M+H]	306.1036	-0.1	8.25
Difenoconazole	$C_{19}H_{17}CI_2N_3O_3$	[M+H]	406.0720	-0.5	9.00
	$C_{20}H_{19}F_{3}N_{2}O_{4}$	[M+H] <sup>+</sup>	409.1370	0.1	8./5
Iriflumizole	C <sub>15</sub> H <sub>15</sub> CIF <sub>3</sub> N <sub>3</sub> O	[M+H]	346.0929	0.0	8.18
Amisulbrom	$C_{13}H_{13}BrFN_5O_4S_2$	[M+H]	467.9628	0.3	9.05
Protenotos	C <sub>11</sub> H <sub>15</sub> BrClO <sub>3</sub> PS	[M+H] <sup>+</sup>	374.9402	-0.2	8.63
Buprofezin	C <sub>16</sub> H <sub>23</sub> N <sub>3</sub> OS	[M+H]⁺	306.1635	0.0	8.80
Piperonyl butoxide	$C_{19}H_{30}O_5$	[M+NH <sub>4</sub> ] <sup>+</sup>	356.2432	-0.2	9.05
Butachlor	C <sub>17</sub> H <sub>26</sub> CINO <sub>2</sub>	[M+H] <sup>+</sup>	312.1725	-0.1	8.90
Quinoxyfen	C <sub>15</sub> H <sub>8</sub> Cl <sub>2</sub> FNO	[M+H] <sup>+</sup>	308.0040	0.0	9.26
Pyridaben	C <sub>19</sub> H <sub>25</sub> CIN <sub>2</sub> OS	[M+H] <sup>+</sup>	365.1449	0.2	10.13
Fenpropimorph	C <sub>20</sub> H <sub>33</sub> NO	[M+H]⁺	304.2635	-0.2	7.10

\*Pesticides standard solution contains 50 ppb each of metalaxyl and metalaxyl-M.



Fig. 5 Total Ion Current Chromatogram (TICC) of Pesticides Standard Solution



Fig. 6 Extracted Ion Chromatograms (XICs) of 79 Pesticide Compounds

#### ■ Linearity of Calibration Curve

Linearity of the calibration curve for each pesticide was evaluated by generating a 7-point calibration curve with the range 0.5-50 ppb in solvent (acetonitrile) and in spinach extract. Both in solvent and in extract, linearity showed very good results (coefficient of determination  $R^2$ : 0.99 or more) for all the

compounds. For 68 out of 79 compounds detection was within 2.5 ppb or less. Calibration curves for fludioxonil, which was detected in negative mode, in solvent and in extract are shown in Fig. 7 as an example, and calibration ranges for all the compounds are shown in Table 3.



Fig. 7 Calibration Curve of Fludioxonil (Left: in Solvent, Right: in Spinach Extract)

	Calibrati	on Range
Compound	(p	pD)
	in solvent	extract
Acetamiprid	1-50	2.5-50
Dimethoate	1-50	0.5-50
Bromacil	5-50	5-50
Propoxur	1-50	1-50
Isouron	0.5-50	0.5-50
Fluometuron	0.5-50	0.5-50
Pyraclonil	0.5-50	0.5-50
Metalaxyl/Metalaxyl-M*	0.5-50	0.5-50
Methidathion	1-50	1-50
Flumioxazin	5-50	5-50
Chlorbufam	25-50	25-50
Ethiprole	1-50	1-50
Paclobutrazol	1-50	1-50
Barban	5-50	5-50
Benthiavalicarb-isopropyl	0.5-50	0.5-50
Tiadinil	1-50	1-50
Triadimenol	2.5-50	2.5-50
Triflumizole Metabolite	0.5-50	0.5-50
Prometryn	0.5-50	0.5-50
Tetraconazole	0.5-50	0.5-50
Flusilazole	0.5-50	0.5-50
Bensulide	1-50	1-50
Flubendiamide	1-50	2.5-50
Kresoxim-methyl	0.5-50	1-50
Pyrazoxyfen	0.5-50	0.5-50
Famoxadone	2.5-50	2.5-50
Phoxim	1-50	0.5-50
Trichlamide	2.5-50	2.5-50
Metconazole	1-50	0.5-50
Pyraclofos	0.5-50	0.5-50
Bitertanol	2.5-50	2.5-50
Pyrazophos	0.5-50	0.5-50
Diflufenican	1-50	1-50
Pentoxazone	10-50	10-50
Tolfenpyrad	0.5-50	1-50
Pyributicarb	0.5-50	0.5-50
Chlorpyrifos	1-50	1-50
Etoxazole	0.5-50	0.5-50
Cyenopyrafen	0.5-50	0.5-50
Spirodiclofen	1-50	1-50

#### Table 3 Linear Ranges of 79 Pesticides

-

Compound	Calibratio (pr	on Range ob)
Compound	in solvent	in spinach extract
3-Hydroxycarbofuran	1-50	1-50
Cymoxanil	5-50	5-50
Phosphamidon	0.5-50	0.5-50
Terbacil	5-50	5-50
XMC (3,5-xylyl methylcarbamate)	2.5-50	2.5-50
Flutriafol	1-50	1-50
Fensulfothion	0.5-50	0.5-50
Triforine (isomer-1)	10-50	10-50
Triforine (isomer-2)	10-50	10-50
Diethofencarb	1-50	1-50
Fludioxonil	0.5-50	0.5-50
Mandipropamid	0.5-50	0.5-50
Pyriminobac-methyl (E)	0.5-50	0.5-50
Malathion	0.5-50	0.5-50
Bromobutide-debromo	1-50	1-50
Fluopicolide	0.5-50	1-50
Triadimefon	1-50	1-50
Flamprop-methyl	0.5-50	0.5-50
Bromobutide	1-50	2.5-50
Carfentrazone-ethyl	1-50	1-50
Dimethametryn	0.5-50	0.5-50
Penthiopyrad	0.5-50	0.5-50
Tebuconazole	1-50	0.5-50
Benalaxyl	0.5-50	0.5-50
Oxadiargyl	2.5-50	2.5-50
Isoxathion	0.5-50	0.5-50
Prochloraz	1-50	1-50
Pirimiphos-methyl	0.5-50	0.5-50
Difenoconazole	1-50	0.5-50
Trifloxystrobin	0.5-50	0.5-50
Triflumizole	0.5-50	0.5-50
Amisulbrom	25-50	25-50
Profenofos	0.5-50	0.5-50
Buprofezin	0.5-50	0.5-50
Piperonyl butoxide	0.5-50	0.5-50
Butachlor	5-50	5-50
Quinoxyfen	0.5-50	1-50
Pyridaben	0.5-50	0.5-50
Fenpropimorph	0.5-50	0.5-50

\*Metalaxyl and metalaxyl-M were not distinguished.

### Spike and Recovery Test

A spike and recovery test was performed using spinach extract to which 79 pesticides standard solution was spiked at 0.01 mg/kg per sample (concentration in pretreated sample solution was 2.5 ppb). The recovery rate was calculated at the concentration of 2.5 ppb, if not possible, at the concentration of 50 ppb. The results of recovery rate calculated by external standard method and reproducibility (n=5) are shown in Table 4, and the breakdown of the recovery rate is shown in Fig. 8.

Recovery rates were 70-120 % for 77 of the 79 compounds and %RSDs were less than 20 % for all the compounds. Good recovery rate and reproducibility were obtained without significant matrix inhibition by using the Revive ILSP column in sample preparation.



Fig. 8 Breakdown of Recovery Rate

Compound	Recovery Rate (%)	%RSD
Acetamiprid	82.6	11.4
Dimethoate	87.1	7.6
Bromacil	73.1	2.5
Propoxur	103.9	9.9
Isouron	107.1	1.5
Fluometuron	105.4	4.7
Pyraclonil	107.8	5.6
Metalaxvl/Metalaxvl-M*	108.6	3.6
Methidathion	110.5	10.3
Flumioxazin	65.5	1.7
Chlorbufam	85.0	12.7
Ethiorole	104 3	79
Paclobutrazol	101.5	8.4
Barban	92.1	4.0
Benthiavalicarh-isonronyl	108.2	1.0
Tiadinil	00.2	7.1
Triadimenol	100.4	18.6
Triflumizale Metabolite	100.4	10.0
	103.9	4.0
Totracopazala	08.1	5.7 E 1
	98.1	5.1
Flushazole	97.2	5.9
Bensulide	102.3	3.6
Flubendiamide	95.3	14.4
Kresoxim-methyl	95.9	9.9
Pyrazoxyten	116.9	4.9
Famoxadone	84.6	18.6
Phoxim	85.2	12.5
Trichlamide	80.7	17.2
Metconazole	95.2	5.3
Pyraclofos	115.7	5.9
Bitertanol	77.4	6.4
Pyrazophos	112.3	2.1
Diflufenican	83.6	5.6
Pentoxazone	97.9	3.1
Tolfenpyrad	101.3	4.6
Pyributicarb	102.9	2.6
Chlorpyrifos	104.2	13.3
Etoxazole	108.0	2.1
Cyenopyrafen	98.8	4.7
Spirodiclofen	110 5	6.7

#### Table 4 Recovery Rate and Reproducibility (%RSD) (n=5)

\*Metalaxyl and metalaxyl-M were not distinguished.

Recovery

Rate (%)

93.4 81.9

96.8

74.9

88.2

996

104.6

97.9

101.8

100.2

88.4

104.5

100.3

105.2

107.0

96.2

96.5

105.6

112.5

77.3

106.5

979

109.9

93.3

66.0

94.3

98.3

100.4

95.7

96.7

83.4

98.1

92.1

94.1

87.6

100.9

100.4

95.0

103.7

%RSD

14.6

2.6

6.5

3.7

7.7

84

4.0

7.5

43

6.8

4.5

3.2

3.5

5.1

12.2

7.2

6.0

8.4

9.9

60

52

5.7

6.8

3.0

12.0

3.9

6.2

5.1

7.7

3.1

3.1

4.4

3.2

4.1

7.6

3.0

3.2

4.9

3.9

#### ■ Conclusion

Screening analysis of pesticides was performed using positive/negative ion polarity switching by the LCMS-9050. The on-line sample preparation method with a Revive ILSP column made it possible to speed up and simplify the preparation process. Screening analysis of pretreated spinach samples with pesticides standard spiked provided good results for recovery rate, reproducibility and linearity. For sensitivity, 68 out of 79 compounds were detected at 2.5 ppb or less. This method is applicable for the analysis of residual pesticides in other crops.

<References>

1) Lupo, S.A., Romesberg, R.L., Lu, X., 2020. Automated inline pigment removal for the analysis of pesticide residues in spinach by liquid chromatography tandem mass spectrometry, J. Chromatogr. A 1629, 461477.

Nexera, Shim-pack, and LabSolutions Insight Explore are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.



Shimadzu Corporation

For Research Use Only, Not for use in diagnostic procedures.

01-00547-EN First Edition: Apr. 2023

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these

The products in your country. The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See http://www.shimadzu.com/about/trademarks/index.html for details. See http://www.shimadzu.com/about/trademarks/index.html for details. Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not

www.shimadzu.com/an/

they are used with trademark symbol "TM" or "@". Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.