

## Application News

Liquid Chromatograph-Mass Spectrometer LCMS-8045

### High Throughput Testing Method with Newborn's DBS for Hemoglobinopathies by LCMS-8045 with ZenTech Targeted MS/MS Hemo Device

Devika Tupe<sup>1</sup>, Bhaumik H Trivedi<sup>1</sup>, Pratap Rasam<sup>1</sup>, Avinash Prabhavalkar<sup>2</sup>, Madeleine Boulanger<sup>3</sup>, Martin Greffe<sup>3</sup>, Luc Voorn<sup>3</sup>

<sup>1</sup>Shimadzu Analytical (India) Pvt. Ltd. <sup>2</sup>Synergy Medical Systems LLP, <sup>3</sup>Zentech

#### User Benefits

- ◆ End-to-end solution provided by Shimadzu and ZenTech
- ◆ High throughput, sensitive & robust methodology for hemoglobinopathies using Newborn's dried blood spots (DBS)
- ◆ Analysed in MRM-mode targeting the unique peptides of hemoglobin variants.

#### ■ Introduction

##### Hemoglobin

Hemoglobin (Hb) (Fig. 1) is a tetramer composed of two  $\alpha$ -globin and two non- $\alpha$ -globin chains working in conjunction with heme to transport oxygen in the blood. Normal adult hemoglobin (HbA) is designated  $\alpha_2\beta_2$ . Hemoglobin Variant is derived from gene abnormalities affecting the  $\alpha$ -globin genes (HBA1 or HBA2) or  $\beta$ -globin (HBB) structural genes (exons). More than a thousand hemoglobin variants have been identified relative to changes in the globin chains. Qualitative changes correspond to amino acid substitutions resulting in hemoglobinopathies. [1,2,3]

##### Hemoglobinopathies

Hemoglobinopathies are genetic defects that result in one of the chains of the hemoglobin molecule. These are usually single-gene disorders, such as sickle cell disease, and they are usually inherited as autosomal codominant traits. To date, over 1000 hemoglobinopathies have been found but fortunately, only a few have clinical significance. Although there is a formal system for naming hemoglobinopathies, most hemoglobinopathies are known by a letter or identifier from the location where the hemoglobin variant was found or a combination of these. [4] It is important to be able to structurally identify and quantitate these variants in order to detect abnormal hemoglobin quickly and precisely. The conventional technologies have reported false (positive/negative) results as they are very sensitive to pH and a minute deviation impact the results.

This application highlights high throughput, sensitive and reproducible testing method with Newborn's DBS for hemoglobinopathies using Shimadzu LCMS-8045 (Fig. 2) with Targeted MS/MS Hemo device (ZenTech s.a., Avenue du Pré-Aily 10, 4031 Angleur, Belgium) (Fig. 3).




Fig. 1 Structure of Normal Hemoglobin

Liquid chromatography mass spectrometer (LC-MS/MS) methods with high throughput, high sensitivity and specificity have been developed for analysis of Hb variants. The application is capable of analysing  $\alpha$ -thalassemia,  $\beta$ -thalassemia and Hb variants like HbS, HbC, HbE, HbO<sup>Arab</sup>, HbD<sup>Punjab</sup> having fast turnaround time required for large numbers of samples. The techniques have produced accurate analysis that would have been missed by the traditional methods.



Fig. 2 Nexera™ X3 with LCMS-8045



Fig. 3 Targeted MS/MS Hemo device

#### ■ Experimental

##### Sample preparation

The LCMS-8045 triple quadrupole mass spectrometer was coupled to a Nexera X3 UHPLC system (Fig. 2). The detection of normal and variant Hb was performed using the commercially available Targeted MS/MS Hemo device kit. Prior to the LC-MS/MS analysis, the Hb is extracted and denatured from dried blood spots samples before being digested by trypsin. The flow of sample preparation is described (Fig. 4). Controls are provided in the device to ensure the proper course of extraction, denaturation and digestion steps as well as the correct functioning of the mass spectrometer for all peptides from normal and variant Hb. The samples were analysed in MRM-mode. Analytical conditions are listed in Table 1. The optimized MRM transitions are listed in Table 2.




Fig. 4 Flow of sample preparation

Table 1 Analytical Conditions

Mass Spectrometer	: LCMS-8045
Ionisation	: Electrospray Ionization (ESI), positive
Interface Voltage	: 3.5 kV
Desolvation Line	: 200 °C
Heating Gas	: 15 L/min
Interface Temp.	: 400 °C
Nebulizing Gas	: 3 L/min
Drying Gas	: 5 L/min
Heat Block	: 200 °C
CID	: 270 kPa
UHPLC	: Nexera X3
Pump A	: 0.1% Formic acid in Water
Pump B	: 0.1% Formic acid in Acetonitrile
Column Oven	: RT
Injection Volume	: 5 µL
Analytical column	: None
Gradient	: Isocratic 50% A / 50% B
Flow program	: 0.0 min 0.4 mL/min 0.5 min 0.4 mL/min 0.6 min 0.6 mL/min 0.9 min 0.6 mL/min 1.0 min 0.4 mL/min
Autosampler temperature	: 5 °C

Table 2 MRM Transitions

Target compound	Parent ion (m/z)	Product ion 1 (m/z)	Product ion 2 (m/z)
α-peptide 1	536.5	446.2	680.3
α-peptide 2	627.0	233.3	261.2
β-peptide 1	477.0	237.2	502.2
β-peptide 2	658.0	313.3	758.3
β-peptide 3	690.0	378.3	501.3
β-peptide S	462.0	237.2	472.1
β-peptide C	694.5	237.1	244.2
β-peptide E	459.0	214.2	360.2
β-peptide O <sup>Arab</sup>	625.6	249.0	501.1
β-peptide D <sup>Punjab</sup>	689.6	377.1	276.1
γ-peptide 1	550.0	251.3	634.2
δ-peptide 1	480.3	390.3	688.4

Table 3 Sample Status with the Targeted MS/MS Hemo Device

Expected Phenotype	Result Obtained	Number of Samples
Normal	Normal	508
HbS Carrier	HbS Carrier	5
HbE Carrier	HbE Carrier	1
HbD <sup>Punjab</sup> Carrier	HbD <sup>Punjab</sup> Carrier	1




Fig. 5 Example of Chromatogram

## ■ Results

Using the Targeted MS/MS Hemo device from ZenTech, the testing method with Newborn's DBS for hemoglobinopathies can be performed quite easily and very fast, in less than one minute, showing a good peak shape (Fig. 5).

All the dried blood spots samples analysed came from routine testing labs. The status of these DBS are known, and their clinical status has been successfully correlated and assigned (Table 3).

## ■ Conclusion

Similar application note is available on LCMS-8050 by Shimadzu Europe but this application note was developed to showcase that Targeted MS/MS Hemo device from ZenTech can be very well reproduced on Shimadzu LCMS-8045 for testing method with Newborn's DBS of hemoglobinopathies. Needless to say, that ease of sample preparation, rugged technology, ambiguous outcome and overall performance, will appeal to laboratories.

## ■ References

- [1] A.D.A.M. Medical Encyclopedia [Internet]. Atlanta (GA): A.D.A.M., Inc.; ©2005 Thalassemia; updated 2014 Feb 24; cited Apr 10 2015. Available from : <http://www.nlm.nih.gov/medlineplus/ency/article/000587.htm>
- [2] Tietz: Fundamentals of Clinical Chemistry, (6th Edition). Burtis, C. A., Ashwood, E. R., Bruns, D.E. (Eds.). Haryana, India: Elsevier. 2008
- [3] National Library of Medicine. HBB 2009. Retrieved from <http://ghr.nlm.nih.gov/gene/HBB>
- [4] AACC, Hemoglobinopathies, Trefor Higgins DEC.29.2014, <https://www.aacc.org/science-and-research/clinical-chemistry-trainee-council/trainee-council-in-english/pearls-of-laboratory-medicine/2014/hemoglobinopathies>

Nexera is a trademark of Shimadzu corporation in Japan and/or other countries.



Shimadzu Corporation

[www.shimadzu.com/an/](http://www.shimadzu.com/an/)

Shimadzu Analytical (India) Pvt.Ltd.

[www.shimadzu.in](http://www.shimadzu.in)

### For Research Use Only. Not for use in diagnostic procedures.

This has not been approved or certified as a medical device under the Pharmaceutical and Medical Device Act of Japan.

It cannot be used for the purpose of medical examination and treatment or related procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu.

See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

06-SAIP-LC-055-EN

First Edition: Mar. 2023

## Related Products

Some products may be updated to newer models.



> LCMS-8060

High Performance Liquid  
Chromatograph Triple Quadr...

## Related Solutions

> Small Molecule  
Pharmaceutical

> Price Inquiry

> Product Inquiry

> Technical Service /  
Support Inquiry

> Other Inquiry