

Gas Chromatograph HS-20 NX USTL/Brevis™ GC-2050

## Efficient Analysis of Residual Solvents in Pharmaceuticals Using the Compact Model, Brevis GC-2050 (1)

### —JP18 and USP467, Water-Soluble Samples—

Yui Higashi, Eisuke Kobayashi

#### User Benefits

- ◆ The slim and compact design of the Brevis GC-2050 enables the expansion of the number of operational units in the laboratory, allowing for efficient high-throughput analysis.
- ◆ Brevis GC-2050 can perform analysis using alternative carrier gases.
- ◆ Analysis can be performed with *tert*-butyl alcohol and cyclopentyl methyl ether, which have been newly added as Class 2 solvents in ICH Q3C (R8).

#### Introduction

Residual solvents in pharmaceuticals are defined as organic volatile chemical substances that are used or produced in the manufacture of drug substances or additives or in the preparation of drug products. In the Japanese Pharmacopoeia 18<sup>th</sup> Edition (JP18) or the United States Pharmacopoeia (USP) General Chapters <467> Residual Solvents, residual solvents are classified as Class 1 to 3, according to their risk to human health, and the headspace GC method is mainly used to analyze them. The carrier gas that is normally used is He. However, He supply shortages have become an issue recently, so there is a demand to perform analysis using alternative carrier gases, such as H<sub>2</sub> or N<sub>2</sub>.

This article introduces the results of analysis of Class 1 and 2 water-soluble samples using the compact design Brevis GC-2050. For the Procedure A, H<sub>2</sub> and N<sub>2</sub> in addition to He were used. When using an alternative carrier gas, it is necessary to first verify the operation based on USP General Chapter <1467>.

#### Instrument Configuration and Analytical Conditions

Table 1 Analytical Conditions of Water-Soluble Samples

GC Analytical Conditions (Procedure A and B)	
Model	: Brevis GC-2050
Detector	: FID (Flame Ionization Detector)
Column	: A) SH-I-624Sil MS (0.32 mm I.D. × 30 m, d.f. = 1.8 μm) B) SH-PolarWax (0.32 mm I.D. × 30 m, d.f. = 0.25 μm)
Column Temp.	: A) 40 °C (20 min) – 10 °C/min – 240 °C (20 min) Total 60 mins B) 50 °C (20 min) – 6 °C/min – 165 °C (20 min) Total 59.17 mins
Injection Mode	: A) Split 1:5 B) Split 1:10
Carrier Gas Controller	: Linear velocity (He, N <sub>2</sub> , H <sub>2</sub> )
Linear Velocity	: 35 cm/sec
Detector Temp.	: 250 °C
FID H <sub>2</sub> Flowrate	: 32 mL/min
FID Makeup Flowrate	: 24 mL/min (N <sub>2</sub> )
FID Air Flowrate	: 200 mL/min
Injection Volume	: 1 mL

HS-20 NX Analytical Conditions (Same for Procedure A and B)	
Model	: HS-20 NX USTL (Ultra Short Transfer Line)
Oven Temp.	: 80 °C
Sample Line Temp.	: 110 °C
Transfer Line Temp.	: 120 °C
Vial Shaking Level	: Off
Vial Volume	: 20 mL
Vial Equilibrating Time	: 60 min
Vial Pressurizing Time	: 1 min
Vial Pressure	: 75 kPa
Loading Time	: 0.5 min
Load Equilib. Time	: 0 min
Needle Flush Time	: 5 min



Fig. 1 HS-20 NX USTL (Ultra Short Transfer Line) + Brevis™ GC-2050

#### Class 1 Standard Solution Analysis (Water-Soluble Samples)

Figs. 2 and 3 show the analysis results for the Class 1 standard solution.




Fig. 2 Class 1 Standard Solution Chromatogram (Water-Soluble Sample) Using Procedure A

- 1, 1,1-Dichloroethane,
- 2, 1,1,1-Trichloroethane,
3. Carbon tetrachloride,
4. Benzene,
5. 1,2-Dichloroethane




Fig. 3 Class 1 Standard Solution Chromatogram (Water-Soluble Sample) Using Procedure B


- 1, 1,1-Dichloroethane,
2. 1,1,1-Trichloroethane+ Carbon tetrachloride,
3. Benzene,
4. 1,2-Dichloroethane

## ■ Class 2 Standard Solution Analysis (Water-Soluble Samples)

Fig. 4 shows the analysis results for Procedure A, and Fig. 5 shows the analysis results for Procedure B (Class 2A: black, Class 2B: pink, TBA, CPME, MiBK: blue). For system suitability, JP18 specifies that “the resolution between acetonitrile and methylene chloride in the Class 2 mixture A standard solution is not less than 1.0” when using Procedure A, and “the resolution between acetonitrile and cis-1,2-dichloroethene in the Class 2 mixture A standard solution is not less than 1.0” when using Procedure B. Satisfactory results were obtained with both procedures.


Note: The resolutions shown in the Figs. 4 and 5 are reference values and not guaranteed.

Note: A mixture of standard samples of *tert*-butyl alcohol (TBA), cyclopentyl methyl ether (CPME), and methyl isobutyl ketone (MiBK) was separately prepared to the prescribed concentration.



- 1: Methanol
- 2: Acetonitrile
- 3: Methylene chloride (DCM)
- 4: *trans*-1,2-Dichloroethylene
- 5: *cis*-1,2- Dichloroethylene
- 6: Tetrahydrofuran
- 7: Cyclohexane
- 8: Methyl cyclohexane
- 9: 1,4-Dioxane
- 10: Toluene
- 11: Chlorobenzene
- 12: Ethylbenzene
- 13: *m,p*-Xylene
- 14: *o*-Xylene
- 15: Cumene
- 16: Hexane
- 17: Nitromethane
- 18: Chloroform
- 19: 1,2-Dimethoxyethane
- 20: Trichloroethene
- 21: Pyridine
- 22: Methyl butyl ketone
- 23: Tetraline
- 24: *tert*-Butyl alcohol(TBA)
- 25: Cyclopentyl methyl ether (CPME)
- 26: Methyl isobutyl ketone (MiBK)

Fig. 4 Class 2 Standard Solution Chromatogram (Water-Soluble Sample) Using Procedure A



- 1: Cyclohexane
- 2: Methyl cyclohexane
- 3: *trans*-1,2-Dichloroethene + Tetrahydrofuran
- 4: Methanol
- 5: Methylene chloride
- 6: *cis*-1,2-Dichloroethene
- 7: Acetonitrile
- 8: Toluene
- 9: 1,4-Dioxane
- 10: Ethylbenzene
- 11: *p*-Xylene
- 12: *m*-Xylene
- 13: Cumene
- 14: *o*-Xylene
- 15: Chlorobenzene
- 16: Hexane
- 17: 1,2-Dimethoxyethane
- 18: Trichloroethene
- 19: Chloroform
- 20: Methyl butyl ketone
- 21: Nitromethane
- 22: Pyridine
- 23: Tetraline
- 24: *tert*-Butyl alcohol(TBA)
- 25: Cyclopentyl methyl ether (CPME)
- 26: Methyl isobutyl ketone (MiBK)

Fig. 5 Class 2 Standard Solution Chromatogram (Water-Soluble Sample) Using Procedure B

## ■ Analysis of Class 1 and 2 Standard Samples (Water-Soluble Samples) by Procedure A Using H<sub>2</sub> or N<sub>2</sub> Carrier Gas

Figs. 6 to 9 show the separation patterns for Procedure A using the alternative gases H<sub>2</sub> and N<sub>2</sub>. When using an alternative carrier gas, the system suitability should be checked before performing the operation.

### ■ Conclusion

Even though the Brevis GC-2050 is small and compact, it is capable of analyzing residual solvents in pharmaceuticals in accordance with the JP18 and USP General Chapters <467> Residual Solvents. That compactness means the number of units installed in a laboratory can be increased compared with high-end models, so residual solvents in pharmaceuticals can be efficiently analyzed.

## Examples of Analysis of Class 1 and 2 Standard Samples (Water-Soluble Samples) Using an Alternative Carrier Gas




Fig. 6 Class 1 Standard Solution Chromatogram (Water-Soluble Sample) by Procedure A Using N<sub>2</sub> Carrier Gas

1. 1,1-Dichloroethane, 2. 1,1,1-Trichloroethane, 3. Carbon tetrachloride, 4. Benzene, 5. 1,2-Dichloroethane




Fig. 7 Class 1 Standard Solution Chromatogram (Water-Soluble Sample) by Procedure A Using H<sub>2</sub> Carrier Gas

1. 1,1-Dichloroethane, 2. 1,1,1-Trichloroethane, 3. Carbon tetrachloride, 4. Benzene, 5. 1,2-Dichloroethane




Fig. 8 Class 2 Standard Solution Chromatogram (Water-Soluble Sample) Using Procedure A with N<sub>2</sub> Carrier Gas




Fig. 9 Class 2 Standard Solution Chromatogram (Water-Soluble Sample) Using Procedure A with H<sub>2</sub> Carrier Gas  
\* The Class 2A reagent used in the H<sub>2</sub> carrier gas analysis shown in Fig. 9 included MiBK.

- 1: Methanol
- 2: Acetonitrile
- 3: Methylene chloride (DCM)
- 4: *trans*-1,2-Dichloroethylene
- 5: *cis*-1,2- Dichloroethylene
- 6: Tetrahydrofuran
- 7: Cyclohexane
- 8: Methyl cyclohexane
- 9: 1,4-Dioxane
- 10: Toluene
- 11: Chlorobenzene
- 12: Ethylbenzene
- 13: *m,p*-Xylene
- 14: *o*-Xylene
- 15: Cumene
- 16: Hexane
- 17: Nitromethane
- 18: Chloroform
- 19: 1,2-Dimethoxyethane
- 20: Trichloroethene
- 21: Pyridine
- 22: Methyl butyl ketone
- 23: Tetraline
- 24: *tert*-Butyl alcohol(TBA)
- 25: Cyclopentyl methyl ether (CPME)
- 26: Methyl isobutyl ketone (MiBK)

Brevis is a trademark of Shimadzu Corporation in Japan and other countries.



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country. The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

01-00516-EN

First Edition: Sep. 2023