

Application News

Measuring CH₄ and CO₂ Concentrations in Cow Breath for Estimating Total Methane Emissions

Minako Tanaka, Yuya Kurata, and Kazuyoshi Nakajima

User Benefits

- ◆ The transportable gas analyzer can be placed near a feed trough and continuously measure the concentrations of CH₄ and CO₂ in the breath of cows.
- ◆ The gas analyzer is simple to operate with an onboard sample pretreatment unit that includes a built-in sample pump, filter, and moisture removal system.
- ◆ The analyzer can save data on a USB flash drive to facilitate data manipulation on a computer or sharing data between groups.

Introduction

The greenhouse effect of methane (CH₄) is approximately 28 times that of carbon dioxide (CO₂)^{*1}. Cow belching is a significant source of CH₄ emissions and releases enteric CH₄ generated by fermentation in the cow's digestive system. Enteric CH₄ released through cow exhalation or eructation accounts for a significant proportion of the greenhouse gas emissions in agriculture, and concerted efforts are underway to reduce these emissions through feed development, changes in living conditions, and selective breeding.

The standard method for measuring enteric CH₄ emissions from cows is to place the animal in a large, specially designed chamber for several days and measure the total enteric CH₄ emitted over that time. While this method provides extremely accurate measurements, equipment and labor requirements make it unsuitable for measuring CH₄ emissions from large numbers of animals. In 2022, Japan's National Agriculture and Food Research Organization (NARO) published a manual that describes a method of estimating CH₄ emissions based on the ratio of CH₄ to CO₂ (CH₄/CO₂) measured in the breath of cows¹⁾. This method has attracted interest for its relative practicality as it does not require large-scale, dedicated facilities and can be used for single, short-term measurements and to collect data from multiple animals.

Experiment Overview

Estimating enteric CH₄ emissions from cows requires a complete analysis system that performs all parts of the analysis, from sampling of cow breath to measuring gas concentrations. As the first step in creating such a system, this Application News describes an example experiment in which the CGT-7100 transportable gas analyzer was used to collect data on CH₄ and CO₂ concentrations over time in cow breath.

CH₄ Generation in the Cow's Digestive System

Fig. 1 illustrates how CH₄ is generated in the digestive system of cows. The digestive system of ruminant animals such as cows and goats includes a stomach compartment called the rumen. In cows, the rumen is the first stomach compartment of the digestive system. The rumen is home to many different types of microorganisms with a variety of roles, one being to help the cow digest food in a process called rumen fermentation. The microorganisms ferment the cow's diet, generating H₂ in this process. Other microorganisms then metabolize this H₂ into CH₄, which is released from the cow by belching.

CGT-7100


The CGT-7100 is a transportable, all-in-one continuous gas analyzer. The compact chassis of the CGT-7100, which contains an onboard sample pretreatment unit with an electronic cooler for dehumidifying, a filter for dust removal, and a pump to draw sample gas into the analyzer, offers a standalone solution for gas analysis^{*3}.

The CGT-7100 can measure up to three gases: two gases from among CO₂, CO, and CH₄ by non-dispersive infrared (NDIR) absorption, and oxygen measured by an O₂ meter (optional) with a zirconia-based limiting current sensor. The concentration of each gas can be monitored over time across a wide concentration range from ppm to volume percentage levels.


With a minimum sample gas flowrate of 100 mL/min, the CGT-7100 can monitor gas concentrations in small sample gas volumes for a wide range of potential applications.

Non-Dispersive Infrared (NDIR) Absorption

The CGT-7100 uses non-dispersive infrared (NDIR) absorption for measurements. Heteroatomic gases, such as CO₂, have unique infrared absorption spectra. The concentration of individual heteroatomic gases in a mixed sample gas can be determined by passing infrared light through the gas sample and then detecting the amount of infrared light transmitted in absorption bands that are unique to each gas species. The NDIR analyzer uses a proven optical system that has been used in on-line gas analyzers and gives robust and stable measurements.




(1) Large numbers of anaerobic microorganisms (●) in the rumen.



(2) Proteins (●) and carbohydrates (●) in feed are decomposed and fermented by the microorganisms.

The microbial decomposition and fermentation process generates volatile fatty acids,^{*2} hydrogen (●), and carbon dioxide.



(3) Other microorganisms metabolize the hydrogen (●), combining it with carbon to form methane (●), which is released from the animal.

Fig. 1 Illustration of CH₄ Generation in Cow Digestive System



Fig. 2 Shimadzu CGT-7100 Transportable Gas Analyzer

*1 Figure taken from IPCC Fifth Assessment Report

*2 Acetic acid, propionic acid, butyric acid, etc.

*3 Dependent on analytical conditions. Contact a Shimadzu representative for details.

Method

The flow of measurement is shown in Fig. 3. An acrylic board enclosure was constructed around the feed trough to reduce the effects of wind and breath dispersion. A filter was also installed near the sampling line inlet to protect the system from airborne feed particles. Air was sampled from one location in the feed trough and the sample gas was aspirated into the analyzer at a flowrate of 400 mL/min.

The CGT-7100 was operated using the conditions shown in Table 1. The measurement range was set to 0 to 500 ppm for CH₄ and 0 to 1 vol % for CO₂. The concentrations of CH₄ and CO₂ were measured continuously for around 15 minutes while the cow was feeding.




Fig. 3 Experimental Setup

Table 1 Analytical Conditions

| | |
|----------------------|--|
| Analyzer: | CGT-7100 |
| Measured Components: | CH ₄ , CO ₂ |
| Measurement Range: | CH ₄ : 0–500 ppm CO ₂ : 0–1 vol % |
| Flowrate Type: | Low flowrate (100–400 mL/min, variable) |

Results

The data collected during the experiment are shown in Fig. 4. The CGT-7100 successfully measured the concentration of CH₄ and CO₂ in cow breath over time. The sudden rises in CH₄ concentration observed around once every minute are believed to coincide with the cow belching. Some measurements exceeded the upper limit of the measurement ranges set for this experiment. However, all concentrations were within the measurable range of the CGT-7100 analyzer and could be captured in future experiments by setting the measurement ranges with higher upper levels.

A fourth degree polynomial approximation of the results was created to reveal trends in the data. The approximation shows that measured CH₄ levels fell in the middle of the feeding time. This fall is believed to be caused by the cow moving its head, which changed how the cow's breath diffused in the feed trough area. Accordingly, the accuracy of estimated CH₄ emissions can be increased by taking measures to reduce the impact of breath diffusion on sampling. This could be achieved by using multiple sampling points located on the top and bottom and left and right sides of the feed trough area.¹⁾ Multiple sampling points would also reduce measured gas levels due to dilution of the sample gas by air collected from points further away from the mouth of the cow, but should also reduce the impact of the position of the cow's head on sampling and offer more stable measurements.

This experiment demonstrates that the CGT-7100 can measure CH₄ and CO₂ concentrations in cow breath and provides valuable insights about the measurement range and sampling method for future experiments intending to estimate CH₄ emissions in cow breath.

<References>

1) National Agriculture and Food Research Organization (2022). Manual for Estimating Methane Emissions from Fermentation in Cow Rumen

https://www.naro.go.jp/publicity_report/publication/pamphlet/tech-pamph/152088.html

NILGS:Development of a new calculation equation for methane emissions from cattle and its manualization (naro.go.jp)




Fig. 4 CH₄ and CO₂ Concentration in Cow Breath



Fig. 5 Experimental Setup

Conclusion

This Application News presents an example experiment that uses the CGT-7100 to measure the concentration of CH₄ and CO₂ in cow breath. The analysis demonstrates that the CGT-7100 can measure CH₄ and CO₂ concentrations in cow breath and provides valuable insights for future experiments aimed at estimating CH₄ emissions in cow breath. The results of this example experiment suggest that with adjustments to the sampling method, CH₄ emissions in cow breath can be estimated based on CH₄ and CO₂ measurements collected using the CGT-7100.

Related links

For more information about the CGT-7100, please see the product website below.

[CGT-7100: SHIMADZU \(Shimadzu Corporation\)](https://www.shimadzu.com/about/trademarks/index.html)

For information about the NOA-7100 transportable NO_x-O₂ analyzer, please see the link below.

[NOA-7100: SHIMADZU \(Shimadzu Corporation\)](https://www.shimadzu.com/about/trademarks/index.html)



SHIMADZU

Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

01-00755-EN

First Edition: Jul. 2024

➤ Please fill out the survey

Related Products

Some products may be updated to newer models.



➤ CGT-7100

Transportable Gas Analyzer

Related Solutions

➤ Environment

➤ Greenhouse Gas &
CCUS

➤ Price Inquiry

➤ Product Inquiry

➤ Technical Service /
Support Inquiry

➤ Other Inquiry