

Application News

GCMS-TQTM8050 NX Gas Chromatography Mass Spectrometer
HS-20 NX (Trap Model) Headspace Autosampler

Ethylene Oxide and 2-Chloroethanol in Various Rice Samples by Dynamic Headspace and GC-MS/MS

Elgin Ting¹, Cynthia Iahey¹,
1 Shimadzu (Asia Pacific) Pte Ltd.

User Benefits

- ◆ Simple sample preparation procedure, direct sample analysis
- ◆ Sample recovery within 70 % to 130 % without matrix-matched calibration curves

■ Introduction

Ethylene oxide (EO) is sometimes used as a fumigant on food products to control pests and insects during storage and transportation. However, residue of EO might remain in the products, raising concerns about the potential health risks associated with it. Therefore, it is essential to monitor and analyze ethylene oxide levels in food product to ensure compliance with safety regulations and protect consumer health. In this study, both ethylene oxide (EO) and 2-chloroethanol (2-CE) in rice were analyzed. 2-CE, a potential byproduct of EO fumigation, is also a compound of interest due to its toxicity and potential health implications. By utilising dynamic headspace and GC-MS/MS, precise and sensitive detection of both EO and 2-CE residues in rice samples can be achieved. The limit of quantitation is 5 ng/g with area %RSD (n=5) for EO and 2-CE at 5.6 % and 4.9 %, respectively. The spike recovery study at concentrations of 10 ng/g, 25 ng/g and 100 ng/g yielded recoveries within 70 %-130 %. This comprehensive analytical approach enables accurate quantification and identification of these compounds, significantly contributing to the overall safety and quality assurance of the food supply chain.

■ Experimental

Analysis condition

The instruments utilised in this experiment were GCMS-TQ8050 NX and HS-20 NX Trap (Shimadzu Corporation, Japan), as shown in Figure 1. HS-20 NX Trap, equipped with an absorbent trap, enables the dynamic headspace technique which enhances sensitivity compared to static headspace. Since the headspace technique was employed, the sample preparation procedure was greatly simplified compared to manual extraction followed by liquid injection. Additionally, HS-20 NX series features a heated, inert and exceptionally short transfer line to prevent analytes adsorption. The highly efficient detector and patented ion source technology in GCMS-TQ8050 NX provides the best sensitivity for ultra trace analysis, making it suitable for this application.

The analytical conditions are presented in Table 1.



Figure 1: HS-20 NX Trap autosampler with GCMS-TQTM8050 NX.

Table 1: Dynamic headspace and GC-MS/MS analytical conditions for analysis of EO and 2-CE analysis in rice.

Headspace Autosampler Parameter (HS-20 NX, Trap Model)	
Oven Temp.	70 °C
Trap Cooling Temp.	-20 °C
Equilibrating Time	5.0 min
Pressurization Gas	Helium
GC Parameter	
Carrier Gas	Helium
Injection Mode	Split mode, split ratio 15
Flow Control Mode	Linear Velocity, 44.0 cm/s
Column	SH-I-624 Sil MS (60 m x 0.25 mm ID x 1.4 µm df) [P/N: 227-36076-01]
MS/MS Parameter	
Ion Source Temp.	250 °C
Interface Temp.	235 °C

Standard and sample preparations

Ethylene oxide-d4 and 2-chloroethanol-d4, were obtained from Cambridge Isotope Laboratories, whereas ethylene oxide and 2-chloroethanol were from Merck. Acetonitrile was purchased from Kanto Chemical.

The procedure for preparing calibration standards, samples and spiked samples for recovery check was carried out by following Shimadzu Application News No. 04-AD-0299-EN [1]. Rice samples were ground before any analysis and spiking. In the current experiment, the recovery study was conducted at different concentrations: 10 ng/g, 25 ng/g, and 100 ng/g. As water content in rice samples is low, magnesium sulfate was not required in this application.

■ Result

Sensitivity, Repeatability and Linearity

Figure 2 shows the MRM mass chromatograms of EO and 2-CE at 5 ng/g of the standard sample. The limit of quantitation (LOQ) for both compounds are 5 ng/g with signal-to-noise ratios (S/N) greater than 10 (Table 2). The %RSD (n=5) of the area count for both compounds at 5 ng/g were less than 6 % (Table 3).




Figure 2: MRM mass chromatograms of: (a) EO and (b) 2-CE at 5 ng/g in the standard sample

Table 2: Signal-to-noise ratios of EO and 2-CE at 5 ng/g.

Compound Name	S/N	
	EO	2-CE
Injection 1	18	64
Injection 2	17	55
Injection 3	18	52
Injection 4	20	65
Injection 5	23	65

Table 3: Area repeatability %RSD (n=5) at 5 ng/g and linearity.

Compound Name	Repeatability	Linearity
	Area %RSD (n=5)	R ²
EO	4.6	0.9998
2-CE	5.6	0.9974


**Figure 3:** Five-point calibration curves of EO and 2-CE

Figure 3 shows 5-point internal standard calibration curves (5, 10, 25, 50 and 100 ng/g). The R² value for EO is >0.999 whereas 2-CE is >0.997 (Figure 3 and Table 3), demonstrating the capability of the system and method to achieve good linearity within the concentration range of 5 ng/g to 100 ng/g.


Quantitation and Recovery Results

Four different types of rice were purchased from local supermarkets and analyzed for the presence of EO and 2-CE. The results indicated that none of the samples contained EO and 2-CE (Table 4).

The recovery for all samples spiked with 10 ng/g, 25 ng/g, and 100 ng/g of each compound, were within 70 %-130 % range (Table 4). Figure 4 displays the MRM mass chromatograms for EO and 2-CE in brown rice spiked with 10 ng/g of each compound. The results indicate that the combination of dynamic headspace with GC-MS/MS and the internal standard method is suitable for analysing EO and 2-CE in rice. These findings suggest that matrix-matched calibration may not be necessary for these specific four rice samples. However, if recovery falls outside the 70 % to 130 % range due to challenging interferences in other rice matrices, employing matrix-matched calibration curves could be considered to improve the result.

Table 4: EO and 2-CE concentration and recovery results in rice samples

Rice sample		Conc Result (ng/g)		% Recovery	
		EO	2-CE	EO	2-CE
Rice Berry	Unspiked	Not detected	Not detected		
	10 ng/g spiked	8.362	10.213	84	102
	25 ng/g spiked	21.843	25.329	87	101
	100 ng/g spiked	87.747	77.651	88	78
Red Cargo Rice	Unspiked	Not detected	Not detected		
	10 ng/g spiked	8.799	11.094	88	111
	25 ng/g spiked	22.853	20.030	91	80
	100 ng/g spiked	87.766	103.334	88	103
Jasmine Rice	Unspiked	Not detected	Not detected		
	10 ng/g spiked	9.112	8.835	91	88
	25 ng/g spiked	22.392	28.905	90	116
	100 ng/g spiked	89.252	91.039	89	91
Brown Rice	Unspiked	Not detected	Not detected		
	10 ng/g spiked	8.680	10.026	87	100
	25 ng/g spiked	21.671	21.465	87	86
	100 ng/g spiked	91.773	108.369	92	108

**Figure 4:** MRM mass chromatograms of: (a) spiked EO (10 ng/g); (b) 2-CE (10 ng/g) in brown rice.

Conclusion

An internal standard method was successfully developed for the analysis of EO and 2-CE in rice using HS-20 NX Trap with GCMS-TQ8050 NX. The dynamic headspace technique and GC-MS/MS were employed, simplifying the sample preparation process and maximizing the sensitivity. Good linearity (R² value > 0.997) and recovery within 70 % to 130 % were achieved for all the samples analyzed in this experiment, even without the use of matrix-matched calibration curves.

Reference

1. Ting, E., *Ethylene Oxide and 2-Chloroethanol Analysis in Instant Noodles by Dynamic Headspace and GC-MS/MS*, Shimadzu Application News No. 04-AD-0299, accessed on 4 September 2024.
https://www.shimadzu.com/an/sites/shimadzu.com.an/files/pim/pim_document_file/applications/application_note/22545/an_04-ad-0299-en.pdf

GCMS-TQ is a trademark of Shimadzu Corporation or its affiliated companies in Japan and/or other countries



SHIMADZU

Shimadzu Corporation

www.shimadzu.com/an/

SHIMADZU (Asia Pacific) Pte. Ltd,
www.shimadzu.com.sg

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.
The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

04-AD-0308-EN

First Edition: Oct. 2024

➤ Please fill out the survey

Related Products

Some products may be updated to newer models.



➤ HS-20 NX series
Headspace Sampler



➤ GCMS-TQ™8050
NX
Triple Quadrupole Gas Chromatograph
Mass Spectrometry

Related Solutions

➤ Food and Beverages

➤ Food Contamination

➤ Price Inquiry

➤ Product Inquiry

➤ Technical Service /
Support Inquiry

➤ Other Inquiry