

Application News

MALDI-TOF mass spectrometry

Solvent-Free Preparation for MALDI-MS Analysis of Virgin and Recycled Polyethylene Terephthalate (PET)

Nina Timmer¹, Matt Openshaw², Andreas Baumeister¹
1 Shimadzu Europa GmbH, 2 Shimadzu Manchester

User Benefits

- ◆ Simple and quick sample preparation
- ◆ Solvent-free preparation for the analysis of poorly soluble polymers
- ◆ Avoidance of hazardous solvents

■ Introduction

Polyethylene terephthalate (PET) is a widely used polymer found in everyday applications, most notably in single-use beverage bottles (1). Since 2021, taxes have been imposed on non-recycled PET, increasing the demand for efficient methods to differentiate between virgin and recycled PET (2). This has made rapid and straightforward analytical techniques for PET more relevant than ever. However, PET presents several analytical challenges due to its poor solubility, complicating analysis using techniques like MALDI (Matrix-Assisted Laser Desorption/Ionization) (3). PET can be dissolved in solvents such as 1,1,1,3,3,3-hexafluoropropan-2-ol, but this solvent belongs to the highly controversial PFAS (per- and polyfluoroalkyl substances) class, known for its environmental persistence and potential health risks (4-5). The development of solvent-free analytical methods for PET is highly beneficial, as it eliminates the need for environmentally harmful solvents and aligns with growing sustainability standards.

■ Sample preparation

In the initial step was similar to (6). The PET sample is ground into smaller particles. Subsequently, the sample is weighed along with the matrix (dithranol) and the cation donor (sodium trifluoroacetate (NaTFA)), followed by thorough grinding with a mortar and pestle. After grinding, the sample is melted at 300°C and applied on a MALDI target. The prepared sample is then analysed on an AXIMA Performance™ instrument.

■ Measurement Conditions

The sample was analyzed on an AXIMA Performance MALDI-TOF mass spectrometer shown in Fig. 1. Table 1 shows the measurement conditions.

Table 1 Analysis conditions for AXIMA Performance

Mode	Reflectron (positive)
Mass range	1000 – 8000
Laser power	120
Profiles	50
Shots per profile	10
Ion Gate (Da)	1000



Fig. 1 AXIMA Performance™ MALDI TOF Mass Spectrometer

■ Results

The spectrum of virgin PET exhibits a series of consecutive peaks with a separation of 192 mass units, corresponding to the repeating monomeric unit of PET (Fig. 2). Additionally, an intense peak is observed in conjunction with a less intense series, separated by 44 mass units. This 44 mass unit difference suggests the presence of an additional carboxyl group within the polymer chain, indicating that the polymer chains possess a carboxyl-functional end group. In the spectrum of virgin PET are also present, suggesting that the primary chemical structure of the polymer is largely preserved during the recycling process. However, the spectrum of recycled PET reveals an additional variation in peak distribution, which indicates the presence of different terminal groups. These variations can be attributed to the recycling process, during which polymer chains may undergo scission, resulting in the formation of new end-group functionalities or alterations in the polymer chain structure, as evidenced by the altered peak distribution.




Fig. 2: Comparison of virgin PET spectrum (red) with recycled PET spectrum (blue)

■ Conclusion

Using solvent-free preparation methods, PET could be successfully analyzed in its solid phase. This approach is particularly advantageous because PET is difficult to dissolve, making solvent-free techniques highly effective. Furthermore, it was possible to distinguish between virgin PET and recycled PET, demonstrating the efficacy of this analytical method in differentiating between these materials.

■ References

- (1) Benyathiar, P. et al. (2022). Polyethylene Terephthalate (PET) Bottle-to-Bottle Recycling for the Beverage Industry: A Review. *Polymers*, 14(12), 2366. <https://doi.org/10.3390/polym14122366>
- (2) Directorate-General for Budget (European Commission). (2021). The EU's 2021-2027 Long-term Budget and NextGenerationEU: Facts and figures. ISBN: 978-92-76-30627-6. <https://doi.org/10.2761/808559>
- (3) Karim, S. S. et al. (2022). Model analysis on effect of temperature on the solubility of recycling of Polyethylene Terephthalate (PET) plastic. *Chemosphere*, 307, 136050. <https://doi.org/10.1016/j.chemosphere.2022.136050>
- (4) Jang, S. et al. (2012). Characterization of Thermal Degradation of Polytrimethylene Terephthalate by MALDI-TOF Mass Spectrometry. *Bulletin Of The Korean Chemical Society*, 33(3), 833–838. <https://doi.org/10.5012/bkcs.2012.33.3.833>
- (5) Abunada, Z., et al. (2020). An Overview of Per- and Polyfluoroalkyl Substances (PFAS) in the Environment: Source, Fate, Risk and Regulations. *Water*, 12(12), 3590. <https://doi.org/10.3390/w12123590>
- (6) Trimpin, S. et al. (2006). Solvent-Free MALDI-MS: Developmental Improvements in the Reliability and the Potential of MALDI in the Analysis of Synthetic Polymers and Giant Organic Molecules. *J. Am. Soc. Mass Spectrom*, 17, 661–671, <https://doi.org/10.1016/j.jasms.2006.01.007>

AXIMA Performance is a trademark of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.



SHIMADZU

Shimadzu Corporation
www.shimadzu.com/an/

SHIMADZU Europa GmbH,
www.shimadzu.eu

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

05-SCA-295-008-EN

First Edition: Dec. 2024

➤ Please fill out the survey

Related Products

Some products may be updated to newer models.



➤ AXIMA Performance

AXIMA Performance - a highly flexible research-grade...

Related Solutions

➤ Environment

➤ Engineering
Materials

➤ Plastics

➤ Price Inquiry

➤ Product Inquiry

➤ Technical Service /
Support Inquiry

➤ Other Inquiry