

Application News

Energy Dispersive X-ray Fluorescence Spectrometer ALTRACE™

Screening Analysis for Hazardous Heavy Metals in Foods and Food Additives

Kotomi Kuroda, Hirokazu Moriya

User Benefits

- ◆ In comparison with the conventional technology, sensitivity was improved by high voltage X-ray tube and optimization of the optical design.
- ◆ Analysis throughput is improved because continuous analysis of a maximum of 48 samples is possible.
- ◆ Screening for hazardous elements in foods and food additives is possible with only simple sample preparation, such as loading the samples into the sample container.

Introduction

Control of hazardous heavy metals contained in foods and food additives is required in order to protect human health and safety. Therefore, various domestic laws and regulations ⁽¹⁾⁻⁽⁴⁾ have been established in Japan, including "Standards and criteria for food and food additives, etc." and "Japan's Specifications and Standards for Food Additives".

From the viewpoint of sensitivity, the main technique used in analyses of toxic heavy metals is atomic absorption spectrophotometry. However, as disadvantages of this method, powder and solid samples must be dissolved in an acid such as nitric acid or hydrochloric acid, and advanced technology and know-how are required in the actual analysis. In contrast, fluorescent X-ray spectroscopy offers excellent convenience because analysis is possible as-is, without dissolving the specimens.

Shimadzu's new product, the ALTRACE, energy dispersive X-ray fluorescence spectrometer (Fig. 1), makes it possible to analyze toxic heavy metals with high sensitivity thanks to the increased output of the X-ray tube and optimization of the optical design. As an additional advantage, continuous analysis of a maximum of 48 samples is possible, contributing to improved analysis throughput.

This Application News article introduces the following:

1. Analysis of cadmium (Cd) in rice
2. Analysis of arsenic (As) and lead (Pb) in food additives



Fig. 1 ALTRACE™

Calibration Curves

Calibration curves were prepared for cadmium (Cd), arsenic (As), and lead (Pb) using standard solutions for atomic absorption spectrophotometry adjusted to the four concentration levels shown below. Fig. 2 shows the calibration curves.

Here, correction of the sample shape and matrix was done with scattered X-ray. In addition, correction for overlap of As by Pb was also carried out (coexisting element correction, dj method).

Concentrations of calibration curve standard solutions:
0 (blank), 0.5, 1, 5 µg/g

Lower Limit of Detection

Table 1 shows the lower limits of detection of each element calculated from the theoretical statistical variation of the calibration curves.

Table 1 Lower Limits of Detection [µg/g]

Element	Cd	As	Pb
Analytical line	CdKα	AsKα	PbLβ1
Lower limit of detection	0.08	0.08	0.10




Fig. 2 Calibration Curves

1. Analysis of Cadmium in Rice

■ Sample

White rice powder, certified reference material NMIJ CRM 7502-a

■ Sample Preparation

The sample was set in a sample cell lined with polypropylene film (thickness: 5 μm) and subjected to simple compression. Fig. 3 shows an image of the sample.



Fig. 3 Image of Sample

■ Cd Profile

Fig. 4 shows the profile of CdK α . A peak can be observed clearly, even with a trace amount of cadmium (0.548 $\mu\text{g/g}$).




Fig. 4 CdK α Profile

■ Results of Repeatability Test

Table 2 shows the results of a 10-time repeatability test.

Element	Cd
Standard value	0.548
Average value	0.532
Standard deviation	0.050
Coefficient of variation [%]	9.4

2. Analysis of Arsenic and Lead in Food Additive

■ Sample

- Adipic acid
- L(+)-ascorbic acid (vitamin C)
- Citric acid

■ Sample Preparation

The sample materials were set in a sample cell lined with polypropylene film (thickness: 5 μm) and subjected to simple compression. Due to the large grain shape of the citric acid, the sample material was crushed with a crushing machine.

■ Results of Quantitative Analysis

Table 3 shows the results of the quantitative analysis. Both arsenic and lead were below the lower limit of quantitation. Because the lower limit of quantitation is much lower than 1 $\mu\text{g/g}$, it can be understood that this analysis is effective for screening analysis of the organic samples listed in Japan's Specifications and Standards for Food Additives.

Element	As	Pb
Adipic acid	<0.21	<0.19
L(+)-ascorbic acid (vitamin C)	<0.28	<0.17
Citric acid	<0.28	<0.25

* The symbol "<" indicates that the result is below the lower limit of quantitation. The lower limit of quantitation is the theoretical 10 σ value.

■ Analysis Conditions

Instrument	: ALTRACE
Elements	: Cd, As, Pb
Analysis group	: Quantitative analysis
Analysis method	: Calibration curve method
Detector	: SDD
X-ray tube	: Rh target
Tube voltage	: 50 [kV] (As, Pb) 65 [kV] (Cd)
Tube current	: Auto [μA]
Primary filter	: #1 (Cd), #5 (As, Pb)
Atmosphere	: Air
Accumulation time	: 300 [s]
Dead time	: Max. 40 %

■ Conclusion

The analytical sensitivity of the new Shimadzu ALTRACE was greatly improved by increasing the power of the X-ray tube and optimization of the optical design. As a result of this improved performance, ALTRACE is the optimum instrument for screening for hazardous elements in foods and food additives. In addition, because a large number of samples can be analyzed in a short time, ALTRACE also contributes to improved analytical throughput for the customer.

<References>

- (1) Standards and criteria for food and food additives, etc. (1959, Ministry of Health and Welfare notification No. 370)
- (2) Public Health Bureau, Food Safety Standards and Evaluation Division notification No.0308-1, dated March 8, 2024
- (3) Japan's Specifications and Standards for Food Additives, 10th Edition
- (4) CODEX GENERAL STANDARD FOR CONTAMINANTS AND TOXINS IN FOOD AND FEED (CODEX STAN 193-1995)

ALTRACE is a trademark of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.



SHIMADZU

Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

01-00775-EN

First Edition: Feb. 2025

➤ Please fill out the survey

Related Products

Some products may be updated to newer models.



➤ ALTRACE

Energy Dispersive X-ray Fluorescence Spectrometer

Related Solutions

➤ Food and Beverages

➤ Food Additive

➤ Food Contact Materials (FCM)

➤ Price Inquiry

➤ Product Inquiry

➤ Technical Service /
Support Inquiry

➤ Other Inquiry