

Application News

DPIMS™ QT Probe Electrospray Ionization Kit
iMScope™ Imaging Mass Microscope

Distribution Analysis of Plant Alkaloids in *Narcissus*, by MS Imaging

Kaoru Nakagawa¹, Tetsuo Iida¹, Shuichi Shimma², Manami Kobayashi¹
1 Shimadzu Corporation, 2 Osaka University

User Benefits

- ◆ Qualitative screening using the DPIMS QT kit can be conducted with measurement time of approximately 0.5 min.
- ◆ Easily switch between DPIMS QT and iMScope QT.
- ◆ Distribution analysis using iMScope QT is expected to improve the efficiency of drug production.

■ Introduction

Plant alkaloids are used in various pharmaceuticals, such as anticancer drugs and analgesics. Among these plant alkaloids, galanthamine is an Amaryllidaceae-type alkaloid with acetylcholinesterase inhibitors used in the treatment of neurological diseases such as Alzheimer's disease. Although the chemical synthesis of galanthamine has been successfully achieved, *Narcissus* is the main source of its production. Research indicates that galanthamine content varies not only with the type of *Narcissus*, but also with the developmental stage and the part of the plant¹⁾. Pharmaceutical companies are pursuing plant species with higher galanthamine content to increase pharmaceutical productivity.

In this study, we were able to quickly confirm the presence of galanthamine in our *Narcissus* sample (Fig. 1) using the DPIMS QT probe electrospray ionization (PESI) kit and the LCMS quadrupole time-of-flight (Q-TOF) mass spectrometer (MS) (Fig. 2). Subsequently, we analyzed the distribution of galanthamine by MS imaging (MSI) using the iMScope QT and the LCMS Q-TOF MS (Fig. 2).

MSI analysis was performed using the iMScope QT atmospheric MALDI equipped with an optical microscope and the LCMS-9030 (Fig. 2 left). Analytical settings of MSI are shown in Table 2. Data analysis was performed using the IMAGEREVEAL™ MS.

Table 1 Analytical Settings of PESI


Mass spectrometer	
System	: DPIMS QT+LCMS-9030
Polarity	: Positive
DL temp	: 250 °C
Heat block temp	: 50 °C
Interface Voltage	: 3.5 kV
MS Range	: <i>m/z</i> 50-2,000
Measurement Time	: 0.5 min

Table 2 Analytical Settings of MSI


Mass spectrometer	
System	: iMScope QT+LCMS-9030
Polarity	: Positive
DL temp	: 250 °C
Heat block temp	: 450 °C
MS Range	: <i>m/z</i> 280-335
Spatial Resolution (Pitch)	: 10 / 25 μm
Laser Diameter Setting	: 1 / 2
Laser Intensity	: 50 / 60
Laser Repetition Frequency	: 20 kHz
Matrix Coating	
System	: iMLayer
Matrix Used	: CHCA
Coating Method	: Deposition with 0.7 μm Thickness




Narcissus tazetta



Galanthamine
C₁₇H₂₁NO₃



Lycorine
C₁₈H₁₇NO₄



Tazettine
C₁₈H₂₁NO₄

Fig. 1 Structural Formula of Alkaloids in *Narcissus tazetta*

■ Sample Preparation and Analysis Conditions

For our sample, we used *Narcissus tazetta* with leaves growing to about 15 cm in length. Galanthamine was extracted with 50% EtOH aq. from freeze-dissolved leaves. The solution was analyzed via the DPIMS QT and the LCMS-9030 (Fig. 2 right). Analytical settings of PESI are shown in Table 1. Data analysis was performed using the LabSolutions Insight Explore™ software. Locations of *Narcissus tazetta* sections are shown in Fig. 3. Frozen leaves and bulbs were sliced to a thickness of 20 μm using a microtome and mounted on indium tin oxide (ITO) coated glass slides. These were coated with α-cyano-4-hydroxycinnamic acid (CHCA) via vapor deposition by using the matrix sublimation apparatus iMLayer™ at a thickness of 0.7 μm.



Fig. 2 iMScope™ QT (left) and DPIMS™ QT (right)




Fig. 3 Locations of *Narcissus tazetta* sections for MS imaging
Narcissus tazetta were grown in soil culture with leaves growing to about 15 cm. Three types of sections were prepared: bulb, root leaf, and leaf cross section.

■ Qualitative Screening of Galanthamine in *Narcissus Tazetta* by DPiMS QT

The mass spectra obtained quickly in only a few tens of seconds from the extracted solution of freeze-crushed *Narcissus tazetta* leaves in positive mode are shown in Fig. 4.

Accurate mass analysis using analysis software confirmed the presence of galanthamine, other plant alkaloids like lycorine and tazettine, the metabolite choline, and sugars composed of hexose were also detected. Plant alkaloids were detected within a mass error of 1 mDa (Table 3).




Fig. 4 Mass Spectra of *Narcissus Tazetta* Leaves

Table 3 Mass Accuracy of Alkaloids Detected in *Narcissus Tazetta*

Compounds	Ions	Theoretical <i>m/z</i>	Measured <i>m/z</i>	Mass Error (mDa)
Galanthamine	[M+H] ⁺	288.1594	288.1600	0.6
Lycorine	[M+H] ⁺	288.1230	288.1235	0.5
Tazettine	[M+H] ⁺	322.1493	322.1496	0.3

■ Distribution Analysis of Galanthamine in *Narcissus Tazetta* by MSI

The distribution of galanthamine, lycorine, and tazettine in sections prepared from the bulb and leaves (two locations) of the *Narcissus tazetta* shown in Fig. 3 was analyzed by MSI. The MS imaging results of a bulb section at 25 μm spatial resolution and 10 μm spatial resolution are shown in Fig. 5, these images were taken with the 5x objective lens. These results showed strong distributions of galanthamine, lycorine, and tazettine from all samples. The MS images of the bulb showed that galanthamine was distributed mostly in the section of the bulb that later grows and becomes the leaves, whereas MS images of lycorine and tazettine did not show such distribution, suggesting that the distribution area of plant alkaloids in the bulb differs from species to species, but no such differences in distribution were observed in the leaf cross sections.

■ Conclusion

The combination of PESI and Q-TOF was able to detect galanthamine in *Narcissus* with a high mass accuracy of within 1 mDa. The time required for analysis was significantly reduced compared to that by LC or LC/MS with a measurement time of 0.5 min, suggesting that rapid screening is possible. The specific distribution of galanthamine in narcissus bulbs was confirmed by MS imaging analysis using a combination of iMScope QT and Q-TOF, with the ionization unit replaced from PESI to MALDI. Such distribution analysis enables identification of regions with high galanthamine content, facilitating efficient isolation and purification in the pharmaceutical manufacturing process, and is expected to contribute to the reduction of extraction costs.




Fig. 5 Mass Images of Alkaloids in *Narcissus Tazetta*

This application news is made and published by rearranging the data from *Mass Spectrom.* 2024;13(1):A0163. For details, please refer to the original paper.

<References>

- 1) Letter... *Planta Med* 2003; 69: 1166-1168

iMScope, iMLayer, IMAGEREVEAL, LabSolutions Insight Explore, and DPiMS are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.



SHIMADZU

Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

01-00810-EN

First Edition: March 2025

➤ Please fill out the survey

Related Products

Some products may be updated to newer models.



➤ iMScope QT
Imaging Mass Microscope



➤ DPiMST™ QT
Kit for Direct Probe Ionization Mass Spectrometer

Related Solutions

➤ Pharmaceutical and Biopharmaceutical

➤ Small Molecule Pharmaceutical

➤ Life Science

➤ Dementia / Alzheimer's disease

➤ Price Inquiry

➤ Product Inquiry

➤ Technical Service / Support Inquiry

➤ Other Inquiry