

Application News

GC-MS GCMS-QP™2050

Analysis of Formic Acid in Acetone Using GC/MS System

Akara Hashimoto¹ and Masahiro Kawazoe²

1 Shimadzu Corporation, 2 Shimadzu Techno-Research, Inc.

User Benefits

- ◆ The GCMS-QP2050 enables highly accurate measurement of formic acid in solution.
- ◆ The SIM mode enables the detection of formic acid at sub-ppm levels.
- ◆ The Scan mode enables the qualitative analysis of unknown components other than formic acid.

■ Introduction

In research on artificial photosynthesis, the analysis of impurities in chemical products and raw materials has increased the demand for high-sensitivity analysis of formic acid. Previous examples of formic acid analysis using GC-BID and GC-FID (Jetanizer™) have been presented, but in this article high-sensitivity analysis of formic acid using the GCMS-QP2050 (see Fig. 1) is introduced.



Fig. 1 GCMS™-QP2050

■ Insert, Column Phosphoric Acid Treatment

Formic acid is known to adsorb onto the injection port and column, making peak detection difficult. However, by performing phosphoric acid treatment on the glass insert and column before measurement, good peaks can be obtained. The procedure for phosphoric acid treatment of the glass insert is shown in Fig. 2, and the procedure for phosphoric acid treatment of the column is shown in Fig. 3.

1. Prepare the split insert (P/N 227-35007-01).
↓
2. Prepare a 0.3 % phosphoric acid acetone solution.
↓
3. Immerse the split insert of step 1 in the prepared solution of step 2 for more than 1 minute (see diagram on right).
↓
4. Dry at 50 °C for 1 hour.



Fig. 2 Phosphoric Acid Treatment Procedure for Glass Insert

1. Prepare a 100 ppm phosphoric acid methanol solution.
↓
2. Measure four times under the analytical conditions shown in Table 1.
↓
3. Measure methanol ten times at 150 °C in the column (with all other conditions the same as in Table 1) to stabilize the column interior.

Fig. 3 Phosphoric Acid Treatment Procedure for Column

■ Analytical Conditions

Table 1 shows the analytical conditions.

Table1 Analytical Conditions	
GC-MS Model:	GCMS-QP2050 (GC-2030)
Autoinjector:	AOC-30i
[GC]	
Column:	SH-PolarWax (0.32 mm I.D. × 30 m, d.f. = 1.0 μm) P/N 227-36252-01
Column Temp.:	80 °C – 5 °C/min – 180 °C (3 min) (Total 23 min)
Injection Temp.:	180 °C
Purge:	10 mL/min
Injection Mode:	Split
Split Ratio:	2
Carrier Gas Controller:	Constant Linear Velocity (He)
Linear Velocity:	60 cm/sec
Injection Volume:	1 μL
[MS]	
Pumping Speed:	255 L/s
Acquisition Mode:	Scan, SIM
Ion Source Temp.:	200 °C
Interface Temp.:	200 °C
SIM monitoring m/z:	Target: 46, For confirmation: 29
m/z range:	20 - 300
Event time:	0.3 sec

■ Analysis of Formic Acid in Acetone Using SIM Mode

The SIM chromatogram for formic acid in acetone at 0.2 ppm is shown in Fig. 4, the calibration curve for 0.2 ppm to 20 ppm is shown in Fig. 5, and the repeatability of the analysis at 0.2 ppm (n = 5) is presented in Table 2.




Fig. 4 SIM Chromatogram of 0.2 ppm Formic Acid




Fig. 5 Calibration Curve from 0.2 ppm to 20 ppm

Table 2 Repeatability of Analysis at 0.2 ppm (n = 5)

Data	Area
Data 1	40,394
Data 2	40,833
Data 3	40,733
Data 4	41,084
Data 5	42,828
Area Average	41,174
%RSD	2.3

■ Analysis of Formic Acid Using Scan Mode

The Total Ion Chromatogram (TIC) chromatogram of a 20 ppm formic acid standard solution is shown in Fig. 6.




Fig. 6 TIC Chromatogram of a 20 ppm Formic Acid Standard Solution

A peak labeled “a” other than formic acid was observed. By conducting a library search using the mass spectrum of peak a, acetic acid was identified. Thus, by performing analysis in Scan mode, it is possible to identify unknown components.




Fig. 7 Library Search Results for Peak a

■ Comparison of Formic Acid Analyses

A comparison of formic acid analyses using GC-BID, GC-FID (Jetanizer), and GC-MS is shown in Table 3.

Table 3 Comparison of Formic Acid Analysis Across Different Analytical Methods

	GC-BID	GC-Jetanizer	GC-MS*2
Phosphoric acid treatment	Needed	Needed	Needed
Conc. Range*1	Approximately 1 ppm to several 1000 ppm	Approximately 1 ppm to several percent	Approximately 0.1 ppm to several ppm
Carrier Gas	He only	He, N ₂ , H ₂	He recommended

*1 The concentration range varies depending on analytical conditions, samples, and the environment.

*2 The values for GC-MS are for the SIM mode.

■ Conclusion

By performing phosphoric acid treatment on the insert and column, it became possible to analyze formic acid in SIM mode down to 0.2 ppm. Furthermore, good results could be obtained for both linearity and repeatability. Additionally, by analyzing in Scan mode, identification of unknown peaks other than formic acid was also achieved.

<Related Applications>

- Improvement of Sensitivity and Repeatability in Analysis of Formic Acid [Application News No. G279](#)
- High-Sensitivity Analysis of Formic Acid in Methanol Solution Using Jetanizer™ [Application News 01-00833-en](#)

GCMS-QP and Jetanizer are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <https://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol “TM” or “®”.

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you “as is” without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

› Please fill out the survey

Related Products

Some products may be updated to newer models.



› GCMS-QP2050

Gas Chromatograph Mass Spectrometer

Related Solutions

› Chemicals

› Food and Beverages

› Environment

› Plastics

› Price Inquiry

› Product Inquiry

› Technical Service /
Support Inquiry

› Other Inquiry