

Application News

Non-Suppressor Ion Chromatograph

Analysis of Ammonia as an Impurity in Fuel Cell Grade Hydrogen According to ISO 14687 Grade D Using a Gas Generator

Hiroyuki Aikawa, Yujing Jiang, Ayano Tanabe

User Benefits

- ◆ The non-suppressor method provides good linearity of ammonium ions and reduces maintenance costs.
- ◆ Since organic solvents are not used for extraction in pretreatment, environmental impact can be reduced.
- ◆ Using a gas generator ensures the validity of the method including pretreatment.

Introduction

Hydrogen is expected to be a new power and heat source to reduce greenhouse gas emissions. Hydrogen is an easily available substance that can be produced from the electrolysis of water in our daily lives. However, when used in fuel cells, it is necessary to control impurities in hydrogen gas because certain impurities introduced during the manufacturing process deteriorate the electrode and electrolyte membrane, reducing the performance of the fuel cell or shortening its life. When ammonia is used as a carrier of hydrogen gas due to its distribution as a chemical, there is a possibility that ammonia may be introduced into hydrogen. This article presents an example of analysis of ammonia in hydrogen using an ion chromatograph (IC) based on the standard values described in ISO 14687 Grade D.¹⁾

Hydrogen Impurity Standard under ISO 14687 Grade D

ISO 14687 Grade D specifies impurities in hydrogen that affect the quality of fuel cell vehicles. Ammonia is specified to be 0.1 μmol/mol or less, and the concentration converted to ammonium ion (NH₄) by IC is 26 μg/L.

Analysis of Ammonium Ion by IC

A non-suppressor IC (Fig. 1) was used to analyze NH₄. Since the IC can selectively measure ions dissolved in the sample, it is less affected by contaminants than other methods and can measure NH₄ with high accuracy. In addition, by using the same device together with a suppressor IC, it is possible to measure inorganic halogens and formic acid, which are other impurities listed in ISO 14687. Table 1 shows the analytical conditions for NH₄, Fig. 2 shows the chromatogram of a standard sample at a concentration near the specified value (25 μg/L), Fig. 3 and Table 2 shows the calibration curve in the range from 2.5 μg/L to 200 μg/L. As shown in Fig. 3, the linearity was greater than 0.999. Table 3 shows the repeatability of the sample at 10 μg/L, which is around 1/2 of the specified concentration, repeated 7 times.



Fig. 1 Non-suppression ion chromatograph HIC-NS

Table 1 Analytical conditions

System	: HIC-NS
Column	: Shim-pack™ IC-C4 (150 mm × 4.6 mm I.D., 7 μm) ^{*1}
Guard Column	: Shim-pack IC-GC4 (10 mm × 4.6 mm I.D., 7 μm) ^{*2}
Mobile Phase	: 2.5 mM Methane sulfonic acid
Flow Rate	: 1.0 mL/min
Column Temp.	: 40 °C
Injection Vol.	: 100 μL
Vial	: 4 mL, Polypropylene ^{*3}
Detection	: Conductivity

*1 P/N : 228-41616-91, *2 P/N : 228-59900-91, *3 P/N : 228-31537-91




Fig. 2 Chromatogram of standard sample (25 μg/L= equivalent to 0.1 μmol/mol)




Fig. 3 Calibration curve

Table 2 Concentration of Calibration Curve and Peak Area

	Concentration of Calibration Curve (μg/L)	Peak Area
STD1	2.5	481
STD2	10	1135
STD3	25	3281
STD4	50	6328
STD5	100	12272
STD6	200	25711
Coefficient of Determination(r ²)	0.9993	

Table 3 Repeatability (10 µg/L n=7)

	Retention time (min)	Peak area	Concentration (µg/L)
1	4.36	1315	9.99
2	4.36	1355	10.3
3	4.37	1380	10.5
4	4.37	1284	9.74
5	4.37	1292	9.80
6	4.37	1246	9.43
7	4.36	1318	10.0
Average	4.36	1313	9.97
%RSD	0.04	3.41	3.61

Verification of Concentration Using a Gas Generator

Next, using a gas generator 491 Flex², we mixed ammonia standard gas near the standard value into hydrogen and measured NH₄ collected in ultrapure water by IC. The gas collection conditions and flow are shown in Fig. 4. Fig. 5 shows the chromatograms of ultrapure water collected from 1st and 2nd impinger, and Table 4 shows the quantitative values and recovery rates of 3 analyses.

When the ammonia concentration was 0.1 µmol/mol, which is the standard value, the concentration of ammonia collected by 2nd impinger was N/D, indicating that ammonia was properly collected by 2nd impinger. Fig. 6 shows the sum of the concentrations of the first and second impingers and the correlation between the quantitative value of IC and the concentration of gas produced by the gas generator.




Fig. 4 Gas collection conditions and flow




Fig. 5 Chromatograms of NH₄ discharged from gas generator (0.1 µmol/mol)

Table 4 At ammonia standard gas concentration of 0.1 µmol/mol Analysis results of NH₄ quantitative value and recovery rate (3 times)

	NH4 Assay value (µmol/mol)		Total NH ₄ Assay value (µmol/mol)	Recovery (%)
	1st Impinger	2nd Impinger		
1st	0.091	N/D *	0.091	91
2nd	0.092	N/D	0.092	92
3rd	0.084	N/D	0.084	84

* S/N=3 or less is defined as N/D.




Fig. 6 Correlation between quantitative value of IC and gas concentration produced by gas generator

Conclusion

Ammonia, an impurity in hydrogen, was analyzed using a non-suppressor IC, and it was confirmed that the method met the standard value specified by ISO 14687 Grade D. In addition, when ammonia standard gas near the standard value was mixed into hydrogen using a gas generator, a good recovery rate was obtained, which confirmed the validity of the method including the collection flow of pretreatment.

< References >

- 1) ISO 14687 Hydrogen fuel quality-product specifications (2025)
- 2) A product of KIN-TEK in the United States. <https://kin-tek.com/>

<Related Applications>

1. Analysis of Inorganic Halogen, Formic Acid, Ammonia, and Formaldehyde as Impurities in Fuel Cell Grade Hydrogen According to ISO14687 Grade D, [Application News No. 01-01071](#)

Shim-pack is a trademark of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <https://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

01-01070A-EN

First Edition: Dec. 2025
Revision A: Feb. 2026

› Please fill out the survey

Related Products

Some products may be updated to newer models.



› Anion Suppressor
Ion Chromatograph

Related Solutions

› Clean Energy

› Fuel Cells

› Price Inquiry

› Product Inquiry

› Technical Service /
Support Inquiry

› Other Inquiry