

# Application News

## Confirmation of Intact Protein Molecular Weight Using a Single Quadrupole Mass Spectrometer

Miho Kawashima, Noriko Kato

### User Benefits

- ◆ The molecular weight of intact proteins can be confirmed by using the simple handling LCMS-2050 single quadrupole mass spectrometer.
- ◆ By using LabSolutions Insight™ Biologics, the results of multiply-charged ion analyses can be understood easily by a simple visual comparison.

### Introduction

Biopharmaceuticals are a class of pharmaceuticals produced using biotechnologies such as genetic modification and cell culture. The proteins that serve as the active components of biopharmaceuticals are macromolecular compounds with multiple ionization sites. When multiply-charged ions are detected, their molecular weights can be estimated by analysis processing (deconvolution) on an appropriate software program. More accurate estimation of molecular weight is possible by using a high-resolution mass spectrometer. However, when the main components can be assumed in advance, for example, in quality control applications, the molecular weight can be confirmed by using a simple handling single quadrupole mass spectrometer.

This article introduces an example of analysis of proteins using a single quadrupole mass spectrometer for synthesis confirmation of the macromolecular compounds.



Fig. 1 Nexera™, LCMS-2050 System

### Analysis Conditions

The measuring instruments used in this experiment were a Nexera X3 high-performance liquid chromatograph and an LCMS-2050 single quadrupole mass spectrometer. Table 2 and Table 3 show the analysis conditions used with the two instruments, respectively.

Table 2 Liquid Chromatograph (LC) Analysis Conditions

System	: Nexera X3
Column	: Shim-pack Scepter™ C4-300 *1 (50 mm × 2.1 mm I.D., 3.0 μm)
Temperature	: 40 °C
Injection volume	: 5 μL
Mobile phases	: 0.1 % formic acid in Water 0.1 % formic acid in Acetonitrile
Flow rate	: 0.4 ml/min
Time program (%B)	: 5 % (0.00 min)→60 % (5.00-15.00 min)→ 5 % (15.01-18.00 min)

\*1 P/N: 227-31177-03

Table 3 Mass Spectrometer (MS) Analysis Conditions

System	: LCMS-2050
Ionization	: ESI/APCI (DUIS), Positive mode
Interface voltage	: 3.0 kV
Mode	: Scan <i>m/z</i> 500 - 2000
Nebulizing gas flow	: 2 L/min
Drying gas flow	: 5 L/min
Heating gas flow	: 7 L/min
Desolvation temp.	: 450 °C
DL temp.	: 200 °C

### Samples

Table 1 shows the proteins that were measured as model samples. These commercially-available standard samples were dissolved and diluted with ultrapure water for use as the analysis samples.

Table 1 List of Measured Compounds

#	Compounds
1	Myoglobin from equine skeletal muscle
2	Carbonic anhydrase from bovine red blood cells
3	β-casein from bovine milk

### Deconvolution Results

Fig. 2 and Fig. 3 show the mass spectra of the detected peaks and the multiply-charged ion analysis mass spectra obtained by deconvolution of the myoglobin and bovine carbonic anhydrase samples, respectively. LabSolutions Insight Biologics includes the multiply-charged ion analysis algorithm "ReSpect" (Positive Probability Limited), which is useful when estimating the molecular weights of proteins with large molecular weights by multiply-charged ion analysis.




Fig. 2 (A) Mass Spectrum and (B) Multiply-Charged Ion Analysis Mass Spectrum of Myoglobin




Fig. 3 (A) Mass Spectrum and (B) Multiply-Charged Ion Analysis Mass Spectrum of Bovine Carbonic Anhydrase




Fig. 4 Results of Multiply-Charged Ion Analysis (A) Overlay of Component Chromatograms, (B) Multiply-Charged Ion Analysis Mass Spectrum, (C) Mass Spectrum

Fig. 4 shows the analysis screen of  $\beta$ -casein by Insight Explore Biologics. The peaks detected from the chromatogram, the mass spectrum of the obtained peaks, and the mass spectrum after multiply-charged ion analysis can be observed on a single screen.

LabSolutions Insight Biologics can generate component chromatograms combining all different valences and isotope spectra. Since it is possible to overlay the component chromatograms and display the molecular weights calculated by deconvolution as peak tops, the analysis results can be understood easily by visual comparison when multiple components are assumed to exist or impurities are detected.

For example, in this experiment, two peaks were detected from  $\beta$ -casein, and from the overlay of the component chromatograms, it can be understood that the estimated molecular weights of the respective peaks are 24088 and 23979.

### Conclusion

Simple conformation of the molecular weights of three intact protein samples was possible by using the LCMS-2050 single quadrupole mass spectrometer and LabSolutions Insight Biologics analysis software.

By using LabSolutions Insight Biologics, analysis results can be understood easily by visual comparison when multiple components are assumed to exist or impurities are detected.

<Related Application News article>

1. Analysis of Oligonucleotide Impurities Using Single Quadrupole Mass Spectrometer, [Application News No. 01-00974-EN](#)

Nexera, LabSolutions Insight, DUIS, and Shim-pack Scepter are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <https://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

01-01023A-EN

First Edition: Jan. 2026  
Revision A: Mar. 2026

› Please fill out the survey

## Related Products

Some products may be updated to newer models.



### › LCMS-2050

High-Performance Liquid  
Chromatograph Mass Spectro...



### › LabSolutions Insight Biologics

Software for Oligonucleotide Sequence  
Characteriza...

## Related Solutions

› Pharmaceutical and  
Biopharmaceutical

› Biopharmaceutical

› Price Inquiry

› Product Inquiry

› Technical Service /  
Support Inquiry

› Other Inquiry