

## Application News

GC Nexis™ GC-2030

# Detailed Hydrocarbon Analysis by Nexis GC-2030 Using ASTM D6730

Kei Usui, Chisato Kanamori

### User Benefits

- ◆ Detailed compositional analysis of gasoline in accordance with ASTM D6730 is feasible using the Nexis GC-2030.
- ◆ By employing PONAsolution™, a PONA/DHA analysis software, laborious identification tasks can be reduced, and quantitation calculations and report generation can be performed automatically.

### ■ Introduction

ASTM D6730 specifies the analysis of hydrocarbons in gasoline and oxygenates such as ethanol, methyl tert-butyl ether (MTBE), and ethyl tert-butyl ether (ETBE) over a boiling point range up to 225 °C. Identification and quantification of individual compounds in gasoline are critically important for refinery process control and product quality assurance.

In this application news, we report compositional analysis of gasoline conforming to ASTM D6730 using the Nexis GC-2030 gas chromatograph. Gasoline typically comprises over 100 hydrocarbons, making identification tasks highly complex; however, PONAsolution, a PONA/DHA analysis software, simplifies this process.

### ■ System Configuration and Analysis Conditions

The instrument configuration and analytical conditions are listed in Tables 1 and 2, respectively. ASTM D6730 requires separation of more than 100 components; therefore, a slightly polar pre-column (SH PONA Tuning Column, P/N: 227-36339-01) and a nonpolar 100 m main column (SH-1 PONA, P/N: 221-76196-00), together with CRG, were employed. The inlet pressure was adjusted so that the retention time of methane at 35 °C fell within  $7.00 \pm 0.02$  min. Data processing was performed using PONAsolution Ver.6, which contains a retention index library tailored for ASTM D6730.

In contrast, gasoline compositional analysis methods defined in ASTM D8071 and D8369 employ vacuum ultraviolet (VUV) detectors. In those methods, identification is performed using retention indices together with absorption spectral information, enabling deconvolution of coeluting components through spectral matching (see [Application News No. 01-00966-EN](#)). Table 3 compares PONA/DHA analyses using FID and VUV detectors.

Table 1 System configuration

|                |                                                                                                                                                                                             |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| GC Model       | : Nexis GC-2030 / AOC-30i                                                                                                                                                                   |
| Injection Port | : SPL                                                                                                                                                                                       |
| Liner          | : Tapered Split/Splitless Liner with fixed wool - "into wool" (P/N : 221-75191)                                                                                                             |
| Column         | : SH PONA Tuning Column*<br>(P/N : 227-36339-01)<br>(5.0 m × 0.25 mm I.D. × 1.0 μm)<br>*cut the column into 3.2 m<br>+ SH-1 PONA<br>(P/N : 221-76196-00)<br>(100 m × 0.25 mm I.D. × 0.5 μm) |
| Detector       | : FID-2030                                                                                                                                                                                  |
| Options        | : CRG-2030 (CO <sub>2</sub> )                                                                                                                                                               |
| Software       | : LabSolutions, PONAsolution                                                                                                                                                                |

Table 2 Analysis conditions

|                           |                                                                               |
|---------------------------|-------------------------------------------------------------------------------|
| Injection Temperature     | : 250 °C                                                                      |
| Flow Control Mode         | : Pressure (He)                                                               |
| Inlet pressure            | : 293.0 kPa                                                                   |
| Purge Flow                | : 3.0 mL/min                                                                  |
| Injection Volume          | : 0.5 μL                                                                      |
| Split ratio               | : 150                                                                         |
| Column Oven Temp. Program | : 5 °C (10 min) → 5.0 °C/min → 48 °C (59 min) → 1.3 °C/min → 200 °C (5.0 min) |
| Detector Temperature      | : 250 °C                                                                      |
| Makeup Gas                | : N <sub>2</sub> 24 mL/min                                                    |
| Detector Gas              | : H <sub>2</sub> 32 mL/min, Air 150 mL/min                                    |

### ■ PONAsolution

PONAsolution Ver.6 provides retention index libraries for each relevant standard (ASTM D6730, D6729, D5134, D6733). Identification of n-alkanes automatically corrects the retention times of the remaining several hundred components, thus streamlining the identification process.

When only a specific retention region requires adjustment, the two-point identification function is useful (Fig. 1). In the two-point identification procedure, reference and sample chromatograms are compared while assigning corresponding peaks. First, click on an identified component mark in the reference data (Fig. 1, a), then click the matching sample peak (Fig. 1, b). Repeat this for a second pair of peaks. The retention times of components between these two anchor points are then automatically corrected.




Fig. 1 Two-point identification using PONAsolution

Table 3 Comparison of ASTM D8369, D8071 and D6730

|                     | D6730                                 | D8369                                 | D8071                                            |
|---------------------|---------------------------------------|---------------------------------------|--------------------------------------------------|
| Report Items        | PIONA group<br>Individual compounds   | PIONA group<br>Individual compounds   | PIONA group<br>select compounds (EtOH, BTX etc.) |
| Detector            | FID                                   | VUV                                   | VUV                                              |
| Analysis Time (min) | ca. 180                               | 49                                    | 35                                               |
| Qualification       | Retention Indices                     | Retention Indices + VUV spectra       | Retention Indices + VUV spectra                  |
| Quantification      | Corrected Percentage Peak Area Method | Corrected Percentage Peak Area Method | Corrected Percentage Peak Area Method            |

## ■ Peak Resolution


Oxy setup blend (RESTEK: 33034) was analyzed to confirm peak shape and resolution. The chromatogram is shown in Fig. 2. Using the ASTM D6730 method with a slightly polar pre-column (SH PONA Tuning Column) connected to a nonpolar main column (SH-1 PONA) allowed separation of pairs that are difficult to resolve on a single column, for example tert-Butanol (3) and 2-Methylbutene-2 (4), or 1-Methylcyclopentene (8) and Benzene (9).

## ■ Quantification Results

To assess quantitative accuracy, a certified component standard (Tokyo Chemical Industry Co., Ltd., S0429) was analyzed. For all analytes evaluated using ASTM D6730, the difference between measured values and certified values was within 0.5%, indicating excellent agreement.

Table 4 Quantification results for reference material of gasoline components

| Compound | a. Conc. of certificate | b. ASTM D6730 | (vol%)<br> a - b |
|----------|-------------------------|---------------|------------------|
| Benzene  | 0.5                     | 0.552         | 0.052            |
| Toluene  | 16.2                    | 16.529        | 0.329            |
| Xylene   | 10.7                    | 10.755        | 0.055            |
| MeOH     | 1.5                     | 1.414         | 0.086            |
| EtOH     | 5.3                     | 5.194         | 0.106            |
| MTBE     | 5.2                     | 5.268         | 0.068            |
| ETBE     | 5.2                     | 5.187         | 0.013            |



- 1: Ethanol
- 2: C5 (n-pentane)
- 3: tert-Butanol
- 4: 2-Methylbutene-2
- 5: 2,3-Dimethylbutane
- 6: Methyl tert-butyl ether (MTBE)
- 7: C6 (n-hexane)
- 8: 1-Methylcyclopentene
- 9: Benzene
- 10: Cyclohexane

- 11: 3-Ethylpentane
- 12: trans-1,2-Dimethylcyclopentane
- 13: C7 (n-heptane)
- 14: 2,3,3-Trimethylpentane
- 15: Toluene
- 16: C8 (n-octane)
- 17: Ethylbenzene
- 18: p-Xylene
- 19: 2,3-Dimethylheptane
- 20: C9 (n-nonane)

- 21: 5-Methylnonane
- 22: 1-Methyl-2-ethylbenzene
- 23: C10 (n-decane)
- 24: C11 (undecane)
- 25: 1,2,3,5-Tetramethylbenzene
- 26: Naphthalene
- 27: C12 (dodecane)
- 28: 1-Methylnaphthalene
- 29: C13 (tridecane)

Fig. 2 Chromatogram of oxy setup blend

## ■ Data Processing of gasoline sample

After system suitability tests were completed, real gasoline samples were analyzed. Quantitation results are presented as two-dimensional tables by PIONA group and carbon number, and as individual component concentration reports (Fig. 3).

In the chromatogram view of PONAsolution, the identified reference data can be displayed above while performing identification of sample peaks in the lower pane (Fig. 4). By using AART for DHA, alkanes can be identified automatically, and the remaining components are subsequently identified automatically, substantially reducing analysis time\*1.

\*1 Once the library is adjusted, re-adjustment is generally unnecessary until retention time shifts occur due to column degradation.

## ■ Conclusion

Compositional analysis of gasoline in accordance with ASTM D6730 was performed. The combination of SH PONA Tuning Column and SH-1 PONA main column enabled separation of compounds with closely spaced retention times. Moreover, use of PONAsolution software greatly simplifies complex identification tasks.

### <Related Applications>

1. Detailed Hydrocarbon Analysis by GC-VUV Using ASTM D8369, [Application News No. 01-00966-EN](#)

| Group         | Area%   | Wt%     | Vol%    | Mol%    |
|---------------|---------|---------|---------|---------|
| Paraffins     | 15.547  | 15.802  | 17.480  | 15.322  |
| Iso-Paraffins | 31.555  | 32.030  | 34.956  | 29.396  |
| Olefins       | 8.781   | 8.719   | 9.145   | 8.105   |
| Naphthenes    | 3.571   | 3.548   | 3.425   | 3.263   |
| Aromatics     | 25.555  | 24.062  | 20.245  | 18.955  |
| Oxygenates    | 14.232  | 15.069  | 14.184  | 24.958  |
| Unknowns      | 0.759   | 0.769   | 0.565   | 0.000   |
| Total         | 100.000 | 100.000 | 100.000 | 100.000 |

| C#       | Area%   | Mass%   | Vol%    | Mol%    |
|----------|---------|---------|---------|---------|
| C1       | 2.802   | 3.184   | 2.955   | 8.649   |
| C2       | 4.129   | 4.398   | 4.094   | 8.307   |
| C3       | 1.451   | 1.511   | 1.384   | 2.188   |
| C4       | 3.409   | 3.503   | 3.810   | 4.664   |
| C5       | 13.076  | 13.317  | 15.397  | 15.904  |
| C6       | 16.121  | 16.220  | 17.089  | 16.063  |
| C7       | 21.727  | 21.529  | 21.857  | 19.133  |
| C8       | 16.258  | 15.881  | 15.012  | 12.528  |
| C9       | 9.279   | 8.956   | 8.113   | 6.355   |
| C10      | 5.614   | 5.448   | 4.909   | 3.461   |
| C11      | 1.803   | 1.751   | 1.543   | 1.010   |
| C12      | 2.175   | 2.131   | 1.923   | 1.117   |
| C13      | 0.587   | 0.588   | 0.573   | 0.278   |
| C14      | 0.322   | 0.323   | 0.310   | 0.142   |
| C15      | 0.488   | 0.489   | 0.466   | 0.200   |
| C16      | 0.000   | 0.000   | 0.000   | 0.000   |
| Unknowns | 0.759   | 0.769   | 0.565   | 0.000   |
| Total    | 100.000 | 100.000 | 100.000 | 100.000 |

| ID No. | Component                          | R.T.   | RI      | Group | C# | Area       | Wt%   | Vol%  | Mol%  |
|--------|------------------------------------|--------|---------|-------|----|------------|-------|-------|-------|
| 6      | Isobutane                          | 8.209  | 354.168 | I     | 4  | 34862.157  | 0.177 | 0.233 | 0.265 |
| 7      | Methanol                           | 8.594  | 368.593 | X     | 1  | 406399.910 | 3.184 | 2.955 | 8.649 |
| 9      | iso-Butylene                       | 8.976  | 386.523 | O     | 4  | 2482.278   | 0.017 | 0.021 | 0.026 |
| 11     | Butane                             | 9.274  | 400.000 | P     | 4  | 192520.754 | 1.367 | 1.734 | 2.045 |
| 13     | trans-2-Butene                     | 9.658  | 408.900 | O     | 4  | 7547.724   | 0.052 | 0.064 | 0.081 |
| 14     | 2,2-Dimethylpropane                | 9.776  | 411.568 | I     | 5  | 2495.253   | 0.018 | 0.022 | 0.021 |
| 15     | cis-2-Butene                       | 10.249 | 421.930 | O     | 4  | 9797.340   | 0.067 | 0.079 | 0.104 |
| 16     | Ethanol                            | 11.465 | 446.522 | X     | 2  | 598700.323 | 4.398 | 4.094 | 8.307 |
| 18     | 3-Methyl-1-butene                  | 11.837 | 453.519 | O     | 5  | 4645.683   | 0.032 | 0.037 | 0.040 |
| 19     | Isopentane                         | 12.875 | 471.955 | I     | 5  | 874250.487 | 6.163 | 7.306 | 7.434 |
| 24     | 1-Pentene                          | 13.871 | 488.297 | O     | 5  | 13582.111  | 0.093 | 0.107 | 0.115 |
| 25     | Isopropanol                        | 14.061 | 491.285 | X     | 3  | 94432.623  | 0.678 | 0.623 | 0.682 |
| 26     | 2-Methyl-1-butene                  | 14.356 | 495.834 | O     | 5  | 36568.783  | 0.251 | 0.283 | 0.311 |
| 27     | Pentane                            | 14.631 | 500.000 | P     | 5  | 642168.562 | 4.527 | 5.310 | 5.461 |
| 29     | trans-2-Pentene                    | 15.232 | 509.343 | O     | 5  | 39714.292  | 0.272 | 0.308 | 0.338 |
| 31     | cis-2-Pentene                      | 15.793 | 517.731 | O     | 5  | 22063.049  | 0.151 | 0.169 | 0.188 |
| 32     | 2-Methyl-2-Propanol (tert-Butanol) | 15.930 | 519.723 | X     | 4  | 138623.456 | 0.984 | 0.916 | 1.155 |
| 33     | 2-Methyl-2-butene                  | 16.128 | 522.599 | O     | 5  | 78572.850  | 0.538 | 0.597 | 0.668 |
| 37     | 2,2-Dimethylbutane                 | 17.194 | 537.206 | I     | 6  | 38315.037  | 0.269 | 0.304 | 0.272 |
| 39     | Cyclopentane                       | 18.459 | 553.907 | O     | 5  | 7440.675   | 0.050 | 0.047 | 0.063 |
| 40     | n-Propanol                         | 18.697 | 556.878 | X     | 3  | 116010.020 | 0.833 | 0.761 | 1.206 |
| 42     | 3-Methyl-1-pentene                 | 18.914 | 559.551 | O     | 6  | 2769.916   | 0.019 | 0.021 | 0.020 |
| 43     | Cyclopentane                       | 19.263 | 563.700 | N     | 5  | 33353.744  | 0.229 | 0.225 | 0.284 |
| 44     | 2,3-Dimethylbutane                 | 19.446 | 565.977 | I     | 6  | 160298.736 | 1.125 | 1.249 | 1.136 |
| 45     | Methyl tert-butyl ether            | 19.588 | 567.667 | X     | 5  | 140932.910 | 0.994 | 0.985 | 0.961 |
| 47     | 2,3-Dimethyl-1-butene              | 19.745 | 569.522 | O     | 6  | 2925.882   | 0.020 | 0.022 | 0.021 |
| 48     | 2-Methylpentane                    | 19.846 | 570.703 | I     | 6  | 451411.524 | 3.168 | 3.562 | 3.199 |
| 49     | trans-4-Methyl-2-pentene           | 20.013 | 572.650 | O     | 6  | 7872.537   | 0.053 | 0.057 | 0.054 |
| 53     | 3-Methylpentane                    | 20.983 | 583.625 | I     | 6  | 279713.239 | 1.963 | 2.170 | 1.962 |
| 54     | 2-Methyl-1-pentene                 | 21.477 | 589.021 | O     | 6  | 12837.827  | 0.088 | 0.094 | 0.091 |
| 55     | 1-Hexene                           | 21.576 | 590.107 | O     | 6  | 5931.440   | 0.034 | 0.037 | 0.036 |
| 56     | 2-Butanol                          | 21.696 | 593.524 | X     | 4  | 118216.776 | 0.839 | 0.763 | 0.985 |
| 58     | Hexane                             | 22.518 | 600.000 | P     | 6  | 474371.965 | 3.329 | 3.708 | 3.362 |

Fig. 3 Two-dimensional PIONA/carbon number table and component concentration report




Fig. 4 PONAsolution chromatogram view

Nexis, Labsolutions, PONAsolution and AOC are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.



Shimadzu Corporation

[www.shimadzu.com/an/](http://www.shimadzu.com/an/)

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <https://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

01-00863-EN

First Edition: Apr. 2026

➤ Please fill out the survey

## Related Products

Some products may be updated to newer models.



➤ **Nexis GC-2030**  
Gas Chromatograph



➤ **PONAsolution**  
Software for PONA Analysis and  
Detailed Hydrocarbon...

## Related Solutions



➤ **Biomass Fuels**

Hydrocarbon  
➤ Processing Industry  
(Petrochemical, Ch

➤ Petroleum refinery

➤ Clean Energy

➤ Price Inquiry

➤ Product Inquiry

➤ Technical Service /  
Support Inquiry

➤ Other Inquiry