

# Classification and Visualization of Beer Quality Using GC-MS and GC-FID

**ASMS 2019** WP 257

Yusuke Takemori<sup>1</sup>; Yui Higashi<sup>1</sup>; Takero Sakai<sup>1</sup>;

Ryo Takechi<sup>2</sup>; Motoki Sasaki<sup>3</sup>; Narihiro Suzuki<sup>3</sup>

1 Shimadzu Corporation, Kyoto, Japan,

2 Shimadzu Scientific Instruments, Columbia, MD,

3 Ise Kadoya Brewery, Ise, Japan

# Classification and Visualization of Beer Quality Using GC-MS and GC-FID

## Overview


Despite the same brand and identical brewing techniques, it is widely known that beer taste and quality is not consistent from plant to plant. Therefore, brewers normally check beer quality mainly by sensory test method, and in turn, adjust beer taste in order to reduce these differences

as much as possible. In this poster, by metabolomics profiling using GC-MS and GC-FID, we introduce a new approach towards classification and visualization of beer quality to identify how specific components influence taste from plant to plant.

## Introduction

In food science, global metabolite analysis, or 'metabolomics', is increasingly applied to a number of value-added food production areas including food-safety assessment, quality control, food authenticity, origin and processing. GC-MS are widely used for measuring the total amount of metabolites in a sample. Whereas GC-FID might

be less utilized in Metabolomics than GC-MS so far. This is because GC-MS is more powerful in peak annotation, and peak annotation is considered necessary for Metabolomics. But GC-FID has a massive feature; low cost, easy to use, and better repeatability.



**GC-MS(MS)**

- Permits comprehensive measurement of several hundred compounds in a single measurement
- Better peak annotation
- First-choice comprehensive measurement

**GC**

- Simple measurement
- Ideal for efficient routine measurement of specific compounds
- Better repeatability
- Low cost

Figure 1 Instrument Features

# Classification and Visualization of Beer Quality Using GC-MS and GC-FID

In order to introduce a new approach towards classification and visualization of beer quality, we tried doing the following two studies.

### Study1: Classification of different brand beers

Five brand beers were analyzed by GC-FID and GC-MS, and then principal component analysis (PCA) was performed using each GC-FID and GC-MS data. The results of PCA from each GC-FID and GC-MS data were compared. From loading plot of GC-MS data, we can successfully identify how specific components influence this classification.

### Study2: Classification of same brand beers

The different plant and lot beers in the same brand were analyzed by GC-FID and GC-MS, and then principal component analysis (PCA) was performed using each GC-FID and GC-MS data. The results of PCA from each GC-FID and GC-MS data were compared. From loading plot of GC-MS data, we can successfully identify how specific components influence this classification.

## Methods

### Samples

After metabolites were extracted from beers, they were derivatized. Table 1 shows pretreatment procedure for the sample of GC-MS and GC-FID, and Table 2 shows types of analyzed beer samples.

Table 1 Pretreatment procedure




Table 2 Types of analyzed beer samples

Types of beer	Study1	Study2
Lager A	○	
Pale ale A (Plant: a)	○	○
Pale ale A (Plant: b, Lot: a)		○
Pale ale A (Plant: b, Lot: b)		○
Pale ale A (Plant: b, Lot: c)		○
Pale ale B	○	
Pale ale C	○	
IPA A	○	

## Classification and Visualization of Beer Quality Using GC-MS and GC-FID

### Analytical conditions and software

Pretreated beers were analyzed by GC-MS triple quadrupole mass spectrometer and GC-FID. For statistical analysis, SIMCA 15 software (INFOCOM CORPORATION) was used.

Table 3 Analytical condition in GC-MS/MS

<b>GC-MS/MS</b>	: GCMS-TQ8040
<b>Software</b>	: Smart Metabolites Database (475 compounds)
<b>GC conditions</b>	
Column	: DB-5 (30 m × 0.25 mm I.D., 1.00 μm)
Carrier Gas	: He
Injection Temperature	: 280 °C
Control Mode	: Linear velocity (39.0 cm/sec)
Injection Method	: Splitless
Sampling time	: 1 min
Oven Temperature	: 100°C (4 min) - 10 °C/min - 320 °C (11 min)
<b>MS conditions</b>	
Ion Source Temperature	: 200 °C
Interface Temperature	: 280 °C
Tuning Mode	: Standard
Measurement Mode	: MRM
Loop Time	: 0.25 seconds
MRM Transitions	: 950 Transitions

Table 4 Analytical condition in GC

<b>GC</b>	: GC-2030
Column	: SH-Rtx-1 (60 m × 0.32 mm I.D., 1.00 μm)
Carrier Gas	: He
Injection Temperature	: 280 °C
Control Mode	: Linear velocity (25.0 cm/sec)
Injection Method	: Split (1:15)
Oven Temperature	: 40°C - 4 °C/min - 320 °C (15min)
Detector	: FID
Detector Temperature	: 330 °C

# Classification and Visualization of Beer Quality Using GC-MS and GC-FID

## Results

Principal component analyses (PCA) for GC-MS/MS and GC-FID data were performed in the following study 1 and 2.

### Study1: Classification of different brand beers

#### Principal component analyses for GC-MS/MS

Figure2 shows Score plot and loading plot from GC-MS/MS data. From Score plot of figure2, we can confirm that five types of beers were successfully classified. Loading plot of figure2 means which compounds are relatively higher concentration among five beers and how specific compounds influence this classification. Table5 shows higher concentration compounds in each beer from loading plot. In IPA A and Pale ale C, a lot of some sugars were remained.

#### Principal component analyses for GC-FID

Figure3 shows the score plot from GC-FID data. We can get almost the same score plot as one of GC-MS/MS. This result suggest that GC-FID can also used for classification of beers quality. Whereas, it is difficult how specific compounds influence this classification from only FID data, because qualitative analyses performance on GC-FID is less than GC-MS.




Figure 2 Score plot and loading plot from GC-MS/MS data





Figure 3 Score plot from GC-FID data

Table 5 Higher concentration compounds in each beer from loading plot

Types of beer	Higher concentration compounds
Lager A	Phenylpyruvic acid, Glutaric acid Lyxose, Xylose, Arabinose Threo-b-hydroxyaspartic acid 2-Ketoglutaric acid 2-Hydroxyglutaric acid
Pale ale A	Maleic acid, Cadaverine, Maltitol 4-Aminobutyric acid, Dopamine Tryptophan, Oxalic acid
Pale ale C	Galactose, Galacturonic acid Glucose, Mannose, Erythrulose Homogentisic acid, Glucuronic acid Asparagine
IPA A	Sebacic acid, Fructose, Sorbose Tagatose, Psicose

# Classification and Visualization of Beer Quality Using GC-MS and GC-FID

## Study2: Classification of same brand beers

### Principal component analyses for GC-MS/MS

Figure4 shows score plot and loading plot from GC-MS/MS data. From Score plot of figure4, we can confirm that the same brand beers brewed in different plants were divided. Plus, different lots of same plant can be successfully classified in the score plot. Loading plot of figure4 means which compounds are relatively higher concentration among the same brand beers and how specific compounds influence this classification. Table6 shows higher concentration compounds in beers brewed in Plant a and b from loading

plot. In the beer brewed in Plant a, a lot of metabolites from sugar were remained. These results suggest that this approach could be a powerful tool for adjusting food quality.

### Principal component analyses for GC-FID

Figure4 shows the score plot from GC-FID data. We can get almost the same score plot as one of GC-MS/MS. This result suggests that GC-FID can also be used for classification of beers quality in same brand.





Figure 4 Score plot and loading plot from GC-MS/MS data




Figure 4 Score plot from GC-FID data

Table 6 Higher concentration compounds in each plant beer from loading plot

Types of plant	Higher concentration compounds
Plant a	3-Phenylactic acid, Trehalose Glyceric acid, Fructose 1-phosphate Nonanoic acid, 2-Hydroxyisobutyric acid Caproic acid, Glucose 6-phosphate Sedoheptulose 7-phosphate Mannose 6-phosphate Glucose 6-phosphate
Plant b	Allose, Lysine, Tyramine, Methionine Glutamic acid, Galactose Phenylpyruvic acid, Tryptamine 2'-Deoxyuridine, Cystamine-d8, Uridine

## Classification and Visualization of Beer Quality Using GC-MS and GC-FID

### Conclusions

- Score plot from GC-MS/MS and GC-FID successfully categorized five brands of beer, different plants and lots in same brand beer.
- From the loading plots in GC-MS/MS, we could identify important compounds which determine differences of beer brands, same brand beers brewed in different plants and product lots.
- This study suggest a new approach towards classification and visualization of beer quality.
- In the future, GC-FID potential for this filed will be tested in more detail.

### Acknowledgments

We have had the great support of CEO Narihiro Suzuki and Mr. Motoki Sasaki of Ise Kadoya Brewery.

Disclaimer: The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures.