

High Mass MALDI-TOF MS Strategy for Antibody Aggregate Characterization



Brian Feild¹, William Hedgepeth¹, Suvajyoti Guha^{2,3}, Ryan J. Wenzel⁴, Michael R. Zacchariah

¹Shimadzu Scientific Instruments (Columbia, MD), ²Mechanical Eng. Dept. University of Maryland (College Park, MD) ³National Institute of Standards and Technology (Gaithersburg, MD) ⁴CovalX (Zurich, Switzerland)


Introduction

Antibody aggregation remains a complex question to address in the development and manufacture of therapeutic antibodies and biosimilars. This primary degradation product can lead to several undesirable consequences such as immunological response and decreased efficacy. Current technologies for analyzing aggregation products include size exclusion chromatography (SEC), light scattering and analytical ultracentrifugation; however each of these techniques have limitations. As a result, regulatory agencies such as the FDA are leaning towards requesting complimentary data to improve characterization of therapeutic proteins^{1, 2}. Recent developments in High Mass Matrix Assisted Laser Desorption Time of Flight Mass Spectrometry (HM MALDI-TOF MS) can fill these gaps by providing higher resolution data with improved mass accuracy to a mass range of 1.5 – 2.0 MDa. This presentation focuses on illustrating the capabilities of an integrated HM MALDI-TOF MS, the MegaTOF which integrates a standard AXIMA MALDI-TOF MS (Shimadzu, Kyoto Japan) with the high mass detection capabilities of the CovalX detector (CovalX, Zurich, Switzerland)


Background





Therapeutic Protein Aggregation



Crosslinking Chemistry



MegaTOF Aggregation Workflow



Summary

The AXIMA MegaTOF provides new opportunities to biopharmaceutical researchers to better understand and characterize protein aggregation. This HM MALDI-TOF MS system extends the mass accuracy, resolution, speed and user friendly interface provided my MALDI instrumentation to a mass range previously unavailable with this technique. Combined with CovalX cross-linking reagents and the Complex Tracker Software, the MegaTOF provides a complete solution for identifying aggregation products up to 1.5 MDa – 2.0 MDa.

References

- 1 Kozlowski S. M.D. *Potential Need for Measurement Standards to Facilitate R&D of Biological Drugs*. Congressional Testimony 09/24/2009 -<http://www.fda.gov/NewsEvents/Testimony/ucm183596.htm>
- 2 Arakawa T., et al., *Aggregation Analysis of Therapeutic Proteins Part 1*. Bioprocess International, 2006, 4(10): 42-43.
- 3 Seyfried B.K., et al., *MALDI linear TOF mass spectrometry of PEGylated (glyco)proteins*. J Mass Spectrom. 2010 Jun;45(6):612-7.