

High Performance Mass Detection

AXIMA Assurance



AXIMA Assurance™ - Flexibility and Quality


The AXIMA Assurance is designed with the general analytical and life science laboratory in mind. Incorporating a variable repetition rate 50Hz N₂ laser, the system provides high quality, rapid MALDI mass spectra and an array of software tools for data processing and reporting.




Positive and negative ion modes are included as standard, allowing greater flexibility and extending the compound categories that may be analyzed. The system also incorporates a patented beam blower to optionally remove unwanted low mass ions and prevent detector saturation.

High sensitivity is achieved using near-axis laser irradiation and advanced ion optics for enhanced ion transmission. Pulsed extraction of ions from the MALDI source improves resolution and enhanced calibration algorithms with easy-to-use software provide more accurate data.

Unparalleled flexibility is achieved by a variety of sample target formats including standard microtitre plate format 96 or 384 well targets. FlexiMass™ microscope slides (plain or 48 well targets) and a wide variety of adaptors for unconventional sample layouts are also available. The standard sample target formats are fully compatible with common laboratory robots.



ACTH 18-39 demonstrating linear resolution of >6000 FWHM



Typical oligonucleotide spectrum

AXIMA Assurance™ - Software Solutions


Intuitive software for achieving the maximum result with minimum user input, ideal for novice and expert users alike.

The AXIMA Assurance is controlled by the Launchpad™ suite of software, common to all AXIMA mass spectrometers, permitting manual or fully automated operation, facilitating the seamless analysis of as few or as many samples as required. Intelligent optimization of acquisition conditions may be employed permitting auto-tuning for specific samples.


Ideally suited for high-throughput QA/QC application areas, for example:

- Oligonucleotides/primers
- Synthetic peptides/proteins
- Polymer analysis
- Small molecules

Application-centric data processing software packages are available to provide solutions to many commonly asked questions.



Synthetic polymer spectrum



Typical protein mixture

Synthetic polymers

Polymers and copolymers can be characterized using our unique polymer software, Polymer Analysis™, providing useful structural information and statistics in a text report format.

Oligonucleotides

Oligo Analysis™ offers fully automated QC analysis of large numbers of oligonucleotides or peptides, complete with a report indicating the presence or absence of the target compound, an estimate of the purity and occurrence of known contaminants, adducts or truncated/expanded analogues.

Biomarker discovery

This exciting area encompassing clinical sample screening and profiling is comprehensively addressed using automated acquisition methods and refined data processing.

Data can be easily exported to third party software packages to allow comparative experiments using a number of standard data formats including ASCII, mzXML and mzData.

System support

All AXIMA systems can be fully supported throughout their lifetime using sophisticated web-based service diagnostics and realtime remote monitoring. Highly trained specialist local service support engineers are available to install and maintain AXIMA mass spectrometers. A wide range of service contracts is available, catering for all budgets and requirements, including IQ/OQ environments and high-throughput QA laboratories.

Full training courses are offered by MALDI experts at our regional corporate training centers or at the customer site and may be tailored for specific requirements and applications.



Kratos Analytical Ltd.

Wharfside, Trafford Wharf Road, Manchester M17 1GP, UK

Phone: +44 161 888 4400 Fax: +44 161 888 4402

URL <http://www.shimadzu.com/an/>



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice.

Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.