

GC/MS/MS Forensic Toxicology Database


# Quick-DB Forensic



# Provides a Total Solution for All Processes of Toxicological Analysis

Quick-DB™ Forensic is a GC/MS/MS database for forensic toxicological analysis that includes information for the sample preparation, analysis and data processing of 68 toxicological drug components common in toxicology cases. In combination with the GCMS-TQ™ series Gas Chromatograph-Mass Spectrometer, it provides a total solution for all analytical processes.


Due to the simple QuEChERS method used for the sample preparation, even first-time users can pretreat samples easily. In addition to analytical conditions, such as for optimized MRM transitions, the database also includes retention indices, calibration curve information, and other data processing conditions. The database also includes an automatic method creation function (Smart MRM™) for creating analytical methods easily. Due to the preexisting calibration curve information, semi-quantitative values can be calculated easily without using standards. Furthermore, it features a performance management function that automatically determines the system status. By using this function, users can quickly identify areas with low performance to ensure optimized performance throughout the analysis.




## Simple-to-Operate Dedicated Software

The dedicated launcher displayed on-screen helps navigate the necessary operations. Therefore, even first-time users can analyze samples easily by simply following the instructions on the launcher.

During autotuning or when measuring performance check (negative control) samples, the software automatically uses a performance check function to check the status of the MS, consumable use times, and columns in order to evaluate the overall status of the system.




Quick-DB Forensic Launcher Screen




## Navigates Appropriate Pretreatment Steps

For toxicological analysis of whole blood samples, the extraction and purification processes have been optimized using the QuEChERS method.

Preparing samples is simplified by using the Sample Preparation navigator to display photographs with instructions for each step of the QuEChERS cleanup process.




Sample Preparation Navigator



Overview of QuEChERS Protocol






## Achieve Analysis without Standards

To improve the accuracy of semi-quantitation without a standard, the internal standard that is selected and the recovery rate of sample preparation need to be considered. For internal standards, the Quick-DB Forensic database uses six types of deuterium-labeled toxicological substances that exhibit the same behavior as registered compounds. By using the calibration curve information from a whole blood sample spiked with a standard and pretreated, the quantitative error due to matrix interference, recovery rates of sample preparation, and other factors, which can have a significant effect on quantitative accuracy, can be reduced and accurate semi-quantitation values can be obtained without using a standard.


### Using a Deuterium-Labeled Internal Standard to Correct for Sensitivity Variations

Deuterium-labeled internal standards can also be used to correct for GC-MS/MS sensitivity variations, even if the variations cause the response to target compounds in calibration curves to be different than in a customer's instrument.



### Using a Calibration Curve for a Pretreated Sample to Correct for Recovery Rates


By using the calibration curve information from a pretreated sample, accurate quantitative values can be obtained. The calibration curve information is used to correct differences in recovery rates for each target compound or correct for effects caused by adsorption, which occurs when creating calibration curves using standards.



## Smart Scheduling Injects Blanks Automatically

If a high concentration of a target compound is detected in a test sample, carryover within the system may cause doubt about a false positive in the next test sample. If that occurs, the test sample is normally measured again to confirm the results. This may be an impossible task if only a small quantity of sample is available and pretreatment must be done.

Therefore, the Quick-DB Forensic database includes a Smart Scheduling function. This function automatically measures a solvent blank if the concentration of a target compound measured in a test sample is higher than a specified value. Smart Scheduling can prevent the risk of carryover causing a false positive.




## Automates Time-Consuming Derivatization Processes\*

Due to the limited amount of mass spectral information available for phenethylamine drugs, such as methamphetamine, MDA, and MDMA, as well as their tendency for adsorption, they are normally derivatized for measurement. The manual derivatization process requires significant time and can be difficult and lengthy. In contrast, the Quick-DB Forensic database enables on-column TFA derivatization using MBTFA (N-Methylbis-trifluoroacetamide). The test samples can be automatically TFA-derivatized and measured by simply placing the test samples and derivatizing agent in the auto-injector.

### Manual TFA Derivatization Method


Dry the sample under the stream of nitrogen gas.  
10 to 20 minutes



Anhydrous trifluoroacetic acid  
15 minutes at 55 °C


Dry the sample and dissolve it using solvent.  
10 minutes

**Total 40 minutes**



### On-Column TFA Derivatization Method

**Fully Automates Derivatization Process**





**Step 1** An MBTFA derivatizing agent and test sample are successively aspirated into a syringe.

**Step 2** The mixture is injected into the sample injection port. The sample is vaporized in the injection port and delivered to the capillary column.


**Step 3** After waiting a few seconds for the target components to enter the capillary column, the MBTFA derivatizing agent is injected.

**Step 4** The target components are derivatized inside the capillary column and then analyzed by GC-MS.





Mass chromatogram of phenethylamine drug (concentration: 10 ng/mL) added to whole blood sample

\*AOC-20i and AOC-20s are required for this automation.



## More Efficient Multianalyte Data Processing

Quantitation results for a series of data sets can be displayed at the same time for data processing using LabSolution Insight. In addition, chromatograms for each test sample data can be displayed side-by-side for each target compound, making it easy to confirm semi-quantitation results.



## Performance Information Obtained Using this System Can Be Confirmed

System performance results, such as the limit of detection (LOD), limit of quantitation (LOQ) and recovery of sample preparation, using this system are registered. The quantitation limit, pretreatment recovery rate, and other values for detected compounds can be confirmed as reference information. If customer wants to increase the number of registered compounds, the information can also be registered.

ID <sup>1</sup>	Compound Name	Classification	Cal. Type	Cal. Range (ng/mL)	LOD <sup>2</sup> (ng/mL)	LOQ <sup>2</sup> (ng/mL)	Accuracy (Low)	Accuracy (High)		
							Spiked Conc. (ng/mL)	%Accuracy	Spiked Conc. (ng/mL)	%Accuracy
8	Barbital	Barbiturates	Linear	1-500	0.5	1.7	25	105.2	250	90.8
9	Allobarbitol	Barbiturates	Linear	5-500	0.7	2.1	25	76.3	250	90.6
10	Amobarbital	Barbiturates	Linear	1-500	0.6	1.6	25	96.8	250	90.8
11	Penobarbital	Barbiturates	Linear	1-500	0.4	1.3	25	97.5	250	92.0
12	Secobarbital	Barbiturates	Linear	1-500	0.4	1.1	25	98.8	250	92.1
13	Phenobarbital	Barbiturates	Linear	1-500	0.5	1.6	25	102.9	250	86.9
14	Primidone	Barbiturates	Linear	1-500	0.8	2.3	25	96.6	250	101.8
15	Medazepam	Benzodiazepines	Linear	1-500	0.6	1.8	25	102.1	250	97.8
16	Oxazepam	Benzodiazepines	Linear	10-500	2.3	7.2	25	79.0	250	107.7
17	Fludazepam	Benzodiazepines	Linear	1-500	0.2	0.7	25	105.1	250	98.7
18	Diazepam	Benzodiazepines	Linear	1-500	0.3	0.8	25	115.5	250	99.6
19	Clobazepam	Benzodiazepines	Linear	1-500	0.3	0.8	25	91.0	250	100.4
20	Clobazam	Benzodiazepines	Linear	1-500	0.2	0.6	25	113.5	250	91.7
21	Flunitrazepam	Benzodiazepines	Linear	1-500	0.7	2.1	25	83.0	250	91.2
22	Bromazepam	Benzodiazepines	Linear	5-500	1.4	4.4	50	98.7	250	104.6
23	Prizepam	Benzodiazepines	Linear	1-500	0.5	1.6	25	81.9	300	88.8
24	Nimetazepam	Benzodiazepines	Linear	1-500	0.9	2.6	25	79.6	250	80.9
25	Flurazepam	Benzodiazepines	Linear	1-500	0.3	0.9	25	92.1	250	97.5
26	Nitrazepam	Benzodiazepines	Linear	5-300	0.5	1.5	25	148.8	250	87.2
27	Bromazepam	Benzodiazepines	Linear	5-300	0.7	1.9	25	146.1	250	85.3
28	Esazolam	Benzodiazepines	Linear	1-300	0.4	1.3	25	146.9	250	84.8
29	Alprazolam	Benzodiazepines	Linear	1-500	0.4	1.2	25	105.6	250	99.9
30	Tofazepam	Benzodiazepines	Linear	5-500	1.5	4.5	50	79.6	300	92.7

## Performance Report for Quickly Confirming System Status

System performance details displayed on the launcher in real time can also be output in specialized evaluation reports. Two performance management methods are used: one to evaluate instrument performance via MS tuning and the other to evaluate system performance by measuring the performance check sample. These two methods ensure any problems can be confirmed in an easy-to-understand manner.

Instrument Performance Report																							
Tuning File Name: C:\GCMSolution\bin\MS Tuning Files\MS Tuning 0214																							
Tuning Date/Time: 2016/02/14																							
<b>Comprehensive Check</b>		<b>Pass</b>																					
<b>Results in Detail</b>																							
<b>1. FWHM</b>	Pass	<table border="1"> <thead> <tr> <th>m/z</th> <th>FWHM (%)</th> <th>Pass Criteria</th> <th>Check Result</th> </tr> </thead> <tbody> <tr> <td>69</td> <td>0.38</td> <td>0.50</td> <td>- 0.12</td> <td>Pass</td> </tr> <tr> <td>87</td> <td>0.38</td> <td>0.50</td> <td>- 0.12</td> <td>Pass</td> </tr> <tr> <td>92</td> <td>0.38</td> <td>0.50</td> <td>- 0.12</td> <td>Pass</td> </tr> </tbody> </table>	m/z	FWHM (%)	Pass Criteria	Check Result	69	0.38	0.50	- 0.12	Pass	87	0.38	0.50	- 0.12	Pass	92	0.38	0.50	- 0.12	Pass	<p>Evaluate the resolution of the mass spectrum by checking the FWHM of each peak.</p>	
m/z	FWHM (%)	Pass Criteria	Check Result																				
69	0.38	0.50	- 0.12	Pass																			
87	0.38	0.50	- 0.12	Pass																			
92	0.38	0.50	- 0.12	Pass																			
<b>2. Detector voltage</b>	Pass	<table border="1"> <thead> <tr> <th>Detector Voltage (kV)</th> <th>Check Result</th> </tr> </thead> <tbody> <tr> <td>0.00</td> <td>1.80 kV Pass</td> </tr> </tbody> </table>	Detector Voltage (kV)	Check Result	0.00	1.80 kV Pass	<p>Evaluate the state of the ion source and operation of the detector by checking the detector voltage.</p>																
Detector Voltage (kV)	Check Result																						
0.00	1.80 kV Pass																						
<b>3. High Mass Relative Intensity</b>	Pass	<table border="1"> <thead> <tr> <th>m/z</th> <th>Intensity (%)</th> <th>Pass Criteria</th> <th>Check Result</th> </tr> </thead> <tbody> <tr> <td>802</td> <td>4.05%</td> <td>2.00%</td> <td>-</td> <td>Pass</td> </tr> </tbody> </table>	m/z	Intensity (%)	Pass Criteria	Check Result	802	4.05%	2.00%	-	Pass	<p>Evaluate whether it is in the right mass range by checking the relative intensity of m/z 802.</p>											
m/z	Intensity (%)	Pass Criteria	Check Result																				
802	4.05%	2.00%	-	Pass																			
<b>4. Vacuum Leakage</b>	Pass	<table border="1"> <thead> <tr> <th>SEC Pressure (Pa)</th> <th>Pass Criteria</th> <th>Check Result</th> </tr> </thead> <tbody> <tr> <td>227.80</td> <td>0.00</td> <td>Pass</td> </tr> </tbody> </table>	SEC Pressure (Pa)	Pass Criteria	Check Result	227.80	0.00	Pass	<p>Evaluate vacuum leakage by checking the SEC pressure of m/z 802.</p>														
SEC Pressure (Pa)	Pass Criteria	Check Result																					
227.80	0.00	Pass																					
<b>5. Ion Source Status</b>	Pass	<table border="1"> <thead> <tr> <th>Sample Voltage (V)</th> <th>Check Result</th> </tr> </thead> <tbody> <tr> <td>0.1</td> <td>0.0 V Pass</td> </tr> </tbody> </table>	Sample Voltage (V)	Check Result	0.1	0.0 V Pass	<p>Evaluate the state of the ion source by checking the voltage of LENS 1.</p>																
Sample Voltage (V)	Check Result																						
0.1	0.0 V Pass																						

Instrument Performance Report

System Management Report																																																																																																																																																																				
Sample Name: C:\GCMSolution\bin\MS Tuning Files\MS Tuning 0214																																																																																																																																																																				
Sample No: 010001100																																																																																																																																																																				
Reference: 010001100																																																																																																																																																																				
<b>Comprehensive Check</b>		<b>Pass</b>																																																																																																																																																																		
<b>1. Ion Ratio</b>																																																																																																																																																																				
Pass	<table border="1"> <thead> <tr> <th>Compound Name</th> <th>Ratio (%)</th> <th>Reference Ratio (%)</th> <th>Pass Criteria</th> <th>Check Result</th> </tr> </thead> <tbody> <tr> <td>Benzo(a)P</td> <td>30.27%</td> <td>40.00%</td> <td>37.51%</td> <td>- 6.73%</td> <td>Pass</td> </tr> <tr> <td>Benzo(b)F</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(k)F</td> <td>17.28%</td> <td>15.00%</td> <td>16.25%</td> <td>- 1.27%</td> <td>Pass</td> </tr> <tr> <td>Benzo(a)A</td> <td>18.20%</td> <td>15.00%</td> <td>16.67%</td> <td>- 1.47%</td> <td>Pass</td> </tr> <tr> <td>Benzo(b)A</td> <td>18.20%</td> <td>15.00%</td> <td>16.67%</td> <td>- 1.53%</td> <td>Pass</td> </tr> <tr> <td>Benzo(e)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(f)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(g)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(h)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(i)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(j)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(l)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(m)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(n)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(o)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(p)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(q)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(r)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(s)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(t)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(u)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(v)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(w)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(x)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(y)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> <tr> <td>Benzo(z)P</td> <td>10.01%</td> <td>10.00%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> </tbody> </table>	Compound Name	Ratio (%)	Reference Ratio (%)	Pass Criteria	Check Result	Benzo(a)P	30.27%	40.00%	37.51%	- 6.73%	Pass	Benzo(b)F	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(k)F	17.28%	15.00%	16.25%	- 1.27%	Pass	Benzo(a)A	18.20%	15.00%	16.67%	- 1.47%	Pass	Benzo(b)A	18.20%	15.00%	16.67%	- 1.53%	Pass	Benzo(e)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(f)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(g)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(h)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(i)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(j)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(l)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(m)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(n)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(o)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(p)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(q)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(r)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(s)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(t)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(u)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(v)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(w)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(x)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(y)P	10.01%	10.00%	10.00%	- 0.01%	Pass	Benzo(z)P	10.01%	10.00%	10.00%	- 0.01%	Pass	<p>Evaluate the ratio of the mass spectrum by comparing the relative ratio of each compound.</p>	
Compound Name	Ratio (%)	Reference Ratio (%)	Pass Criteria	Check Result																																																																																																																																																																
Benzo(a)P	30.27%	40.00%	37.51%	- 6.73%	Pass																																																																																																																																																															
Benzo(b)F	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(k)F	17.28%	15.00%	16.25%	- 1.27%	Pass																																																																																																																																																															
Benzo(a)A	18.20%	15.00%	16.67%	- 1.47%	Pass																																																																																																																																																															
Benzo(b)A	18.20%	15.00%	16.67%	- 1.53%	Pass																																																																																																																																																															
Benzo(e)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(f)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(g)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(h)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(i)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(j)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(l)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(m)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(n)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(o)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(p)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(q)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(r)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(s)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(t)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(u)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(v)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(w)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(x)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(y)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
Benzo(z)P	10.01%	10.00%	10.00%	- 0.01%	Pass																																																																																																																																																															
<b>2. Retention Time</b>																																																																																																																																																																				
Pass	<table border="1"> <thead> <tr> <th>Compound Name</th> <th>Retention Time (min)</th> <th>Pass Criteria</th> <th>Check Result</th> </tr> </thead> <tbody> <tr> <td>Benzo(a)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(b)F</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(k)F</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(a)A</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(b)A</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(e)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(f)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(g)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(h)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(i)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(j)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(l)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(m)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(n)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(o)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(p)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(q)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(r)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(s)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(t)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(u)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(v)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(w)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(x)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(y)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> <tr> <td>Benzo(z)P</td> <td>11.1</td> <td>11.000</td> <td>Pass</td> </tr> </tbody> </table>	Compound Name	Retention Time (min)	Pass Criteria	Check Result	Benzo(a)P	11.1	11.000	Pass	Benzo(b)F	11.1	11.000	Pass	Benzo(k)F	11.1	11.000	Pass	Benzo(a)A	11.1	11.000	Pass	Benzo(b)A	11.1	11.000	Pass	Benzo(e)P	11.1	11.000	Pass	Benzo(f)P	11.1	11.000	Pass	Benzo(g)P	11.1	11.000	Pass	Benzo(h)P	11.1	11.000	Pass	Benzo(i)P	11.1	11.000	Pass	Benzo(j)P	11.1	11.000	Pass	Benzo(l)P	11.1	11.000	Pass	Benzo(m)P	11.1	11.000	Pass	Benzo(n)P	11.1	11.000	Pass	Benzo(o)P	11.1	11.000	Pass	Benzo(p)P	11.1	11.000	Pass	Benzo(q)P	11.1	11.000	Pass	Benzo(r)P	11.1	11.000	Pass	Benzo(s)P	11.1	11.000	Pass	Benzo(t)P	11.1	11.000	Pass	Benzo(u)P	11.1	11.000	Pass	Benzo(v)P	11.1	11.000	Pass	Benzo(w)P	11.1	11.000	Pass	Benzo(x)P	11.1	11.000	Pass	Benzo(y)P	11.1	11.000	Pass	Benzo(z)P	11.1	11.000	Pass	<p>Evaluate the accuracy of retention time comparison by comparing the retention time ratio of each compound.</p>																																																						
Compound Name	Retention Time (min)	Pass Criteria	Check Result																																																																																																																																																																	
Benzo(a)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(b)F	11.1	11.000	Pass																																																																																																																																																																	
Benzo(k)F	11.1	11.000	Pass																																																																																																																																																																	
Benzo(a)A	11.1	11.000	Pass																																																																																																																																																																	
Benzo(b)A	11.1	11.000	Pass																																																																																																																																																																	
Benzo(e)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(f)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(g)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(h)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(i)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(j)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(l)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(m)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(n)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(o)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(p)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(q)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(r)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(s)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(t)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(u)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(v)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(w)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(x)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(y)P	11.1	11.000	Pass																																																																																																																																																																	
Benzo(z)P	11.1	11.000	Pass																																																																																																																																																																	
<b>3. Column</b>																																																																																																																																																																				
Pass	<table border="1"> <thead> <tr> <th>Compound Name</th> <th>Peak Area</th> <th>Pass Criteria</th> <th>Check Result</th> </tr> </thead> <tbody> <tr> <td>Benzo(a)P</td> <td>1.0</td> <td>0.000</td> <td>Pass</td> </tr> </tbody> </table>	Compound Name	Peak Area	Pass Criteria	Check Result	Benzo(a)P	1.0	0.000	Pass	<p>Evaluate the operation of the column by comparing the peak area.</p>																																																																																																																																																										
Compound Name	Peak Area	Pass Criteria	Check Result																																																																																																																																																																	
Benzo(a)P	1.0	0.000	Pass																																																																																																																																																																	
<b>4. Glass Liner/Column</b>																																																																																																																																																																				
Pass	<table border="1"> <thead> <tr> <th>Compound Name</th> <th>Peak Area</th> <th>Pass Criteria</th> <th>Check Result</th> </tr> </thead> <tbody> <tr> <td>Benzo(a)P</td> <td>0.0</td> <td>0.000</td> <td>Pass</td> </tr> </tbody> </table>	Compound Name	Peak Area	Pass Criteria	Check Result	Benzo(a)P	0.0	0.000	Pass	<p>Evaluate the operation of the glass liner and the column by checking the peak area of the mass spectrum (Benzo(a)P) measured with the regular sample (Phenacetin).</p>																																																																																																																																																										
Compound Name	Peak Area	Pass Criteria	Check Result																																																																																																																																																																	
Benzo(a)P	0.0	0.000	Pass																																																																																																																																																																	
<b>5. Recovery of Sample Preparation</b>																																																																																																																																																																				
Pass	<table border="1"> <thead> <tr> <th>Compound Name</th> <th>Recovery Ratio (%)</th> <th>Pass Criteria</th> <th>Check Result</th> </tr> </thead> <tbody> <tr> <td>Benzo(a)P</td> <td>10.01%</td> <td>10.00%</td> <td>- 0.01%</td> <td>Pass</td> </tr> </tbody> </table>	Compound Name	Recovery Ratio (%)	Pass Criteria	Check Result	Benzo(a)P	10.01%	10.00%	- 0.01%	Pass	<p>Evaluate the ratio of each mass spectrum by comparing the relative ratio of the sample and standard.</p>																																																																																																																																																									
Compound Name	Recovery Ratio (%)	Pass Criteria	Check Result																																																																																																																																																																	
Benzo(a)P	10.01%	10.00%	- 0.01%	Pass																																																																																																																																																																

System Performance Report

## Specifications

### Product Contents

Quick-DB Forensic launcher program, database files, method files, batch files, report format files, and 2-step AOC control software

### Applicable Instruments

GC/MS : GCMS-TQ series  
 Auto-injector : AOC™-20i+s  
 Workstation : GCMSsolution™ (Ver. 4.44 or later) + LabSolutions Insight™ (Ver. 2.0 SP1 or later)

### Required Parts and Reagents

Column : SH-Rxi™-5Si1 MS 30 m × 0.25 mm I.D., df = 0.25 μm (P/N 221-75954-30)  
 Pretreatment : Q-sep QuEChERS extraction salts packet (AOAC 2007.01) (Restek Corporation P/N: 26238)  
 : Q-sep QuEChERS cleanup kit (Restek Corporation P/N: 26424)

Internal standard substance : Alprazolam-D5 (SIGMA-ALDRICH, Cat: A-902-1ML)  
 Diazepam-D5 (SIGMA-ALDRICH, Cat:D-902-1ML)  
 Carbamazepine-D10 (SIGMA-ALDRICH, Cat:C-094-1ML)  
 Secobarbital-D5 (SIGMA-ALDRICH, Cat:S-001-1ML)  
 Diazinon-D10 (SIGMA-ALDRICH, Cat:74332-5MG)  
 Phenobarbital-D5 (SIGMA-ALDRICH, Cat:P-018-1ML)  
 1-Methyl-3-phenylpropylamine (SIGMA-ALDRICH, Cat:M70533-25G)

n-Alkane Standard : C8-C40 Alkanes Calibration Standard (SIGMA-ALDRICH, Cat: 40147-U)  
 TFA derivatizing agent : MBTFA

In addition, general pretreatment equipment, such as a centrifuge, is also required.

This product was developed jointly with the Department of Legal Medicine & Bioethics, Nagoya University Graduate School of Medicine.

## List of Compounds (68 target compounds, six internal standard compounds, and two QC compounds)

Compound	Classification	Compound	Classification
Allobarbitol	Barbiturates	Coumaphos	Pesticides
Amobarbital	Barbiturates	Cyanophos	Pesticides
Barbital	Barbiturates	Diazinon	Pesticides
Pentobarbital	Barbiturates	Edifenphos	Pesticides
Phenobarbital	Barbiturates	EPN	Pesticides
Primidone	Barbiturates	Etrimfos	Pesticides
Secobarbital	Barbiturates	Fenitrothion	Pesticides
Alprazolam	Benzodiazepines	Fenthion	Pesticides
Bromazepam	Benzodiazepines	Iprobenfos	Pesticides
Brotizolam	Benzodiazepines	Isofenphos	Pesticides
Clobazam	Benzodiazepines	Isoxathion	Pesticides
Clonazepam	Benzodiazepines	Malathion	Pesticides
Clotiazepam	Benzodiazepines	Methidathion	Pesticides
Diazepam	Benzodiazepines	Phosalone	Pesticides
Estazolam	Benzodiazepines	Terbufos	Pesticides
Etizolam	Benzodiazepines	Amphetamine-TFA	Phenethylamines
Fludiazepam	Benzodiazepines	MDA-TFA	Phenethylamines
Flunitrazepam	Benzodiazepines	MDMA-TFA	Phenethylamines
Flurazepam	Benzodiazepines	Methamphetamine-TFA	Phenethylamines
Medazepam	Benzodiazepines	Sertraline	SSRI
Nimetazepam	Benzodiazepines	Amitriptyline	Tri/Tetracyclics
Nitrazepam	Benzodiazepines	Carbamazepine	Tri/Tetracyclics
Oxazepam	Benzodiazepines	Chlorpromazine	Tri/Tetracyclics
Prazepam	Benzodiazepines	Clomipramine	Tri/Tetracyclics
Tofisopam	Benzodiazepines	Imipramine	Tri/Tetracyclics
Triazolam	Benzodiazepines	Mianserin	Tri/Tetracyclics
Acetylfentanyl	Drugs of Abuse	Nortriptyline	Tri/Tetracyclics
Cocaine	Drugs of Abuse	Promethazine	Tri/Tetracyclics
Diphenidine	Drugs of Abuse	Thioridazine	Tri/Tetracyclics
Fentanyl	Drugs of Abuse	Trimipramine	Tri/Tetracyclics
JWH-018	Drugs of Abuse	1-Methyl-3-phenylpropylamine-TFA	ISTD
JWH-122	Drugs of Abuse	Alprazolam-d5	ISTD
JWH-210	Drugs of Abuse	Carbamazepine-d10	ISTD
Ketamine	Drugs of Abuse	Diazepam-d5	ISTD
MAM2201	Drugs of Abuse	Diazinon-d10	ISTD
THC	Drugs of Abuse	Secobarbital-d5	ISTD
Butamifos	Pesticides	Phenobarbital-d5 (QC)	QC
Chlorpyrifos	Pesticides	Secobarbital-d5 (QC)	QC

## Note

- 1.The accuracy of the information contained in the database and the usefulness of information obtained as a result of the use of this information is not guaranteed.
- 2.Quantitative information obtained using this system are quantitative values determined without using standard samples. To calculate more accurate quantitative values, be sure to perform tests using a standard sample.
- 3.To reliably identify substances registered with this database, perform measurement using the system requirements of the method file included with the product.
- 4.This database is intended for research purposes. It cannot be used for clinical diagnostic applications.

Quick-DB, GCMS-TQ, Smart MRM, AOC, GCMSsolution, LabSolutions Insight are trademarks of Shimadzu Corporation.  
Q-Sep and Rxi are registered trademarks of Restek Corporation.

**For Research Use Only. Not for use in diagnostic procedures.**

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.

Shimadzu Corporation  
www.shimadzu.com/an/