

Enables the Brevis™ GC-2050 to be Controlled by Thermo Scientific™ Dionex™ Chromeleon™ 7

# Shimadzu GC Driver for Chromeleon 7




Brevis™ GC-2050



Nexis™ GC-2030

## ■ Enables Control of the Brevis GC-2050, a System that Offers Uncompromised Analytical Performance in a Small Footprint

The compact Brevis GC-2050, with a system width of 350 mm including the autoinjector (AOC™-30i), enables even more efficient use of laboratory space. Compared to the same configuration of the Nexis GC-2030, the system width has been reduced by approximately 35 %. Despite its compact design, the instrument can accommodate two industry-standard 7-inch capillary columns, and it also supports dual-line simultaneous analysis as well as seamless switching between AOC and HS-20 analysis.




## ■ Facilitates Easy Daily Maintenance with Support for the Easy sTop Function



The GC-2050 is equipped with a function (Easy sTop) that automatically lowers the GC inlet temperature and simplifies liner (consumable) replacement by simply pressing the MAINT. button on the upper right of the GC front panel. It can also be operated from the Chromeleon window, helping to ensure a convenient workflow in the lab.

Note: Operating the Easy sTop function from Chromeleon is possible with both the GC-2030 and GC-2050.




## ■ Eco Idling for Eco-Friendly Analysis – Now with Automatic Eco Mode Support



The system learns the analysis patterns and visualizes the workload, while also suggesting an eco-friendly operation schedule. During analysis, it delivers high performance, and when idle, it automatically switches to an energy-saving ECO idling mode, balancing both analytical tasks and eco-efficiency.

Note: Compatible with GC-2030 and GC-2050.



## Product Lineup

Description	Version
<b>Shimadzu GC Driver for Chromeleon 7</b>	<b>2.50</b>

Note: This product supports Nexis GC-2030, Brevis GC-2050 and HS-20/HS-10 system control and HS-20 standalone control.  
GC-2010 (Plus/Pro) and GC-2014 (c) can be controlled using the Thermo Fisher Scientific Inc. driver provided with Chromeleon 7.

## Controllable Hardware

The Shimadzu GC Driver for Chromeleon 7 supports control of the following units.

**GC Unit** Nexis GC-2030, Brevis GC-2050

**Options** AOC-30i, AOC-20i (Plus) autoinjector, AOC-20s (U) autosampler, HS-20 (NX)/HS-10 headspace sampler, dual injection system

Unit	Device name
<b>Sample Injector</b>	<b>GC-2030</b> : SPL-2030, WBI-2030, OCI-2030 (NX), PTV-2030, SINJ-2030, SPI-U(1.0)
	<b>GC-2050</b> : SPL-U(1.0), PTV-U(1.0), WBI-U(1.0), OCI NX-U(1.0), SINJ-U(1.0), SPI-U(1.0)
<b>Detector</b>	<b>GC-2030</b> : FID-2030, TCD-2030, ECD-2010 Exceed, FPD-2030, FTD-2030, BID-2030, SCD-2030, PTC-2030, AD BOARD
	<b>GC-2050</b> : FID-U(1.0), FPD-U(1.0), ECD-2010 Exceed U, FTD-U(1.0), BID-U(1.0), STCD-U(1.0), AD BOARD
<b>Advanced Flow Technology</b>	<b>GC-2030</b> : Backflush, detector splitting, detector switching, heart-cut system
	<b>GC-2050</b> : Backflush, detector splitting
<b>Additional Temperature Controller</b>	<b>GC-2030</b> : Auxiliary temperature control unit
	<b>GC-2050</b> : Auxiliary temperature control unit
<b>Additional Flow Controller</b>	<b>GC-2030</b> : APC (3 auxiliary channels), APC (1 auxiliary channel)
	<b>GC-2050</b> : APC (3 auxiliary channels), APC (1 auxiliary channel)
<b>Options</b>	<b>GC-2030</b> : Gas selector, Low-temperature control solenoid valve set CRG-2030, External equipment control relay PGC-2030, PRG Box
	<b>GC-2050</b> : Gas selector, Low-temperature control solenoid valve set CRG-2030, PRG-2030, PRG BOX

The following unit can be controlled using the Thermo Fisher Scientific driver included with Chromeleon 7.

By using the Shimadzu GC Driver for Chromeleon in conjunction with standalone control of the HS-20 headspace, the system consisting of both the GC and HS-20 (NX) can be controlled.

**GC Unit** GC-2010 (Plus) and GC-2014

**Options** HS-20 (NX) headspace sampler, AOC-20i (Plus) autoinjector, AOC-20s (U) autosampler

Unit	Device name
<b>Sample Injection Port</b>	<b>GC-2010 (Plus)</b> : SPL-2010 (Plus), OCI / PTV-2010 (Plus) <b>GC-2014</b> : SPL-2014, DINJ-2014
<b>Detector</b>	<b>GC-2010 (Plus)</b> : FID-2010 (Plus), TCD-2010 (Plus), ECD-2010 (Plus), FPD-2010 (Plus), FTD-2010 (Plus) <b>GC-2014</b> : DFID-2014, SFID-2014, TCD-2014, ECD-2014
<b>Additional Flow Controller</b>	<b>GC-2010 (Plus)</b> : APC (3 auxiliary channels)
<b>Low-Temperature Oven Controller</b>	<b>GC-2010 (Plus)</b> : CRG-2010 low-temperature control solenoid valve unit

Note: Up to six Shimadzu GC units can be controlled from a single instrument server.

It is not recommended to connect other vendor instruments to the same instrument server. Please prepare a dedicated server for controlling Shimadzu instruments.

Brevis, AOC and Nexis are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.

Thermo Scientific, Dionex, Chromeleon and Thermo Fisher Scientific are trademarks of Thermo Fisher Scientific Inc. and its affiliated entities.



Shimadzu Corporation

[www.shimadzu.com/an/](http://www.shimadzu.com/an/)

### For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.