

Software for PONA/DHA Analysis

# PONAsolution for LabSolutions

ASTM Standards: D5134, D6729, D6730, D6733.

PONA/DHA analysis determines the individual chemical components in petroleum products and groups them into PONA types based on chemical structure. Understanding these compositions is essential for process control and quality management in refineries. The PONAsolution is an analytical software for LabSolutions that effectively supports PONA/DHA analysis.




PONA: Paraffins, Olefins, Naphthenes, and Aromatics  
DHA: Detailed Hydrocarbon Analysis

## Simple Data Management and Workflow

PONAsolution does not require data file conversion. After analysis, the software can automatically read and analyze data. All settings and results information are stored in a single data file, simplifying data management.

### Workflow for PONA/DHA analysis



- Support for dual line analysis
- Simultaneous connection of 4 GCs and automatic data analysis on a single PC
- Simple Results View
- PDF Report Output/Print Function

Results View

Report

## Rich Identification Support Functions


### AART for DHA |

#### Automatic library correction using retention index

The libraries used for each ASTM standards are fully registered with retention index. With the help of n-alkanes for multiple point correction, high-precision retention time correction can be achieved, allowing for accommodation in the following situations:


- When the column is trimmed or changed
- When the library is installed to a new instrument

AART: Automatic Adjustment of Retention Time



## Two Point Identification | Partial Library Correction

If you want to modify only a specific interval, PONAsolution can automatically modify the retention time of components by specifying a range.




“Two point identification” and “PIONA filter” can be used in combination.

For example, you can modify the retention time of the library with two-point identification, focusing only on the aromatic components (A).

## PIONA Filter

If you only want to check the position of olefin (O), switch the marker display on the chromatogram for each PIONA type.



## Reliable quality and total support

### Application | Detailed Hydrocarbon Analysis (DHA)

Controlling the concentration of components in petroleum products is extremely important for process control and quality control at refineries. Gasoline contains hundreds of components, and DHA/PONA analysis requires high column separation performance. By using the GC-2030 and SH-PONA columns, optimal separation and reliable results can be achieved. For more details, please refer to the application notes.

#### Application

#### System Configuration

GC	GC-2030 (SPL,FID)
Cooling System	CRG-2030 (LN <sub>2</sub> /CO <sub>2</sub> )
Auto Injector	AOC-30i/AOC-20i Series

#### Software requirements

GC Control software	LabSolutions LC/GC after Ver. 5.96 SP3
Software for PONA/DHA analysis	PONAsolution Ver.6
OS	Windows 10, 11
PC	Please check the PC requirements for LabSolutions LC/GC.

### Comprehensive support system

Shimadzu offers high-quality analysis systems combining GC equipment, software and consumables. Leave installation, inspection and after-sales support to us.



AOC, LabSolutions and PONAsolution are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries. Windows is either a registered trademark or a trademark of Microsoft Corporation in the United States and/or other countries.



Shimadzu Corporation

[www.shimadzu.com/an/](http://www.shimadzu.com/an/)

#### For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol “TM” or “®”.

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol “TM” or “®”.

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you “as is” without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.