

Collective analysis

The collective analysis function has been added in version 1.1

The screenshot displays the IMAGEREVEAL software interface, which is used for differential analysis of mass spectrometry data. The interface is divided into several panels:

- Left Sidebar:** Contains a vertical menu of functions. The 'Collectively Analyze' button, located under the 'Differential Analysis' section, is highlighted with a yellow rectangle.
- ROI List:** A table listing regions of interest (ROIs) for analysis.
- Data Matrix Table:** A table showing the data matrix for the selected ROIs.
- MS Image:** A panel displaying a mass spectrometry image (TIC) with a color scale and a 'Copy Information' button.
- Graph:** A panel showing a mass spectrum plot with peaks labeled with their m/z values.
- MS Image List:** A panel displaying a list of mass images (TIC) for the selected ROIs.

The 'Collectively Analyze' button is located in the left sidebar, under the 'Differential Analysis' section, and is highlighted with a yellow rectangle.


No.	Use	File Name	ROI Na...	Attribute
1		Testicle_9A...	All	Group A

No.	Use	Tag	Label	m/z	Formula	Adduct Ion	Matrix	Polarity
-----	-----	-----	-------	-----	---------	------------	--------	----------


No.	Name	Value
1	Normalize	None

m/z	Intensity
721.48186	
767.49182	
795.32084	
796.32363	
797.32374	
798.32545	
837.53900	
885.53782	

Switch on the “collective analysis” and the contents of the menu bar changes



Collective analysis settings screen



(If image registration and ROI settings are necessary, complete these beforehand)

Collectively Analyze

Pre-processing
None

Data Matrix
Non-target
Auto
0.2Da

Test

PCA

PLS

Save project file

Normalize
None TIC XIC

Import Export

No.	Use	m/z	Tolerance
-----	-----	-----	-----------

It is possible to change settings for pre-processing, data matrix settings, test, PCA, PLS, and project file saving.

It is possible to select testing, PCA, PLS, and project file saving.

Reference Value Setting

Minimum Threshold Value (%) 0.00

Specified Method
☐ Range
☒ Center + Tolerance

Back Next Confirm Cancel

1. Pre-processing (normalization)

Collectively Analyze

Pre-processing
None

Data Matrix
Non-target
Auto
0.2Da

Test

PCA

PLS

Save project file

Normalize
None TIC XIC

Import Export

No.	Use	m/z	Tolerance
Spectrum normalization method can be set.			

☐ Reference Value Setting

Minimum Threshold Value (%) 0.00

Specified Method
☐ Range
☒ Center + Tolerance

Back Next Confirm Cancel

2. Data matrix settings

Collectively Analyze

Pre-processing
None

Data Matrix
Non-target
Auto
0.2Da

Test

PCA

PLS

Save project file

Analysis Method
Target Non-target

Threshold Value
0.000 %

m/z Range
Auto
Manual 10.00000 - 1000.00000 Da Set the file range

Bin Size
0.2000 Da

Labeling
Matrix Clusters

Specified Peak Exclusion
Exclude Specified m/z

For the data matrix settings, you can choose between target and non-target settings.

If using peak picking results, first carry out peak picking from the main window. If the ROI spectrum is necessary at this point, data matrix calculations should be carried out beforehand in the normal procedure.

Back Next Confirm Cancel

2. Data matrix settings

If peak picking is necessary:

The screenshot displays the IMAGEREVEAL software interface with several panels. The 'Data Matrix Table' panel is empty. The 'MS Image' panel shows a color-coded mass spectrum image. The 'Graph' panel displays a mass spectrum plot with peaks labeled at m/z 721.48186, 767.49182, 796.43363, 797.42374, 837.53900, and 885.53782. The 'Analysis Parameters' panel shows a table with one row: '1' under 'No.', 'Normalize' under 'Name', and 'None' under 'Value'. The 'MS Image List' panel shows a table with one row: 'Testicle_9AA_PL_5x_1_AREA01.i.mdx' under 'File Name' and 'TIC' under 'Type'. The 'Peak Picking' button in the 'Graph' panel is highlighted with a yellow box and a circled '2'. The 'Testicle_9AA...' checkbox in the 'MS Image List' panel is highlighted with a yellow box and a circled '1'.

ROI List

No.	Use	File Name	ROI Name	Attribute
1		Testicle_9A...	All	Group A

Data Matrix Table

No.	Use	Tag	Label	m/z	Formula	Adduct Ion	Matrix	Polarity
-----	-----	-----	-------	-----	---------	------------	--------	----------

MS Image

Compound Name/Comment:
TIC

File Name:
Testicle_9AA_PL_5x_1_AREA01.i.mdx

Type:
TIC

Copy Information

Graph

Spectrum Box Plot

Intensity

m/z

Peak Picking

MS Image List

D...	File Name	Sp
<input checked="" type="checkbox"/>	Testicle_9AA...	Wh

Testicle_9AA_PL_...

TIC

Analysis Parameters

No.	Name	Value
1	Normalize	None

2. Data matrix settings

Peak picking

Peak Picking

Parameter Settings

m/z Range: 699.98492 - 900.01906 Da Set Default

Smoothing: Savitzky-Golay

Number of Data Points: 9

Number of peaks to detect: 1000

☐ Threshold Value: 0.000 %

☐ Detect Monoisotopic

Minimum Peak Number for Isotope Cluster: 1

Matching Tolerance (ppm): 1

☐ Specified Peak Exclusion: Exclude Specified m/z

Tolerance: 0.2000 Da

Execute

Spectrum Graph

☒ Show Peaks Add to Peak List

Intensity

m/z

Peak List 702 Peak

No.	m/z	Intensity
1	795.52230	1443185.71229
2	796.52490	681252.41170
3	797.52445	261516.97704
4	809.50950	180146.24418
5	767.49232	119581.51946
6	885.53868	114637.12529
7	810.51225	87232.64859
8	796.01916	80876.03661
9	798.52494	75169.20472
10	795.78487	63858.99839
11	886.54092	61671.66389
12	768.49495	54027.83907
13	797.04571	48932.13258
14	796.78374	47286.68287
15	837.53880	47282.07643
16	857.50787	45927.95820
17	793.50711	44473.40080
18	721.47935	39547.56227
19	823.54383	39259.81788
20	821.53421	37386.97173
21	811.51235	34631.03988
22	796.29718	34286.77276
23	883.51835	28751.02937
24	794.51308	27390.47551
25	887.54315	25609.46170
26	748.51008	25004.12189

Save As Compound Template Use As Target List Close

1

2

2. Data matrix settings

The peak picking results are automatically input and become target peaks

Collectively Analyze

Analysis Method: **Target** Non-target ☐ Threshold Value: 0.000 %

Pre-processing: ☒ None

Data Matrix ☒ Target Peak List 0.2Da

Test: ☒

PCA: ☒

PLS: ☒

Save project file: ☒ !

Compound List

Used Compound Template: Peak List

Excluded Compound Template:

	<input checked="" type="checkbox"/>	m/z	Compound Name	Formula	Matrix	Polarity	Adduct Ion
1	<input checked="" type="checkbox"/>	795.52230	795.52230		Any	Bipolar	
2	<input checked="" type="checkbox"/>	796.52490	796.52490		Any	Bipolar	
3	<input checked="" type="checkbox"/>	797.52445	797.52445		Any	Bipolar	
4	<input checked="" type="checkbox"/>	809.50950	809.50950		Any	Bipolar	
5	<input checked="" type="checkbox"/>	767.49232	767.49232		Any	Bipolar	
6	<input checked="" type="checkbox"/>	885.53868	885.53868		Any	Bipolar	
7	<input checked="" type="checkbox"/>	810.51225	810.51225		Any	Bipolar	
8	<input checked="" type="checkbox"/>	796.01916	796.01916		Any	Bipolar	
9	<input checked="" type="checkbox"/>	798.52494	798.52494		Any	Bipolar	
10	<input checked="" type="checkbox"/>	795.78487	795.78487		Any	Bipolar	
11	<input checked="" type="checkbox"/>	886.54092	886.54092		Any	Bipolar	
12	<input checked="" type="checkbox"/>	768.49495	768.49495		Any	Bipolar	
13	<input checked="" type="checkbox"/>	797.04571	797.04571		Any	Bipolar	
14	<input checked="" type="checkbox"/>	796.78374	796.78374		Any	Bipolar	
15	<input checked="" type="checkbox"/>	837.53880	837.53880		Any	Bipolar	
16	<input checked="" type="checkbox"/>	857.50787	857.50787		Any	Bipolar	
	<input checked="" type="checkbox"/>	793.50711	793.50711		Any	Bipolar	

Tolerance: 0.0200 Da

Back Next Confirm Cancel

3. Test

Collectively Analyze

☒ Pre-processing
None

☒ Data Matrix
Non-target
Auto
0.2D

☒ Test

☒ PCA

☒ PLS

☒ Save project file !

There are no setting items.

Back Next Confirm Cancel

No configuration items are displayed, but ROI settings with two or more attributes are required

PCA

Collectively Analyze

Pre-processing
None

Data Matrix
Target
Peak List
0.2Da

Test

PCA

PLS

Save project file


Number of Principal Components
☒ Auto
☐ Manual 5


Pre-processing
Pareto Scale


Back Next Confirm Cancel


PLS


Collectively Analyze


☒  Pre-processing
None

☒  Data Matrix
Target
Peak List
0.2Da

☒  Test

☒  PCA

☒  PLS

☒  Save project file

Number of Latent Variables
☐ Auto
☒ Manual 5

Pre-processing
Pareto Scale

ROI List

Import

Export

No.	File Name	ROI Name	Attribute	Y value
1	Testicle_9AA_PL_SL_5x_1...	ROI001	Group A	1.00000
2	Testicle_9AA_PL_SL_5x_1...	ROI002	Group B	2.00000
3	Testicle_9AA_PL_SL_5x_1...	ROI003	Group C	3.00000

Y values must be entered to use PLS.

Back

Next

Confirm

Cancel

Saving a project file

Collectively Analyze

Save As

Input is required.

To save the file, select the destination folder and enter the name of the project file to be saved.

Uncheck the box if you do not want to save the file.

Pre-processing
None

Data Matrix
Target
Peak List
0.2Da

Test


PCA

PLS


Save project file

Back Next Confirm Cancel

Once all settings are complete, click “Confirm”















A confirmation window appears



The execution details are displayed.

Collectively Analyze

	Pre-processing None	 Executing	Processing Time : The pre-processing time is included in data matrix creation.
	Data Matrix Target Peak List 0.02Da	 Executing	Processing Time : --:--
	Test	 Not executed	Processing Time : --:--
	PCA	 Not executed	Processing Time : --:--
	PLS	 Not executed	Processing Time : --:--
	Save project file	 Not executed	Processing Time : --:--
			Total processing time :--:--

Create Data Matrix

Calculating...25792

Cancel

Close

Once the process is complete, analysis results appear on the main window

