Shimadzu Packed Column for HPLC

Shim-pack™ SCR-101N

Instruction Manual

■ PRINCIPLE FOR SEPARATION

AND PROCEDURE OF APPLICATION

● Principle for Separation:

Shim-pack SCR-101N is a packed column for sugar analysis which separation mode is a size exclusion chromatography. In this case, water is used as the mobile phase and the degree of retention depends on the size of solute molecule relative to the size of pore of packings.

Small molecules will permeate the pores and eluate slowly. On the other hand, relatively large molecules are excluding from pores and eluate faster.

In case of sugars, they will be eluted in the order of trisaccharides (such as raffinose, etc.), disaccharides (such as cane sugar, etc.), and monosaccharides (such as grape sugar, fruit sugar, etc.).

● Column Packing Material:

This Column is packed with a cation exchange resin which is sulfonated polystyrene-divinylbenzene copolymer.

● Analyzable Samples:

The COLUMN is available exclusive to the analyses of sugars. However, it is also available for analyses of any neutral materials which are water-soluble and have several hundreds of molecular weights. For example, it is also available for separation of sugar-oilgomer, alcohols, etc. As basic materials may be absorbed into the packings preventing the elution from the COLUMN, take care not to inject them into the COLUMN too much.

It is recommended to use the COLUMN for analyses of sugars contained in fruits, juice, Japanese "SAKE", wine, syrup, and beverages, etc.

● Mobile Phase:

For the COLUMN, use distilled water or deionized water as the mobile phase.

CAUTION: Take care of such that changing the pH of water or addition of any organic solvents to the water will swell the particles of the packings more or less, thereby possibly causing any change of the bed of COLUMN, resulting in deteriorating the performance of COLUMN.

In case the COLUMN is to be mounted on a System used another mobile phase before, it is required to clean the flowing channel with water completely and replaces the mobile phase with water completely.

As letting any outside air get into the COLUMN forming bubbles may deteriorate the performance of COLUMN, it is necessary to completely remove the air from the flowing channel of System and then mount the COLUMN on the System.

● Pressure and Flow Rate:

The most adequate flow rate of the mobile phase in case of the COLUMN is approximately 0.5 to 1 mL/min. Too high flow rate will raise the pressure at the inlet of COLUMN abnormally, deteriorating the performance of Separation. On the other hand, too low flow rate will lengthen the analyzing time required. Therefore, it is necessary to keep the flow rate within the range of 0.5 to 1 mL/min mentioned above so as maintain the due performance of COLUMN for a long period.

NOTE: The flow rate of 0.5 mL/min will make the pressure at the inlet of COLUMN of approximately 2 MPa; the flow rate of 1 mL/min will make the pressure of approximately 4 MPa, varying more or less depending on the temperature of COLUMN and each COLUMN.

● Temperature:

Although the COLUMN generally tends to exhibit its performance for separation better in higher temperature, resulting in capability of proceeding with samples more, it is recommended to use it below 60 °C in temperature for prevention of any deterioration of the COLUMN.

■ MOUNTING AND DISMOUNTING OF
COLUMN ON AND FROM SYSTEM


Clean the flowing channel of System well with water and replace with water completely prior to mounting the Shim-pack SCR-101N on the System.

As the proper FLOW direction of the mobile phase is indicated on the COLUMN, never to flow it in the reverse direction.

The column is connected with the following parts. Endure that the fittings are connected properly to avoid creating dead volume between the tubing and the column interface.

Item Name	P/N	Comments
Ferrule 1.6F 316L	228-16000-10	1/pkg
Male nut 1.6MN	228-16001	1/pkg

Observe the appropriate connection method as illustrated in the Figure.

column connection

Do not overtighten as this will result in damage to the fittings.

Never to change the pressure abruptly. It is required to change the pressure gradually.

Take care not to allow the mobile phase kept in the storage bath to be used up during the operation for analysis.


Prior to stopping the operation of System, make sure Than all the ingredients constituting the sample introduced into the COLUMN have been eluted, If the analysis has been made at the temperature of COLUMN above room temperature, it is necessary to deliver the mobile phase continuously in the same condition as that during the analysis until the temperature of COLUMN has returned to room temperature.

In case of making the analysis successively tomorrow, too, leave the System as it is.

In case the analysis will not be made for coming 2 to 3 days, disconnect the connection between the outlet of COLUMN and the detector and then close up to the outlet of COLUMN with the stop plug to prevent the inside of COLUMN from drying. In case the COLUMN will not be used for coming fairly long period, it is required to flow the mobile phase through it for about one hour at least once every 4 weeks to prevent for decomposing of filled water. The column is delivered filled with water. If the column will not use for long time, it is also required periodic cleaning with fresh water mentioned above.

In case the dismantling the COLUMN from the System, close up the both ends of it with the Stop plug immediately after its removing to prevent its inside from drying.

Store the COLUMN in a place which is kept from direct sunshine and is little subject to temperature changes. In cold season, it is recommended to store it indoors, taking care not to let the packings kept in the COLUMN freeze. In hot season, it is recommended to store it in the place which temperature is under 40 °C.

Typical Example of Chromatogram showing the Performance of COLUMN

■ TECHNICAL SUPPORT

Shim-pack SCR-101N columns are manufactured, inspected, packaged and shipped under strict standards of quality control. Should you find any defect in performance, please contact your local Shimadzu representative, who will ensure your complete satisfaction.

We regret that we cannot guarantee the lifetime of columns, also that we cannot accept any claim when performance has deteriorated due to noncompliance with the operation procedures elucidated above, or as a result of normal aging.

※ The contents of this instruction sheet are subject to change without notice.

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