High Performance Packed Column for HPLC

CoreFocus

Shim-pack SR-C30 [Reversed phase column]

INSTRUCTION MANUAL

■ Introduction

This is the instruction manual for the Shim-pack SR-C30. Shim-pack SR-C30 is a fully porous silica column. To ensure stable and long-term use of the Shim-pack SR-C30, which has excellent column performance, please read this instruction manual carefully and use it correctly.

■ Operating Precautions

- Please check for any abnormalities in the appearance and packaging box of the column.
- Verify that there are no errors in the name of the packing material, particle size, column size, etc., by checking the labels attached to the packaging box and the column.
- The inspection report includes the column lot number, column serial number, and chromatographic test conditions. Please keep the report for future reference.
- The number of theoretical plate stated in the inspection report is for quality checks at the time of shipment and does not indicate the performance (specifications) of the column.

■ Column Performance

Shim-pack SR-C30 undergoes rigorous QC testing to ensure quality and stability. The shipping solvent is Acetonitrile / Water = 60 : 40 (v : v). Refer to the table below for the specification of Shim-pack SR-C30.

Item	Specification		
Functional Group	Triacetyl group		
Pore Size	12 nm		
Surface Area	200 m ² /g		
Carbon Content	19%		
Column hardware	Stainless (SUS)		
Usable pH range	2-7.5		
Common temperature	20 - 40°C		
Maximum temperature	50°C		

P/N	Particle Size	Column I.D.	Column Length
227-31218-01	3 μm	4.6 mm	100 mm
227-31218-02			150 mm
227-31218-03	5 μm	4.6 mm	150 mm
227-31218-04			250 mm

NOTE

Near the pH limit of use, please use a mobile phase containing at least 10% organic solvent. At the pH limit, column lifetime may be shortened depending on conditions such as temperature and mobile phase composition. Additionally, continuing to use the column at the maximum temperature limit may also reduce its lifetime.

■ Column Installation

- The flow direction of the column is indicated by an arrow on the column label (→). When installing the column, ensure that the flow direction arrow matches the mobile phase flow direction.
- Use PEEK or SUS tubes with an inner diameter of 0.25 – 0.3 mm and an outer diameter of 1.6 mm. However, do not use PEEK tubes when working with fluorinated organic solvents such as tetrahydrofuran (THF), chloroform, concentrated sulfuric acid, concentrated nitric acid, dichloroacetic acid, acetone, dichloromethane, dimethyl sulfoxide (DMSO), and hexafluoroisopropanol (HFIP).
- Please use a male nut or fitting for connecting the column. When making connections, be careful to avoid creating any unnecessary gaps. If there are gaps at the connection points, it may lead to liquid leaks and a decrease in column performance (theoretical plates, peak symmetry). For information on connecting the column, please also refer to the system guide and the operating manual of the LC system.
- Male nuts or fittings can be obtained with the following product names and part numbers (P/N).

Item Name	P/N	Remarks	Pressure Limit
Male Nut, PEEK	228-18565-84	5 pcs	20 MPa
Male Nut Fitting Kit	228-45717-01	2 pcs	35 MPa
UHPLC Fitting 2 S	228-56867-41	1 pc	130 MPa

NOTE

Particulates or air in the system flow line may deteriorate the column. Before connecting the column, please flush the flow line with the mobile phase to clean it.

- If peaks are tailing more on the early eluting compounds than later eluting compounds, it's possible that there is void volume. Make sure the connection tube is fully inserted into the column joint. As a guideline, tighten the connection when the tube extends about 5 mm from the tip of the ferrule, and when removed, the tubes should extend 2-3 mm from the tip of the ferrule.

- Make sure to use appropriate internal diameter and tube length from the injector and to the detector to reduce system dead volume and peak broadening.

■ Column Handling Precautions

- To avoid damage to the column and prevent deterioration in performance, do not drop or bump the column.
- To maximize column lifetime, use the column within the pressure limit shown in the following table.

Particle Size	Column I.D.	Column Length	Maximum Pressure
3 μm	4.6 mm	100 mm	20 MPa
		150 mm	30 MPa
5 μm	4.6 mm	150 mm	20 MPa
		250 mm	25 MPa

*Continuous use near the maximum pressure limit can shorten the column's lifetime. Pressure varies depending on factors such as column length, column temperature, and the type of organic solvent, so please adjust the flow rate accordingly.

- The column should be disconnected from the system only after the pressure shows 0 MPa.
- Sudden pressure fluctuations can lead to early degradation of the column. Please avoid slow operation of the sample injection valve.
- If there are issues such as poor reproducibility of compound retention times, baseline drift, or noise at the start of the analysis, it may indicate that the column has not been adequately equilibrated. It is recommended to flush the column with at least five times its volume using the mobile phase intended for the analysis until improvements are observed. If no improvements are seen, please consider reviewing the equipment and analytical conditions.

NOTE

The elution order, retention time, peak shape, etc. may change significantly when ion-pairing reagents are used in the mobile phase. These reagents are often difficult to completely remove, so we recommend that columns with a history of using ion-pair reagents be used exclusively for analyses with ion-pair reagents.

■ Mobile Phase Selection

- Both aqueous and non-aqueous solvents can be used as mobile phases; however, repeatedly switching between solvents with extremely different polarities may lead to a decline in column performance. Common organic solvents that can be used include acetonitrile, methanol, and tetrahydrofuran (THF). However, when using THF, please pay attention to the solvent resistance of the materials used for the tubes and other components in the flow path.

NOTE

When switching from aqueous solvents to non-aqueous solvents, it is recommended to flush with an intermediate polarity solvent, such as isopropanol or dioxane, at least ten times the column volume (approximately 25 mL for a 150 $\times$ 4.6 mm I.D. column) before making the substitution. Additionally, since intermediate polarity solvents may result in higher pressure, please be cautious and reduce the flow rate by half during the flushing process to manage the pressure.

- Conditions such as pH and temperature extremes affect the column lifetime. Typical operating temperature is between 20°C and 40°C.

■ Sample Considerations

Samples should be dissolved in a solvent with a composition as similar as possible to that of the mobile phase (initial solvent during gradient elution). Injecting a large volume of samples dissolved in a solvent with higher solubility than the mobile phase can lead to decreased separation efficiency or precipitation of the sample at the column inlet. Additionally, to prevent the precipitation of salts contained in the sample or the sample solvent within the column, please confirm their miscibility with the mobile phase before injection.

■ Column Clogging and Ghost Peaks

The causes of pressure rise and peak splitting may include clogging or contamination of the filter at the column inlet.

- Filter the mobile phase using a 0.45 μm membrane filter before flushing the flow line and connecting the column.
- Please filter the sample solution using a membrane filter (0.2 - 0.45 μm) before injection. Using a guard column* can help prevent clogging of this column.

*The guard column is a cartridge-type column sold separately from the analytical column. Additionally, a separate cartridge holder is required when using the cartridge column. Cartridge columns and holders can be obtained using the following product numbers.

P/N	Item	Particle size	Column I.D.	Column Length
227-31218-05	Cartridge columns	3 μm	4.0 mm	10 mm
227-31218-06				20 mm
227-31218-07		5 μm	4.0 mm	10 mm
227-31218-08				20 mm
227-31172-03	Holders	-	-	10 mm
227-31172-04		-	-	20 mm

- Possible causes of baseline drift and noise include pump malfunction due to dissolved air, decreased light intensity when using a UV detector, bubble formation during high-temperature analysis, and solvent purity. If baseline drift and noise occur and do not improve after flushing the column, please consider reviewing the equipment and analytical conditions.
- If ghost peaks caused by system-derived contaminants such as mobile-phase-related impurities are detected, installing a Ghost Trap DS* between the pump and the injector may help remove the ghost peaks. The Ghost Trap DS is available under the following product name and P/N.

Item	P/N	Description	Size	Pressure Limit
Ghost Trap DS	228-59921-92	Two cartridges and one holder	7.6 mmID ×30 mm	35 MPa
	228-59921-94	Two cartridges and one holder	4.0 mmID ×20 mm	
Ghost Trap DS-HP	228-59931-91	Packed type (1 pc)	2.1 mmID ×30 mm	100 MPa

*When using a mass spectrometer as a detector, bleed noise originating from this product may occur; therefore, please do not use it. Additionally, in analyses using ion-pair reagents, the ion-pair reagents may be retained by this product, potentially affecting retention time and peak shape.

■ Cleaning and Storing the Column

In general, columns should be cleaned using the following methods. The flushing volume during cleaning and replacement should be at least five times the column volume.

- If the mobile phase does not contain buffers or salts, increase the concentration of the organic solvent that makes up the mobile phase to wash away substances with significant retention remaining in the column. Organic solvents can be used up to 100%. If highly lipophilic components are adsorbed, adding THF may be effective.
- If the mobile phase contains buffers, salts, or ion-pairing agents, first replace it with a water/organic solvent mixture (in the same ratio as the mobile phase) that does not contain these components, and then clean as described above. For buffers or salts around 50 mmol/L, direct replacement with a 60% acetonitrile aqueous solution is possible. During this process, take care to ensure that salts do not precipitate.
- After using the column near the pH limits, washing with water alone may cause column degradation. It is recommended to clean with the aforementioned water/organic solvent mixture or a 60% acetonitrile aqueous solution.
- If macromolecules such as proteins or polysaccharides are adsorbed onto the column, they are generally difficult to remove by cleaning. For samples containing these substances or samples with many impurities, it is recommended to perform pre-treatment in advance.
- During long-term storage, replace the shipping solvent with Acetonitrile/Water = 60:40 (v:v). When storing, seal with the provided plug and keep it in a cool, clean, dark place with minimal temperature fluctuations and low humidity.

■ Disposal Precautions

When disposing of the column, do so in accordance with the processing standards determined by law, separately from general industrial waste and household garbage.

■ Technical Support

Shim-pack SR-C30 columns are manufactured, inspected, packaged and shipped under strict standards of quality control. Should you find any defect in performance, please contact your local Shimadzu representative, who will ensure your complete satisfaction.

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