


Fundamentals of HPLC

2 _ Role of column and Separation Modes in HPLC

Shimadzu Corporation
Analysis and Measurement Division
Global Application Development Center
Content : Natsuki Iwata (May,2021)


Column and mobile phase in the HPLC configuration



Factors affecting separation in HPLC

- **Separation mode**

- Separation of target compounds from co-existing matrix are accomplished using various interactions provided from the combination of the column and the mobile phase.
- Separation mode mainly depends on the column chemistry
- Mobile phase basically consist of two solvents that have strong and weak elution power to target compounds respectively
(Except size exclusion mode)



Typical separation modes


- **Reversed phase chromatography (Reversed phase mode)**
- **Normal phase chromatography (Normal phase mode)**
- **Ion exchange chromatography (Ion exchange mode)**
- **Size exclusion chromatography (Size exclusion mode)**

Polarity and hydrophobicity of the compounds affecting separation


What is polarity ?

- • • In physical chemistry, it's defined by the dipole moment of compound.

e.g. Water is a **polar** compound,
Methane is a **non-polar** compound.



- • • In chromatography, compounds with highly hydrophilic functional groups are called as polar compounds (**Hydrophilic**).



e.g. **Hydrophilic groups** : hydroxyl group, carboxyl group, amino group (**High polarity**).
Hydrophobic groups : alkyl chains (methyl groups etc.), benzene rings (**Low polarity**)

In reversed phase analysis, **compounds with higher hydrophilic groups elute faster** (Compounds with higher hydrophobic groups elute slower).

-OH Hydroxy group (Hydrophilic group)
-CH₃ Methyl group (Hydrophobic group)

Normal phase and Reversed phase

	Stationary phase	Mobile phase
Normal phase	High polarity (hydrophilic)	Low polarity (hydrophobic)
Reverse phase	Low polarity (hydrophobic)	High polarity (hydrophilic)

Reversed phase chromatography

● Mechanism of retention

- Separation mode using hydrophobic interaction

● Features


- Wide range of applications
 - >Applicable to most of the compounds
- High separation efficiency
 - >Easy to get sharp peaks
 - >Analysis time can be short.
- Neutral hydrophilic compounds are poorly retained
 - >Not suitable for carbohydrates such as sugar

● Typical mobile phase

- Weak elution power: Water or Buffer
 - Adjust the retention of ionic compounds by pH
- Strong elution power : Organic solvents such as methanol and acetonitrile
 - Adjust the retention by concentration of organic solvent

● Typical columns

- ODS (Octadecyl: C18), C8 (Octyl), TMS (Trimethyl), Phenyl, etc.



LAAN-E-LC 215



7

7

Typical analysis of reversed phase chromatography

● Separation by difference in hydrophobicity

- Lower hydrophobic compounds elute faster



Analytical conditions

Column: Shim-pack VP-ODS (150 mm x 4.6 mm I.D.)
 Mobile phase: Water/acetonitrile = 30/70
 Flow rate: 1 mL/min
 Column temperature: 40 °C
 Detection: UV at 245 nm

LAAN-E-LC 215

8

8


Normal phase chromatography

● Mechanism of retention

- Separation mode using hydrophilic interaction

● Features

- Normally water is not used as mobile phase
 - > Suitable for hydrolyzable compounds
 - > Convenient for preparative purification
- Saturated hydrocarbon compounds can not be retained
- Basic compounds adsorbed strongly
- High structural recognition ability
 - > Suitable for the separation of geometrical isomers



● Typical mobile phase

- Solvents with weak elution power: Low Polarity solvents : hexane
- Solvents with strong elution power: Highly polarity solvents : ethanol and 2-propanol
 - > Adjust retention by concentration of polarity organic solvent
 - > Adjust retention of ionic compounds by adding acid/base
- Mixture of acetonitrile and water may be used in particular cases (hydrophilic interaction chromatography, HILIC).

● Typical columns

- SIL (silica gel), CN (Cyano), NH₂ (Amino)


Ion exchange chromatography

● Mechanism of retention

- Separation mode using ion interaction based on ionic bond

● Features

- Suitable for all of the ionic compounds
 - > Suitable for both organic and inorganic compounds
- Theoretically neutral compounds can not be retained
- Cations and anions can not be separated in one analysis (Columns can retain one but can not retain the other)



● Typical mobile phase

- Aqueous solution (Mixing with acetonitrile, etc. if necessary) containing salts such as a buffer solution
 - > Higher concentration (ionic strength) buffer results in faster elution
 - > In cation exchange mode, increase the pH to elute faster
 - In anion exchange mode, decrease the pH to elute faster.

● Typical columns

- Strong cation exchange (SCX), weak cation exchange (WCX)
- Strong anion exchange (SAX), weak anion exchange (WAX)

Size exclusion chromatography (SEC)

- **Mechanism of retention**

- Separate molecules by their solvated size (molecular weight)

- **Features**

- Separation is only depends on the size of molecules
->Larger molecules elute faster
- Theoretically, all molecules should elute between the exclusion limit and penetration limit of column
-> Allowable molecular weight range is determined by the characteristics of column

- **Typical eluent**

<Non-aqueous mode, GPC>


- Tetrahydrofuran, chloroform, dimethylformamide, etc.

< Aqueous mode, GFC>

- Water, buffer

- **Typical columns**

- Elution time can be controlled by changing columns that have particular molecular weight range respectively
- Non-aqueous type (organic solvents) and aqueous type (water or aqueous buffer solution)



LAAN-E-LC 215

11

11

Thank you for your attention.

The content of this video is based on the latest information available at the time of publication.

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, Whether or not they are used with trademark symbol "TM" or "®".

12